Instructions for Use

Leukocyte Alkaline Phosphatase Kit

Procedure No. 85

Intended Use

Sigma-Aldrich Alkaline Phosphatase kits are intended as a general laboratory use histological stain. Alkaline Phosphatase reagents are for professional "In Vitro Diagnostic Use" only. This manual, qualitative, histochemical procedure demonstrates alkaline phosphatase activity in leukocytes of human blood or bone marrow films. Histological visualization of leukocyte associated alkaline phosphatase is a unique technique currently commonly used in medicine.

In 1929, Kay first suggested the presence of alkaline phosphatase in leukocytes.¹ However, not until many years later did Kaplow introduce a practical staining method for demonstrating the leukocyte enzyme.² Kaplow used sodium α -naphthyl phosphate as substrate, fast blue RR as the diazonium salt and recommended a scoring method for semiquantitative comparisons.³ Also suggested were improvements resulting in good preservation of cellular morphology. Subsequently, α -naphthyl phosphate was replaced by naphthol AS phosphate. It is noteworthy that a good correlation has been found between cytochemical and biochemical techniques for assaying leukocyte alkaline phosphatase activity (LAPA).³

With the Sigma-Aldrich method, leukocyte alkaline phosphatase is determined in blood or bone marrow films after they are gently fixed to slides. The fixed films are then incubated in a solution containing naphthol AS-MX phosphate. As a result of phosphatase activity, naphthol AS-MX is liberated and immediately coupled with a diazonium salt, forming an insoluble, visible pigment at sites of phosphatase activity.

Reagents

Fast Blue RR Salt (Cat. No. FBS25-10CAP)

Fast Blue, C.I. 37155, Preweighed capsules. Actual weight per capsule will vary with dye lot and has been optimized by assay.

Naphthol AS-MX Phosphate Alkaline Solution (Cat. No. 855-20ML)

Naphthol AS-MX phosphate, 0.25% (w/v), buffered at pH 8.6, 25°C. Danger. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapors/spray. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.

Mayer's Hematoxylin Solution (Cat. No. MHS1-100ML)

Hematoxylin, certified, 1 g/L, C.I. 75290, sodium iodate, 0.2 g/L, aluminum ammonium hydroxide, 5 g/L, and stabilizers

Citrate Concentrated Solution (Cat. No. 854C-20ML)

Citric acid-sodium citrate, 1.5 mol/L. Danger. Causes skin irritation. Causes serious eye damage. Avoid breathing dust/fume/gas/mist/vapors/spray. May cause respiratory irritation. Wear protective gloves/eye protection/face protection.

Special Materials Required but Not Provided

- Scott's tap water substitute concentrate, Catalog No. S5134
- Acetone, ACS reagent

Storage and Stability

Store Fast Blue RR Salt below 0°C. Store Naphthol AS-MX Phosphate Alkaline Solution in refrigerator (2–8°C). Reagents are stable until expiration date.

Store Mayer's Hematoxylin Solution tightly capped at room temperature (18–26°C) protected from light. Do not return solution to original container after use in Coplin jar. Discard when the time required for suitable staining exceeds the time recommended in the procedure by more than 5 minutes.

Store Citrate Concentrated Solution at room temperature (18–26°C). Store Citrate Working Solution in refrigerator (2–8°C). Citrate solutions are suitable for use in the absence of microbial growth.

Preparation

Fast Blue RR Salt, Naphthol AS-MX Phosphate Alkaline Solution and Mayer's Hematoxylin Solution are provided ready for use in the procedure.

Prepare Citrate Working Solution by diluting 2 mL Citrate Concentrate Solution to 100 mL with deionized water.

To prepare Fixative Solution (citrate buffered acetone, 60%), warm Citrate Working Solution to room temperature (18–26°C). With constant stirring, add 2 volumes of Citrate Working Solution to 3 volumes of acetone. Discard after use.

Scott's Tap Water Substitute Working solution is prepared by diluting 1 volume of Scott's Tap Water Substitute Concentrate with 9 volumes deionized water.

Precautions

These IVD's are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

Fresh blood or bone marrow films or samples anticoagulated with heparin may be used.⁴ AVOID EDTA.⁴ Blood smears should be stained for enzyme activity within 8 hours after preparation. However, if this is not possible, gradual loss of alkaline phosphatase activity may be delayed by fixation and storage overnight in freezer.⁴ Films should be dried at least 1 hour prior to fixation and 3 hours post-fixation before freezing.

Notes

There is a scarcity of data concerning compounds which may interfere with leukocyte alkaline phosphatase activity (LAPA). Certain drugs and other substances are known to influence circulating alkaline phosphatase activity.⁵ Oral contraceptives, cortisol and stress may result in elevated leukocyte alkaline phosphatase scores.⁶

Perform procedure using positive controls. These can be obtained from patients with pyogenic leukocytosis, or women in the third trimester of pregnancy or during the first several days postpartum. Leukocyte alkaline phosphatase scores from these persons usually exceed 100. A negative control can be prepared from a normal fixed smear by immersing it in boiling water for 1 minute to inactivate the enzyme. Control films can be preserved up to 1 year if stored fixed, wrapped in Parafilm® at –70°C. These films should be air dried at least 1 hour prior to fixation and 3 hours post-fixation before freezing.

It is strongly recommended that each laboratory establish its own expected range, characteristic for the local population.

The procedure depends upon subjective rating of stained cells. This can result in a wide variation of ratings obtained. The temperature of the reaction mixture must be kept between 18–26°C. Lower temperatures will result in significantly lower scores. Above 30°C, marked increases in activity will occur. Eosinophils do not stain but can be recognized by bilobate nuclei and large refractile granules.

Procedure

1. Measure 48 mL distilled water into a suitable container and adjust temperature to 18–26°C.
2. Prepare diazonium salt solution:
Dissolve contents of one Fast Blue RR Salt capsule in distilled water (from Step 1). A magnetic stirrer may be helpful.
3. Add 2 mL Naphthol AS-MX Phosphate Alkaline Solution to diluted diazonium salt solution (from Step 2). Mix.
4. Bring Fixative Solution to room temperature (18–26°C). Fix slides by immersing in Fixative Solution for 30 seconds. Rinse gently in deionized water for 45 seconds. Do not allow slides to dry.
5. Add slides to alkaline-dye mixture (from Step 3) and incubate at 18–26°C for 30 minutes. Protect immersed slides from direct light. Discard alkaline-dye mixture after use.
6. After 30 minutes, remove slides and rinse thoroughly in deionized water for 2 minutes. Do not allow slides to dry.
7. Place slides in Mayer's Hematoxylin Solution for 10 minutes.
8. Rinse counterstained slides for 3 minutes in deionized water. This will result in red violet nuclear staining.
9. Evaluate microscopically. If coverslipping is required use only an aqueous mounting media.

Performance Characteristics

Method of Scoring

Scan the film (900X) and select a thin area where erythrocytes are barely touching. Sites of phosphatase activity will appear as blue granules. Select 100 consecutive segmented and band form neutrophilic granulocytes. Rate from 0 to 4+ on the basis of quantity and intensity of precipitated dye within the cytoplasm of these cells. For characteristics of scoring, refer to Table 1. The sum of the ratings of 100 cells is regarded as the score.

Table 1. Characteristics of Scoring*

Cell Rating	Precipitated Azo Dye in Cytoplasm			Background of Cytoplasm
	Amount** (%)	Granule Size	Staining Intensity	
0+	None	-	None	None
1+	50	Small	Faint to Moderate	Colorless to very pale blue
2+	50-80	Small	Moderate to Strong	Colorless to pale blue
3+	80-100	Medium to Large	Strong	Colorless to blue
4+	100	Medium to Large	Brilliant	Not visible

* Table 1 represents modification of observations made by Kaplow.^{2,3}

** Percentage of volume of cytoplasm occupied by azo dye precipitate.

To obtain the leukocyte alkaline phosphatase activity (LAPA) score, the number of cells counted is multiplied by the value for cell rating. These figures are added to obtain the LAPA score as shown in the following examples:

Cell Rating	Counted Score	LAPA Score
0	60	0
1+	20	20
2+	14	28
3+	5	15
4+	1	4
Total	100	67

Expected Values

The following scores were obtained at Sigma-Aldrich for 20 normal individuals. Blood was drawn into heparinized tubes and films were prepared within one hour.

Fast Blue RR	
LAPA Mean $\pm$ 1SD	96 $\pm$ 44
LAPA Range	52-140

The following range of scores were obtained from those 20 normal individuals:

Fast Blue RR	
	32-182

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Expected Observations

In humans, alkaline phosphatase activity is restricted to mature and band-form granulocytes. Occasionally weak staining may be observed in lymphocytes. Bone marrow osteoblasts and endothelial cells stain strongly.

A series of slides was prepared using blood obtained from normal subjects, females in the last trimester of pregnancy, early postpartum patients, and individuals exhibiting a leukemoid reaction. Negative controls were prepared by heat inactivation as described. Various combinations of Fast Blue RR and counterstains (Methylene Blue Solution and Mayer's Hematoxylin Solution) were performed in replicate and the LAPA scored. Values obtained on multiple slides prepared on each subject were in close agreement and did not appear influenced by the staining/counterstaining technique employed. For example, 7 smears made using blood from a prenatal patient yielded scores ranging from 235–269 with a mean, standard deviation and coefficient of variation of 255, 11.3 and 4.4% respectively. Negative controls revealed a complete absence of stained material.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
FBS25	Fast Blue RR Salt	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
855	Naphthol AS-MX Phosphate Alkaline Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
MHS1	Mayer's Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
854C	Citrate Concentrated Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

85L1:

- H290 May be corrosive to metals.
- H331 Toxic if inhaled.
- P234 Keep only in original packaging.
- P261 Avoid breathing mist or vapors.
- P271 Use only outdoors or in a well-ventilated area.
- P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.
- P390 Absorb spillage to prevent material damage.
- P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

85L2:

- H290 May be corrosive to metals.
- H319 Causes serious eye irritation.
- H331 Toxic if inhaled.
- P234 Keep only in original packaging.
- P261 Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.
- P264 Wash skin thoroughly after handling.
- P271 Use only outdoors or in a well-ventilated area.
- P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No 14, National Committee for Clinical Laboratory Standards
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248–250

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 2.0	2014
Rev. 3.0	2016
Rev. 4.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Removed reference #4 and renumbered the references. Added a revision history table. Added additional danger statement for 854C-20ML and 855-20ML. Added Warnings and Hazards. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Vorgeschlagene Norm: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No 14, National Committee for Clinical Laboratory Standards
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248-250

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Références

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No 14, National Committee for Clinical Laboratory Standards
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248-250

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Bibliografia

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929.
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955.
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968.
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No 14, National Committee for Clinical Laboratory Standards.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581.
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248-250.

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Referenser

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, Nr. 14, National Committee for Clinical Laboratory Standards
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, sid. 248–250

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Referências

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No. 14, National Committee for Clinical Laboratory Standards
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248-250

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Reference

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No 14, National Committee for Clinical Laboratory Standards
- Young DS. Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6): 579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248-250

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Referanser

- Kay HD: Plasma phosphatase in osteitis deformans and in other diseases of the bone. Br J Exp Pathol 10:253, 1929
- Kaplow LS: A histochemical procedure for localizing and evaluating leukocyte alkaline phosphatase activity in smears of blood and marrow. Blood 10:1023, 1955
- Kaplow LS: Leukocyte alkaline phosphatase cytochemistry: Applications and methods. Ann NY Acad Sci 155:911, 1968
- Proposed Standard: Histochemical Method for Leukocyte Alkaline Phosphatase H22-P, Vol 4, No 14, National Committee for Clinical Laboratory Standards
- Young DS: Effects of Drugs on Clinical Laboratory Tests. Annals of Clinical Biochemistry. 1997;34(6):579-581
- Elias JM: Principles and Techniques in Diagnostic Histopathology, Noyes Publications, New Jersey, 1982, pp 248–250

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MDSS CH GmbH
MDSS GmbH