

User Guide

### JumpStart™ *Taq* DNA Polymerase

With 10X reaction buffer and a vial of 25 mM MqCl<sub>2</sub>

D4184

### **Product Description**

JumpStart<sup>™</sup> *Taq* DNA Polymerase is an optimized blend of our high-performance *Taq* DNA Polymerase and JumpStart<sup>™</sup> *Taq* antibody. The neutralizing monoclonal antibody binds to *Taq* DNA polymerase, inactivating it until the complex is dissociated at  $\geq$ 70 °C, providing simple and efficient hot start PCR. Because JumpStart<sup>™</sup> *Taq* DNA Polymerase is inactive at room temperature, reaction mixtures can be prepared on the bench and transported directly to the thermocycler as the enzyme will become active in the first denaturation step of the cycling process.

The enzyme is supplied at 2.5 units/ $\mu$ L and comes with an optimized 10X reaction buffer free of magnesium chloride (MgCl<sub>2</sub>). A separate vial of MgCl<sub>2</sub> is provided to allow optimization of Mg<sup>2+</sup> concentration for best reaction efficiency.

### **Features**

- Reduces nonspecific amplification and primer-dimer artifacts
- · Increases product yield
- Allows for room-temperature PCR set up
- Suitable for large number of thermal cycles (>35)
- Average amplification length is 0.1 to 3 kb

### **Applications**

JumpStart™ *Taq* can be used for amplification of:

- Genomic DNA
- cDNA
- Low copy number targets
- Multiple targets (multiplex PCR)

### Unit Definition

One unit incorporates 10 nmol of total deoxyribonucleoside triphosphates into acid precipitable DNA in 30 minutes at 74 °C.

### Reagents Provided

- JumpStart™ Taq DNA Polymerase,
   2.5 units/µL in 50% glycerol, 87.5 mM KCl,
   17.5 mM Tris-HCl, pH 8.0, 75 µM EDTA,
   7.5 mM DTT, and stabilizers
- 10X PCR Buffer without MgCl<sub>2</sub>, 100 mM Tris-HCl, pH 8.3, 500 mM KCl
- MgCl<sub>2</sub> solution, 25 mM

### Materials and Reagents Required

(Not included, see Product Ordering)

- Deoxynucleotide (dNTP) Mix, containing 10 mM each of dATP, dCTP, dGTP, and dTTP sodium salts
- Nuclease free water
- Custom ordered primers specific to gene target
- PCR tubes or plates
- Sample containing template DNA
- Thermal cycler

### Precautions and Disclaimer

This product is for R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

### Storage/Stability

1

Store at -20 °C. JumpStart<sup>TM</sup> Taq DNA Polymerase is stable for 2 years in the supplied storage buffer and at the supplied concentration. Storing JumpStart<sup>TM</sup> Taq DNA Polymerase below -20 °C is not recommended, as repeated freeze-thaw cycles of this product may adversely affect its function.



### Directions for Use

For best reproducibility, assemble a master mix of PCR reagents by multiplying the number of reactions needed (plus 10% to account for pipetting error) by the suggested volumes in the table below. Aliquot reaction mixture into PCR tubes.



Reagent	<b>Final Concentration</b>	Amount per 20 µL reaction
10X PCR buffer	1X	2 μL
$MgCl_2$	1.5-7 mM	Variable
JumpStart™ Taq	0.05 U/μL	0.4 μL
dNTP Mix, 10 mM	200 μΜ	0.4 μL
Primers	0.1-0.5 μΜ	Variable
Template	Variable	Variable

**Note:** JumpStart<sup>TM</sup> Taq DNA Polymerase is a magnesium ion-dependent enzyme, optimal concentrations of template DNA, primers, and MgCl<sub>2</sub> will be target-specific. For more information on ideal concentrations based on application see <u>Technical Guide</u>.



Recommended input template is 10 ng DNA; however, JumpStart™ *Taq* may amplify as little as a single copy of non-complex template or 10-100 copies of complex genomic template. For cDNA templates use a 1:10 reaction dilution for medium to highly expressed targets, or a 1:2 to 1:5 dilution for low expression targets.

A suggested thermocycling protocol using JumpStart<sup>™</sup> *Taq* DNA Polymerase is provided below:



Initia	I denaturation	94 °C	2 min
S	Denaturation	94 °C	15 seconds
40 cycles	Annealing	60 °C or 5 °C below lowest primer $T_{\text{m}}$	30 seconds
4	Extension	72 °C	1 min/kb
Final	extension	72 °C	1 min
Hold		4 °C	<sub>∞</sub>

Amplification parameters will vary depending on primers, template, and instrument used. For tips on optimizing PCR conditions as well as a 2-step cycling protocol please see Technical Guide.



D4184ug Rev 11/22

### Technical Guide

### Considerations for Primer Design

Thoughtful primer design is essential for PCR efficiency and specificity. For successful amplification consider the following:

- Select an 18-30 nucleotide-long sequence with 40-60% G/C content and even distribution of all 4 bases.
- Avoid repetitive elements or self-complementary sequences >3 bp.
- Primer pairs should not differ in length by >3 bp and should not contain complementarity to one another.
- Maintain calculated primer T<sub>m</sub> between 55-60 °C, permitting only 2-3 °C variation between primer pairs.
- Priming efficiency can be increased by including a terminal G at the 3' end; however, the number of Gs or Cs in the last 5 bases of the primer sequence should be no more than 3.
- Ensure each primer sequence is unique to the gene of interest and is absent in other genes in the gDNA sample or within the vector.

### Handling gDNA Templates

To prevent genomic DNA (gDNA) shearing, add template last and mix gently using a wide pore pipet tip. DO NOT VORTEX!

#### Optimization of PCR Conditions

PCR involves the cycling of denaturing, annealing, and extension steps for DNA synthesis by a polymerase enzyme. To obtain the best product yield and accuracy, each step must be optimized.

The **denaturing** step (94-96 °C) activates the JumpStart™ *Taq* DNA Polymerase and separates double-stranded DNA strands, making it accessible to primers.¹ The duration of this step should be long enough to denature DNA but not so long that it compromises *Taq* DNA polymerase integrity. High salt conditions, GC-rich (>55%) templates, and gDNA templates may require longer denaturation times and/or higher temperatures. For maximum retention of JumpStart™ *Taq* activity during thermocycling, use 94 °C for denaturation.

- The annealing temperature can be calculated by subtracting 5 °C from the lowest reaction primer T<sub>m</sub>. The annealing time should be long enough for the primer to anneal to the template but not too long for non-specific annealing to occur.<sup>1</sup>
- The optimal **extension** temperature for *Taq* DNA Polymerase is 72 °C; however, lower temperatures may be used for some reactions. Extension time depends on length and complexity of the target sequence. For complex templates, use 1 minute/kb, with 15 seconds added if the PCR product is >2 kb. Short or non-complex templates may be amplified with extension times of 30 seconds/kb.
- The number of cycles needed for amplification depends on the amount of template input, with higher amount of input requiring less cycling.¹ Generally, 25-30 cycles are sufficient to produce detectable product; however, low concentration templates may require up to 45 cycles.
- To maintain enzyme fidelity, or accuracy of nucleotide incorporation, limit the number of PCR cycles and use an equimolar concentration of each dNTP.<sup>1</sup> A low magnesium concentration is also important to maintain enzyme fidelity.<sup>1</sup>

### Multiplex PCR

When performing multiplex PCR, competition between products for reagents may occur. Consider adjusting the following for optimization:<sup>2</sup>

- Proportion of primer pair concentration: if a target sequence produces a relatively "weaker" signal, the amount of primer used may be increased to compensate. For sequences with low copy numbers, or high-complexity, primer concentration can be used at 0.3-0.5 μM.
- Primer concentration can also be decreased for target sequences producing "stronger" signal to achieve balance. For high copy number or low-complexity sequences, primer concentration can be used at 0.04-0.4 μM.
- dNTP: Perform a stepwise increase of dNTP to a concentration ≤400 μM. Keep MgCl<sub>2</sub> concentration constant for this optimization.
- PCR Buffer: Use 2X buffer concentration for the reaction instead of 1X.
- Supplement the reaction with additional MqCl<sub>2</sub>, or PCR-enhancing additives.

### Recommended MgCl<sub>2</sub> Concentrations By Application

Because JumpStart<sup>TM</sup> Taq DNA Polymerase is a magnesium ion-dependent enzyme, the optimal concentrations of template DNA, primers, and MgCl<sub>2</sub> will be target-dependent. The optimal MgCl<sub>2</sub> concentration is also dependent upon the intended application. See table below for recommended ranges of MgCl<sub>2</sub> to use in reactions containing JumpStart<sup>TM</sup> Taq DNA Polymerase.

Application	Recommended MgCl <sub>2</sub> Concentration range
Endpoint PCR	1.5-3.5 mM
SYBR® green-based qPCR	3-5 mM
Probe-based qPCR	4-7 mM

### Two-Step PCR Amplification

Application of a two-step PCR process is possible when the annealing and extension temperatures are similar.

Initia	al denaturation	94 °C	3 minutes
Cycles	Denaturation	94 °C	3 seconds
40 Cyc	Annealing/ extension	60 °C*	15-30 seconds
Fina	Extension**	72 °C	1 min
Hold		4 °C	∞

<sup>\*</sup>Consult primer  $T_m$  regarding temperature selection. Extension time is target dependent, with larger targets requiring more than the recommended time.

### PCR-Enhancing Additives

When optimizing PCR conditions for a new experiment, the following can be added to the reaction mix individually. After performing PCR amplification, samples with and without additive can be compared using agarose gel electrophoresis or other standard methods to look for improved product specificity and yield.

Additive	Purpose
BSA (10-100 μg/mL)	<i>Taq</i> DNA polymerase stabilization <sup>3</sup>
Formamide (1.25-10%)	Increases specificity in G/C rich regions <sup>4</sup>
DMSO (Up to 5%)	Accelerates strand renaturation <sup>5</sup> Nucleic acid thermal stability against depurination <sup>5</sup>
Glycerol (Up to 10%)	Increases thermal stability of the polymerase and lowers the temperature necessary for strand separation <sup>5</sup>
Ammonium sulfate (15-30 mM)	Affects the denaturing and annealing temperatures of the DNA <sup>6</sup>
Single strand binding protein (0.7-1.5 µg)	Inhibits formation of secondary structures, improving fidelity and <i>Taq</i> processivity <sup>7</sup>
Betaine (0.8-1.6 M)	Reduces base pair composition dependence of DNA melting <sup>8</sup>

D4184ug Rev 11/22 4

<sup>\*\*</sup>The final extension step is needed when the objective is to use the PCR product for downstream applications (for example, cloning, or agarose gel electrophoresis), and is otherwise optional (for example, for qPCR).

### Troubleshooting Guide

#### **Problem**

### Suggestions

- Titrate MgCl<sub>2</sub> concentration in 0.5 mM increments using the provided supplemental MgCl<sub>2</sub>. See "Recommended MgCl<sub>2</sub> concentrations by application" section in the <u>Technical Guide</u> for expected concentration ranges based on application. Each amplicon target must be optimized individually.
- Adjust the annealing temperature in 2–3 °C increments or use a gradient PCR to find the optimal annealing temperature.
- Increase the number of amplification cycles. If currently using 25-30 cycles, increase the cycle number to 35-40.
- For complex templates like human genomic DNA, increase the initial denaturation time by 1-2 minutes and/or increase the denaturation temperature to 95 °C to overcome denaturation difficulties.

## No or low product amplification

- Check concentration of input template. For complex templates like intact eukaryotic genomic DNA, 1000 genome copies may be required for amplification of difficult targets. For highly concentrated templates, such as purified plasmid, consider diluting 1:1000 to improve amplification.
- Assess DNA quality to ensure absence of PCR inhibitors in sample. If presence of inhibitors is suspected, DNA can be diluted 1:10-1:100. Alternatively, lysis and DNA purification can be performed using the GenElute™ genomic DNA miniprep kits.
- Refer to "PCR-Enhancing Additives" section of the <u>Technical Guide</u> to improve amplification.
- If yield is too low for downstream applications, increase the reaction volume to 50–75  $\mu L$ .

# Amplification of nonspecific product(s)

- Raise the annealing temperature in 2–3 °C increments or use a gradient PCR to find
  the optimal annealing temperature. Raising the temperature improves the specificity
  of binding by the primers; however, it may also result in reduced binding and
  extension of the primers.¹ If raising the annealing temperature causes reduced yield of
  the specific product without eliminating side reaction products, it may be necessary to
  redesign the primers to improve specificity.
- Take precautions to avoid crossover contamination of PCR with both specific and nonspecific PCR products, including primer-dimer artifacts.<sup>9</sup>
- Titration of JumpStart<sup>™</sup> *Taq* may be necessary to optimize PCR efficiency, especially if the reaction conditions vary from those recommended in this document. In this case, increase the concentration of JumpStart<sup>™</sup> *Taq* by two- or four-fold. Increasing the concentration of JumpStart<sup>™</sup> *Taq* beyond this level may inhibit PCR.
- The use of more than 5% v/v DMSO with JumpStart™ *Taq* is not recommended as it may interfere with the enzyme-antibody complex. Other co-solvents, salts, and extremes in pH can also reduce the affinity of the JumpStart™ *Taq* antibody for the *Taq* DNA Polymerase and compromise its effectiveness for hot start PCR.

### References

- 1. Green MR, Sambrook J., Cold Spring Harb Protoc. 6, 436-456(2019).
- Markoulatos, P., et al., J. Clin. Lab. Anal. 16, 47-51 (2002).
- Paabo, S., et al., Nucleic Acids Res., 16, 9775-9787 (1988).
- Sarkar, G., et al., Nucleic Acids Res., 18, 7465 (1990).
- Cheng, S., et al., Proc. Natl. Acad. Sci. USA, 91, 5695-5699 (1994).
- Grunenwald H. PCR Protocols. Methods in Molecular Biology™. 226, 89-99 (2003).
- 7. Chou, Q., Nucleic Acids Res., 20, 4371 (1992).
- 8. Rees, W. A., et al., Biochemistry, 32, 137-144 (1993).
- 9. Kwok, S., and Higuchi, R., Nature 339, 237-238 (1989).

### **Notice**

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

### **Technical Assistance**

Visit the tech service page at SigmaAldrich.com/techservice.

### Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at <u>SigmaAldrich.com/terms</u>.

### **Contact Information**

For the location of the office nearest you, go to SigmaAldrich.com/offices.

### **Product Ordering**

Description	Number
Deoxynucleotide (dNTP) Mix, containing 10 mM each of dATP, dCTP, dGTP, and dTTP sodium salts	D7295
Nuclease-free water	W1754
Custom ordered primers specific to gene target	OLIGO
GenElute™-E Single Spin DNA Cleanup Kit	EC600
GenElute™ Bacterial Genomic DNA Kit	NA2120
GenElute™ Mammalian Genomic DNA Miniprep Kit	G1N10
GenElute™ Plant Genomic DNA Miniprep Kit	G2N70
GenElute™ PCR Clean-Up Kit	NA1020
GenElute™ Gel Extraction Kit	NA1111
Precast Agarose Gels	P6222 P5472 P6097 P5972 P5722
1 kb DNA Ladder	D0428
Water, Microbial DNA-free	MBD0025
Nuclease-Free Water, for Molecular Biology	W4502
JumpStart™ <i>Taq</i> Ready Mix	P2893
RED <i>Taq</i> ® Ready Mix	P0982
Glycerol-free JumpStart™ <i>Taq</i> DNA Polymerase	D9310
DMSO	D8418
Single strand binding protein	S3917
Betaine solution	B0300
Mineral Oil	M5904
Bovine Serum Albumin solution	B8667

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck, JumpStart, GenElute, REDTaq and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

