1.15973.0002 **REF****Microscopy****Weigert's iron hematoxylin kit**

for nuclear staining in histology

For professional use only

In Vitro Diagnostic Medical Device

**Intended purpose**

This "Weigert's iron hematoxylin kit - for nuclear staining in histology" is used for human-medical cell diagnosis and serves the purpose of the histological investigation of sample material of human origin. It is a ready-to-use staining kit that when used together with other *in vitro* diagnostic products from our portfolio makes target structures evaluable for diagnostic purposes (by fixing, embedding, staining, counterstaining, mounting) in histological specimen materials, for example histological sections of e.g. the liver, the kidney, the intestine, the placenta and similar.

Unstained structures are relatively low in contrast and are extremely difficult to distinguish under the light microscope. The images created using the staining solutions help the authorized and qualified investigator to better define the form and structure in such cases. Further examinations may be necessary to reach a definitive diagnosis.

Principle

Aluminium hematoxylin (e.g. nuclear fast red aluminium sulfate) are washed out of the tissue sections by acidic staining solutions. Iron hematoxylin solutions (e.g. Weigert's iron hematoxylin), on the other hand, are acid-resistant and are preferentially used for staining nuclei in histological tissue sections in which counterstaining is done with acidic staining solutions.

Examples here include the visualization of elastic fibers according to the van Gieson method (e.g. Elastica van Gieson) and the various trichrome connective-tissue staining methods (e.g. acc. to Masson-Goldner).

Commonly used hemalum solutions such as Mayer's Hemalum, Delafield, Gill and Harris tend to procedure also weak nucleic staining only when used with acidic counter-stains.

Sample material

Starting materials are sections of formalin- or Bouin-fixed tissue embedded in paraffin (3 - 5 µm thick paraffin sections).

Reagents

Cat. No. 1.15973.0002

Weigert's iron hematoxylin kit
for nuclear staining in histology**Package components:**

The staining kit contains

Reagent 1: Weigert's solution A - alcoholic hematoxylin solution	500 ml
Reagent 2: Weigert's solution B - hydrochloric acid iron(III)nitrate solution	500 ml

Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

When using the corresponding auxiliary reagents, the corresponding instructions for use must be observed.

Deparaffinize and rehydrate sections in the conventional manner.

Reagent preparation**Weigert's iron hematoxylin staining solution**

Mix reagent 1 and 2 (Weigert's solution A and Weigert's solution B) in the ratio 1 + 1.

The prepared staining solution remains stable for approx. one working week.

The solution must be exchanged as soon as the cell nuclei appear brown.

Procedure**Staining in the staining cell**

Deparaffinize histological slides in the conventional manner and rehydrate in a descending alcohol series.

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with histological specimen	
Weigert's iron hematoxylin staining solution	5 min
Running tap water	3 min
Counterstaining corresponding to the method used	
Ethanol 70%	1 min
Ethanol 70%	1 min
Ethanol 96%	1 min
Ethanol 96%	1 min
Ethanol 100%	1 min
Ethanol 100%	1 min
Xylene or Neo-Clear™	5 min
Xylene or Neo-Clear™	5 min
Mount the Neo-Clear™-wet slides with Neo-Mount™ or the xylene-wet slides with e.g. Entellan™ new and cover glass.	

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, histological slides can be covered with non-aqueous mounting agents (e.g. DPX new, Entellan™ new, Neo-Mount™) and a cover glass and can then be stored.

Result

Nuclei	blue-black
--------	------------

Trouble-shooting**Weak staining**

Weigert's iron hematoxylin staining solution (mixture of Reagents 1 and 2) can be used for approx. 1 working week after preparation; a new solution should be prepared after this time has elapsed (over-oxidation).

In addition, as an alternative it is also possible to fix the tissue specimen using the Bouin fixing method, since the formalin fixation may result in an insoluble protein network that in some tissues is capable of preventing the penetration of large staining particles (e.g. fuchsin), in particular when the tissue has been fixed for a very long time.

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using histoprocessor systems or automatic staining systems, please follow the instructions for use supplied by the supplier of the system and software.

Analytical performance characteristics

"Weigert's iron hematoxylin kit" stains and thereby visualizes biological structures, as described in the "Result" chapter of this IFU. The use of the product is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, histoprocessing, decisions regarding suitable controls and more.

The analytical performance of the product is confirmed by testing each production batch. The successful participation in international interlaboratory tests on a regular basis provide an additional and unaffiliated confirmation of analytical specificity.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100%:

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Nuclear staining				
Nuclei	8/8	8/8	6/6	6/6

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

The results of this Performance Evaluation confirms that the product is suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used. This method can be supplementarily used in human diagnostics. Further tests must be selected and implemented according to recognized methods. Suitable controls should be conducted with each application in order to avoid an incorrect result.

Storage

Store the Weigert’s iron hematoxylin kit - for nuclear staining in histology at +15 °C to +25 °C.

Shelf-life

The Weigert’s iron hematoxylin kit - for nuclear staining in histology can be used until the stated expiry date. After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C. The bottles must be kept tightly closed at all times. If stored at +15°C to +25 °C, the freshly prepared Weigert's iron hematoxylin staining solution can be used for minimum one working week. The solution must be exchanged as soon as the cell nuclei appear brown. However, the solutions should be discarded when contaminations (e.g. bacteria, fungi), that occur at times, are observed.

Capacity

The package is sufficient for 400 - 500 applications.

Additional instructions

For professional use only. In order to avoid errors, the application must be carried out by qualified personnel only. National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines. Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link “Hints for Disposal of Microscopy Products” at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No.	1.00199	Picrofuchsin solution acc. van Gieson for microscopy	500 ml
Cat. No.	1.00579	DPX new non-aqueous mounting medium for microscopy	500 ml
Cat. No.	1.00591	ELASTIN color solution acc. to Weigert for microscopy	500 ml
Cat. No.	1.00974	Ethanol denatured with about 1% methyl ethyl ketone for analysis EMSURE®	1 l, 2.5 l
Cat. No.	1.04699	Immersion oil for microscopy	100-ml dropping bottle, 100 ml, 500 ml
Cat. No.	1.07961	Entellan™ new rapid mounting medium for microscopy	100 ml, 500 ml, 1 l
Cat. No.	1.08298	Xylene (isomeric mixture) for histology	4 l
Cat. No.	1.09016	Neo-Mount™ anhydrous mounting medium for microscopy	100-ml dropping bottle, 500 ml
Cat. No.	1.09843	Neo-Clear™ (xylene substitute) for microscopy	5 l
Cat. No.	1.15974	Elastica van Gieson staining kit for connective tissue	4x 500 ml

Hazard classification

Cat. No. 1.15973.0002 Please observe the hazard classification printed on the label and the information given in the safety data sheet. The safety data sheet is available on the website and on request.

Main components of the products

Cat. No.	1.15973.0002
Reagent 1	
C.I. 75290	20 g/l
Reagent 2	
Fe(NO ₃) ₃ x 9 H ₂ O	5 g/l
HCl 25%	11.2 g/l

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Literature

1. Romeis - Mikroskopische Technik, Editors: Maria Mulisch, Ulrich Welsch, 2015, Springer Spektrum, 19. Auflage
2. Welsch Sobotta - Lehrbuch Histologie, Editor: Ulrich Welsch, 2006, ELSEVIER Urban&Fischer, 2. Auflage
3. Histotechnik, Gudrun Lang, 2013 Springer Verlag, 2. Auflage
4. Theory and Practice of Histological Techniques, John D Bancroft, Marilyn Gamble, 2008, Churchill Livingstone ELSEVIER, 6th Edition
5. Laboratory Manual of Histochemistry, Linda L. Vacca, 1985, Raven Press
6. Staining Procedures, George Clark, 1981, Williams&Wilkins, 4th Edition
7. Basiswissen Histologie und Zytologie, Karl Heinz Stein, Hellmut Flenker, 1998, uZv, 2. Auflage
8. Histological & Histochemical Methods: Theory & Practice, J. A. Kiernan, 1990, Pergamon Press, 2nd Edition
9. Histological and Histochemical Methods, Theory and practice, J. A. Kiernan, 2015, Scion Publishing Ltd, 5th Edition
10. Conn’s Biological Stains, R.W. Horobin, J.A. Kiernan, 2002, Biological Stain Commission Publication, 10th Edition

H225: Highly flammable liquid and vapor.

H290: May be corrosive to metals.

H319: Causes serious eye irritation.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P234: Keep only in original packaging.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent 1:

H225: Highly flammable liquid and vapor.

H319: Causes serious eye irritation.

Reagent 2:

H290: May be corrosive to metals.

Revision History

Version	Modification Comment
2024-Jul-01	Initial version with the introduction of Revision History

Consult instructions
for use

Manufacturer

Catalog number

Batch code

Caution, consult
accompanying documents

Use by
YYYY-MM-DD

Temperature
limitation

Status: 2024-Jul-01

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.sigmaaldrich.com

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