

# 40839 Atto 647 goat anti-rabbit IgG (whole Molecule)

### **Product Description**

Anti-rabbit IgG (whole molecule) (Sigma, R4880) is developed in goat using purified rabbit IgG as the immunogen. Affinity isolated antigen specific antibody is purified from goat anti-rabbit IgG antiserum to remove essentially all goat serum proteins, including immunoglobulins, which do not specifically bind to rabbit IgG.

Goat anti-rabbit is conjugated to Atto 647 NHS (Abs.max. 645 nm; Em.max. 673 nm), (Cat. No. 18373), then further purified via gel permeation chromatography and dialysis to remove unbound Atto-Dye.

## Reagents

Atto 647 goat anti-rabbit IgG (whole molecule) is provided in unit sizes of 1 ml as 1 mg/ml solution in 0.1 M sodium phosphate, 0.1 M NaCl, pH 7.5, containing 5 mM sodium azide as a preservative.

#### **Working Concentrations**

It is recommended that each individual user determine the optimum working dilution empirically for their systems.

Generally concentrations of 1-10  $\mu$ g /ml are sufficient for many applications.

**Fluorophore / Protein (F/P) Ratio:** 2 to 9 (The exact F/P ratio of each lot is mentioned on the certificate of analysis)

The F/P molar ratio is determined spectrophotometrically as follows:

$$F/P = \frac{A_{644} \times 203,000}{[A_{280} - (A_{644} \times 0.05)] \times 150,000}$$

 $A_{644}=$  absorbance at 644 nm measured in a cuvette with a pathlength of 1 cm  $A_{280}=$  absorbance at 280 nm both measured in a cuvette with a pathlength of 1 cm 203,000 = molar extinction coefficient (\$\epsilon\$) of a IgG at 280 nm [cm\$^-1\$M\$^-1] 150,000 = molar extinction coefficient (\$\epsilon\$) of the Atto 647N dye at 644 nm [cm\$^-1\$M\$^-1] 0.05 = correction factor due to the fluorophore's absorbance at 280 nm

#### Storage / Stability

For continuous use, store at 2-8 °C for up to three months. For extended storage, the solution may be frozen in working aliquots at -20 °C. Frozen aliquots are stable for at least six months. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Protect fluorescent conjugates from light.

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources. © 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

