

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Schneider's Insect Media

Many insect tissue culture media were formulated to mimic the main physico-chemical properties of the body fluid of specific insects. A cursory survey of the formulas of culture media designed for insect tissues reveals great qualitative and quantitative disparities in composition. Different media designed for the same species may exhibit less similarity than two media designed for insects of different orders. Various media have been devised for the *in vitro* culture of *Drosophila* cells and tissues. The most extensively used are Schneider's medium and Echalier and Ohanessian's D-22 medium. *Drosophila* cells have been employed to study a variety of biological processes including genetics, endocrinology, physiology, and cell biology as well as recombinant protein expression.

When supplemented with 5–20% heat-inactivated fetal bovine serum, Schneider's medium has been found to support the rapid growth of both primary and established cultures of cells derived from *Drosophila melanogaster*.^{6.7} The medium has been used for the growth and maintenance of the cell lines originally derived by Schneider⁸ from *Drosophila* embryos as well as the culture of cells from other dipteran species.⁹

	S9895	S0146
	[Powder]	[1×]
COMPONENT	g/L	g/L
Inorganic Salts		
CaCl ₂ (anhydrous)	_	0.6
MgSO₄ (anhydrous)	1.807221	1.807221
KCI	1.6	1.6
NaHCO₃	_	0.4
NaCl	2.1	2.1
Na₂HPO₄	0.7	0.7
Amino acids		
β-Alanine	0.5	0.5
L-Arginine	0.6	0.6
L-Aspartic Acid	0.4	0.4
L-Cystine · 2HCl	0.026732	0.026732
L-Cysteine	0.06	0.06
L-Glutamic Acid	0.8	0.8
L-Glutamine	1.8	1.8
Glycine	0.25	0.25
L-Histidine	0.4	0.4
L-Isoleucine	0.15	0.15
L-Leucine	0.15	0.15
L-Lysine	1.65	1.65
L-Methionine	0.15	0.15
L-Proline	1.7	1.7
L-Serine	0.25	0.25
L-Threonine	0.35	0.35
L-Tryptophan	0.1	0.1
L-Tyrosine · 2Na · 2H₂O	0.07202	0.7202
L-Valine	0.3	0.3

Vitamins and others		
Fumaric acid	0.06	0.06
D(+)-Glucose	2.0	2.0
α-Ketoglutaric acid	0.35	0.35
L-(-)-Malic acid	0.6	0.6
Succinic acid	0.06	0.06
D(+)-Trehalose	2.0	2.0
Yeast Extract	2.0	2.0
Add		
CaCl ₂ (anhydrous)	0.6	_
NaHCO₃	0.4	_

References

- Mitsuhashi, J., Media for insect cell cultures. in: Advances in Cell Culture vol. 2, Maramorosch, K., ed., Academic Press (New York, NY: 1982) pp. 133-196.
- Echalier, G., In vitro culture of Drosophila cells and applications in physiological genetics. in: Invertebrate Tissue Culture, Kurstak, E., and Maramorosch, K., eds., Academic Press (New York, NY: 1976) pp. 131-150.
- Schneider, I., and Blumenthal, A. Drosophila cell and tissue culture. in: Biology and Genetics of Drosophila vol. 2B, Ashburner, M., and Wright, T.R.F., eds., Academic Press (New York, NY: 1978) pp. 266-315.
- Kuroda, Y., *Drosophila* tissue culture: Retrospect and prospect. in: Invertebrate Cell Culture Applications, Maramorosch, K., and Mitsuhashi, J., eds., Academic Press (New York, NY: 1982) pp. 53-104.
- Van der Straten, A. et al., Efficient expression of foreign genes in cultured *Drosphila melanogaster* cells using hygromycin B selection. in: Invertebrate Cell System Applications vol. 1, Mitsuhasi, J., ed., CRC Press (Boca Raton, FL: 1989) pp 183-195.
- 6. Schneider, I., Differentiation of larval *Drosophila* eyeantennal discs *in vitro*. J. Exp. Zool., **156**, 91-104 (1964).
- Schneider, I., Histology of larval eye-antennal disc and cephalic ganglia of *Drosophila* cultured *in vitro*. J. Embryol. Exp. Morphol., **15**, 271-279 (1966).
- 8. Schneider, I., Cell line derived from late embryonic stages of *Drosophila melanogaster*. J. Embryol. Exp. Morphol., **27**, 353-365 (1972).
- Schneider, I., in: Handbook of Practical Tissue Culture Applications, Maramorosch, K., and Hirumi, H. eds., Academic Press (New York, NY: 1979).

JG,ISP,PD,JF,ALF,MAM 04/14-1