

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

SILu™Lite SigmaMAb Adalimumab Monoclonal Antibody Standard

recombinant, expressed in CHO cells

Catalog Number **MSQC16** Storage Temperature –20 °C

Product Description

SILu™Lite SigmaMAb Adalimumab is a recombinant monoclonal antibody with a molecular mass of 150 kDa expressed in CHO cells. SigmaMAb Adalimumab is designed to be used as a standard for optimization of bioanalytical assays of Adalimumab.

Each vial of SigmaMAb Adalimumab contains 500 μg of antibody lyophilized from a solution of phosphate buffered saline. Vial content was determined by measuring A_{280} and using an extinction coefficient ($E^{0.1\%}$) of 1.4.

Sequence Information

SigmaMAb Adalimumab Heavy Chain:

EVQLVESGGGLVQPGRSLRLSCAASGFTFDDYAMHWVRQAP
GKGLEWVSAITWNSGHIDYADSVEGRFTISRDNAKNSLYLQ
MNSLRAEDTAVYYCAKVSYLSTASSLDYWGQGTLVTVSSAS
TKGPSVFPLAPSSKSTSGGTAALGCLVKDYFPEPVTVSWNS
GALTSGVHTFPAVLQSSGLYSLSSVVTVPSSSLGTQTYICN
VNHKPSNTKVDKKVEPKSCDKTHTCPPCPAPELLGGPSVFL
FPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVE
VHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVS
NKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSL
TCLVKGFYPSDIAVEWESNGQPENNYKTTPPVLDSDGSFFL
YSKLTVDKSRWQQGNVFSCSVMHEALHNHYTQKSLSLSPGK

SigmaMAb Adalimumab Light Chain:

DIQMTQSPSSLSASVGDRVTITCRASQGIRNYLAWYQQKPG KAPKLLIYAASTLQSGVPSRFSGSGSGTDFTLTISSLQPED VATYYCQRYNRAPYTFGQGTKVEIKRTVAAPSVFIFPPSDE QLKSGTASVVCLLNNFYPREAKVQWKVDNALQSGNSQESVT EQDSKDSTYSLSSTLTLSKADYEKHKVYACEVTHQGLSSPV TKSFNRGEC

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses.

Preparation Instructions

Reconstitute the contents of the vial by adding 500 μ L of ultrapure water or phosphate buffer and mixing vigorously for a 1mg/mL solution.

Note: Avoid PBS for reconstitution.

If the lyophilized powder does not dissolve completely, make the solution slightly acidic by adding 0.1% formic acid until complete dissolution is achieved. The resulting acidic solution should be neutralized to pH 6-7 by addition of a base or dilution into suitable buffer.


Storage/Stability

Store the lyophilized product at -20 °C.

SILu is a trademark of Sigma-Aldrich Co. LLC


SG, JK, MAM, SM 04/20-1

Appendices Figure 1. Mass spectra



(a) Light Chain, calculated mass: 23,412 Da

(b) Heavy Chain, calculated mass: 50,644 Da



(c) Intact of SigmaMAb Adalimumab, calculated mass: 148,080 Da

Deconvoluted mass spectra of (a) light chain, (b) heavy chain and (c) Intact of SigmaMAb Adalimumab. The reduction was performed in non-denaturing conditions, where the interchain disulfide bonds (which are more susceptible to reduction) will break and produce the light chain and heavy chains, while the intrachain disulfide bonds within each individual domain may remain intact.

Table 1.The calculated molecular mass of intact, light chains, heavy chains of SigmaMAb Adalimumab with the most abundant Glycoforms.

Description	Composition	Modification	Avg. Mass (Da)*	Disulfide Bond
Light chain, Reduced	$C_{1027}H_{1610}N_{282}O_{332}S_6$	NA	23411.91	2 intra-chain
Heavy chain, Reduced	$C_{2191}H_{3392}N_{582}O_{677}S_{15} \\ C_{2247}H_{3484}N_{586}O_{716}S_{15} \\ C_{2253}H_{3494}N_{586}O_{721}S_{15} \\ C_{2259}H_{3504}N_{586}O_{726}S_{15}$	(-Lys) G0F (-Lys) G1F (-Lys) G2F (-Lys)	49198.97 50644.31 50806.45 50968.59	4 intra-chain
Native. Intact, Non-reduced	$\begin{array}{c} C_{6436}H_{9972}N_{1728}O_{2018}S_{42} \\ C_{6548}H_{10156}N_{1736}O_{2096}S_{42} \\ C_{6554}H_{10166}N_{1736}O_{2101}S_{42} \\ C_{6560}H_{10176}N_{1736}O_{2106}S_{42} \\ C_{6566}H_{10186}N_{1736}O_{2111}S_{42} \\ C_{6572}H_{10196}N_{1736}O_{2116}S_{42} \end{array}$	(-2Lys) G0F+G0F (-2Lys) G0F+G1F (-2Lys) G1F+G1F (-2Lys) G1F+G2F (-2Lys) G2F+G2F (-2Lys)	145189.5 148080.2 148242.3 148404.5 148566.6 148728.7	16 (12 intra chain and 4 inter chain)

G0F: GlcNAc2Man3GlcNAc2Fuc G1F: GalGlcNAc2Man3GlcNAc2Fuc G2F: Gal2GlcNAc2Man3GlcNAc2Fuc

^{*}Masses based on NIST Physical Reference Data