

New Phosphorylation Assay Kits

PROTEIN PHOSPHORYLATION ASSAYS

Sigma-RBI is pleased to introduce a series of phosphospecific assay kits for the precise quantitation of phosphorylated or non-phosphorylated cell signaling proteins involved in phosphorylation cascades. Phosphospecific ELISAs offer researchers a sensitive, quantitative and economical alternative to traditional immunoblotting or functional assays.

Assay Platform

The format is a solid phase sandwich Enzyme-Linked Immunosorbent Assay (ELISA). A capture antibody specific for a protein (regardless of phosphorylation state) has been coated onto the wells of the microtiter plate. Antigen binds to the capture antibody. A detection antibody specific for the phosphorylated or non-phosphorylated protein is added and binds to the immobilized protein. An anti-rabbit IgG-HRP completes the four-member sandwich. The reaction is visualized by tetramethylbenzidine (TMB) substrate and the intensity of the color is directly proportional to the concentration of protein present in the original sample. The reaction is read at 450 nm in a microtiter plate reader.

Samples

Cell extracts are the samples of choice, because the proteins involved in the phosphorylation are present in intracellular compartments. Researchers may use the cell extraction procedure provided in the Technical Bulletin or utilize the **Sigma Mammalian Cell Lysis Kit** (Prod. No. **MCL-1**).

Advantages

Sensitive

- Picogram/mL quantities for non-phosphorylated ELISAs
- < 1.0 unit/mL for phosphorylated ELISAs
- 2 to 10 times more sensitive than traditional immunoblotting

Specific

- Highly specific for target protein and phosphorylation site
- Validated against immunoblotting performed with the same antigen and antibodies
- Validated by peptide competition

Quantitative

- Standard curve is run in each assay

Simple

- Precoated plates eliminate the coating step
- All incubations performed at room temperature
- Small sample size – only 10 μ L (diluted to a total volume of 100 μ L)
- Use any microtiter reader available on the market
- No additional equipment required

Stable

- Entire kit is stable in the refrigerator for 12 months
- Components are all stored at 2-8°C
- Stable liquid chromogen and stop reagents
- Two vials of standard extend the life-span of the kit

Fast

- Results in 4 hours

Economical

- Reagents may be used for multiple kit runs

Performance Characteristics

- Low intra- and inter-assay variability
- Linear standard curve
- Recovery ranges > 93%
- Each phosphospecific kit is normalized against non-phosphorylated ELISA

Sample Data using Sigma-RBI Phosphospecific ELISAs

Western Blot:
(38 kDa)



ELISA:
OD 450 nm


0.160	0.212	0.258	0.361	0.554	0.887	1.508	2.978
-------	-------	-------	-------	-------	-------	-------	-------

p38 MAPK
pTpY180/182
(Units/test)


0	0.15	0.31	0.63	1.25	2.5	5	10
---	------	------	------	------	-----	---	----

Detection of p38 MAPK [pThr¹⁸⁰/pTyr¹⁸²] by ELISA vs. Immunoblotting.

The sensitivity of this ELISA was compared to immunoblotting using known quantities of p38 MAPK (pThr¹⁸⁰/pTyr¹⁸²). The data presented show that the sensitivity of the ELISA is approximately 10X greater than that of immunoblotting. The bands shown in the immunoblotting data were developed using rabbit anti-p38 MAPK (pThr¹⁸⁰/pTyr¹⁸²), an alkaline phosphatase conjugated anti-rabbit IgG, followed by chemiluminescent substrate and autoradiography.



Recombinant p38 MAPK was phosphorylated using MKK6 enzyme *in vitro*. Non-phosphorylated p38 MAPK was used as control. The phosphorylated and non-phosphorylated p38 MAPK were analyzed with phospho-p38 MAPK (pThr¹⁸⁰/pTyr¹⁸²) ELISA and p38 MAPK ELISA. This phospho-p38 MAPK (pThr¹⁸⁰/pTyr¹⁸²) ELISA kit is specific for measurement of phosphorylated p38 MAPK at threonine 180 and tyrosine 182. The kit does not detect non-phosphorylated p38 MAPK protein.



Phospho-p38 MAPK (pThr¹⁸⁰/pTyr¹⁸²) ELISA Peptide Blocking. The specificity of this assay for p38 MAPK phosphorylated at threonine 180 and tyrosine 182 was confirmed by peptide competition. The data presented show that only the phosphorylated peptide containing phosphorylated threonine and tyrosine could block the ELISA signal. The same sequence containing non-phosphorylated threonine and tyrosine did not block the signal.

New Phosphorylation Assay Kits

Akt/PKB ELISA
Prod. No. **CS0160**

Phospho-Akt/PKB (pThr³⁰⁸) ELISA
Prod. No. **CS0110**

Phospho-Akt/PKB (pSer⁴⁷³) ELISA
Prod. No. **CS0120**

Epidermal Growth Factor Receptor (EGFR) ELISA
Prod. No. **CS0080**

Phospho-Epidermal Growth Factor Receptor (EGFR)
(pTyr¹⁰⁶⁸) ELISA
Prod. No. **CS0150**

Phospho-Epidermal Growth Factor Receptor (EGFR)
(pTyr¹¹⁷³) ELISA
Prod. No. **CS0140**

ERK 1&2 ELISA
Prod. No. **EK0100**

Phospho-ERK 1&2 (pThr¹⁸⁵/pTyr¹⁸⁷) ELISA
Prod. No. **PE0100**

Insulin Receptor β Subunit ELISA
Prod. No. **CS0090**

Phospho-Insulin Receptor β Subunit (pTyr¹¹⁵⁸) ELISA
Prod. No. **PI0100**

Phospho-Insulin Receptor β Subunit (pTyr¹¹⁶²/pTyr¹¹⁶³) ELISA
Prod. No. **PI0200**


JNK 1&2 ELISA Human
Prod. No. **CS0100**

Phospho-JNK 1&2 (pThr¹⁸³/pTyr¹⁸⁵) ELISA
Prod. No. **CS0130**

p38 MAPK ELISA
Prod. No. **PM0100**

Phospho-p38 MAPK (pThr¹⁸⁰/pTyr¹⁸²) ELISA
Prod. No. **CS0020**

New Cell Signaling & Neuroscience Literature



Assay Kits for Cell Signaling Brochure

- Over 90 assays kits in the areas of apoptosis, cyclic nucleotides, lipid signaling, neuroscience, phosphorylation and more!
- Featured products include: phosphospecific ELISAs, ELISAs, Starbright® Fluorescent technology, enzyme immunoassays and more!




Neuropsychiatric Disorders Brochure

- Featuring overviews on schizophrenia, depression, anxiety, bipolar disorder, attention-deficit hyperactivity disorder
- Over 350 products are featured, including over 50 new products



Apoptosis Brochure

- Overviews and 10 detailed apoptosis pathways outlining Fas signaling, caspase cascades, activation and inhibition of apoptosis and much more!
- 600 kits and reagents for apoptosis research are featured, including over 50 new products



Apoptosis Poster

- Over 300 products
- Numerous pathway illustrations

Fill out & return the Business Reply Card to receive your new Sigma-RBI Literature!