1.14546.0001

Spectroquant® Phosphate Cell Test

P

for the determination of orthophosphate
USEPA equivalent for drinking water and wastewater

1. Method

In sulfuric solution orthophosphate ions react with ammonium vanadate and ammonium heptamolybdate to form orange-yellow molybdovanadophosphoric acid that is determined photometrically ("VM" method).

The method is analogous to APHA 4500-P C.

2. Measuring range and number of determinations

Measuring range	Number of determinations
0.5 - 25.0 mg/l PO ₄ -P	25
1.5 - 76.7 mg/l PO ₄ ³⁻	
1.1 - 57.3 mg/l P ₂ O ₅	

For programming data for selected photometers / spectrophotometers see www.sigmaaldrich.com/photometry.

3. Applications

This test measures only orthophosphate.

Sample material:

Groundwater and surface water, seawater
Wastewater
Industrial water
Boiler water
Nutrient solutions for fertilization
Soils after appropriate sample pretreatment

4. Influence of foreign substances

This was checked individually in solutions containing 10 and 0 mg/l PO₄-P. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %							
Ag ⁺	1000	Cu ²⁺	1000	Pb ²⁺	10	EDTA	0.2%
AsO ₄ ²⁻	50	Fe ³⁺	10	S ²⁻	10	Hydrazine	10
Ca ²⁺	1000	Hg ²⁺	1000	SCN ⁻	1000	Na-acetate	10%
Cd ²⁺	1000	Mg ²⁺	1000	SiO ₃ ²⁻	500	NaCl	20%
CN ⁻	1000	Mn ²⁺	1000	SO ₃ ²⁻	1000	NaNO ₃	20%
Co ²⁺	100	NH ₄ ⁺	1000	Zn ²⁺	1000	Na ₂ SO ₄	20%
Cr ³⁺	50	Ni ²⁺	500				
Cr ₂ O ₇ ²⁻	5	NO ₂ ⁻	1000				

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

25 reaction cells
1 sheet of round stickers for numbering the cells

Other reagents and accessories:

MQuant® Phosphate Test, Cat. No. 1.10428,
measuring range 10 - 500 mg/l PO₄³⁻ (3.3 - 163 mg/l PO₄-P)
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
Phosphate standard solution CRM Certipur®, 1000 mg/l PO₄³⁻,
Cat. No. 1.04690
Hydrochloric acid 25% for analysis EMSURE®, Cat. No. 1.00316
Pipette for a pipetting volume of 5.0 ml

6. Preparation

- Use only phosphate-free detergents to rinse glassware. Otherwise fill with hydrochloric acid (approx. 10%) and leave to stand for several hours.
- Analyze immediately after sampling.
- Check the phosphate content with the MQuant® Phosphate Test. Samples containing more than 25.0 mg/l PO₄-P must be diluted with distilled water.
- The pH must be within the range 0 - 10.** Adjust, if necessary, with sulfuric acid.
- Filter turbid samples.

7. Procedure

Pretreated sample (10 - 40 °C)	5.0 ml	Pipette into a reaction cell, close the cell tightly, and mix.
Measure the sample in the photometer.		

Notes on the measurement:

- For photometric measurement the cells must be clean. Wipe, if necessary, with a clean dry cloth.
- Measurement of turbid solutions yields false-high readings.
- The pH of the measurement solution must be within the range 0.5 - 1.0.
- The color of the measurement solution remains stable for at least 60 min after the end of the reaction time stated above.

8. Analytical quality assurance

recommended before each measurement series

To check the photometric measurement system (test reagents, measurement device, handling) and the mode of working, a diluted phosphate standard solution containing 15.0 mg/l PO₄-P (46.0 mg/l PO₄³⁻) can be used.

Sample-dependent interferences (matrix effects) can be determined by means of standard addition.

Additional notes see under www.sigmaaldrich.com/qa-test-kits. For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Note

Dispose of chemical waste in accordance with the local regulations. Information on disposal can also be found at www.sigmaaldrich.com/spectroquant-retrologistic.

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