Instructions for Use

Leukocyte Acid Phosphatase (TRAP) Kit

Procedure No. 386

Intended Use

Sigma-Aldrich Acid Phosphatase kits are intended as a general laboratory use histological stain. Acid Phosphatase reagents are for professional "In Vitro Diagnostic Use" only. This manual, qualitative, histochemical procedure demonstrates acid phosphatase activity in leukocytes of human blood, bone marrow, and tissue touch preparations. Histochemical visualization of leukocyte associated acid phosphatase is a unique technique currently commonly used in medicine.

Use of substituted naphthol AS phosphates in conjunction with diazonium salts for detection of acid phosphatase in human leukocytes was first reported by Goldberg and Barka.¹ The Sigma-Aldrich procedure for demonstrating this leukocyte enzyme is based upon use of naphthol AS-BI phosphate and fast garnet GBC salt. The latter diazo compound was selected because it couples rapidly at acid pH, forming highly insoluble dye deposits. Tartaric acid is provided in order to demonstrate the presence of tartrate-resistant acid phosphatase.

According to the Sigma-Aldrich procedure, blood films are incubated in a solution containing naphthol AS-BI phosphoric acid and fast garnet GBC. Duplicate films are treated with L(+)-tartrate-containing substrate. Naphthol AS-BI, released by enzymatic hydrolysis, couples immediately with fast garnet GBC forming in soluble maroon dye deposits at sites of activity. Cells containing tartaric acid-sensitive acid phosphatase are devoid of activity, whereas those mononuclear cells containing tartaric acid-resistant phosphatase are not affected by such treatment.

Reagents

Citrate Concentrate (Cat. No. 3861-20ML)

Citrate buffer, 0.38 mol/L, pH 5.4 when diluted according to procedure

Tartrate Solution (Cat. No. 3862-20ML)

L(+)-Tartrate buffer, 0.67 mol/L, pH 5.2 when diluted according to procedure

Acetate Solution (Cat. No. 3863-50ML)

Acetate buffer, 2.5 mol/L, pH 5.2 when diluted according to procedure

Naphthol AS-BI Phosphoric Acid Solution (Cat. No. 3864-40ML)

Naphthol AS-BI phosphoric acid, 12.5 mg/mL, in N,N'-dimethyl formamide. Danger. Harmful in contact with skin. Causes serious eye irritation. Harmful if inhaled. May damage the unborn child. Wear protective gloves/protective clothing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Fast Garnet GBC Salt (Cat. No. 38615-10CAP)

Fast garnet GBC salt, purified, 15 mg/capsule, C.I. 37210. Danger. May cause cancer. Obtain special instructions for use. IF exposed or concerned: Get medical advice/attention.

Acid Hematoxylin Solution (Cat. No. 2852-100ML)

Hematoxylin, certified, 1 g/L, C.I. 75290, sodium iodate, aluminum ammonium sulfate, and stabilizers, pH 3.3 at 25°C.

Special Materials Required but Not Provided

- Acetone, ACS Reagent
- 37°C water bath

Storage and Stability

Store Citrate Concentrate, Tartrate Solution and Acetate Solution in refrigerator (2–8°C). Discard if turbidity develops.

Store Fast Garnet GBC Salt and Naphthol AS-BI Phosphoric Acid Solution in refrigerator (2–8°C). Color of Naphthol AS-BI Phosphoric Acid Solution may vary from pale to dark yellow. Reagents are stable until expiration date shown on bottles.

Store Acid Hematoxylin Solution at room temperature (18–26°C). Discard when time for suitable staining exceeds time recommended in procedure by more than 5 minutes.

Dilute citrate solution is suitable for use in the absence of visible microbial growth. Citrate/Acetone solution should be kept tightly capped when not in use. Prepare fresh daily.

Preparation

Dilute Citrate solution is prepared by adding 2.0 mL Citrate Concentrate to 18.0 mL water.

The Citrate/Acetone solution is prepared by adding 20.0 mL dilute citrate solution to 30.0 mL acetone.

NOTE: Keep Citrate/Acetone Solution tightly capped when not in use. Prepare fresh daily.

Precautions

These IVD's are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Blood, bone marrow films, tissue touch preparations and cytocentrifuge preparations may be used. Either EDTA or heparin will serve as anticoagulants. Blood or bone marrow films may be stored fixed at room temperature (18–26°C) for several weeks or unfixed for several days without appreciable change in activity.² Do not ship whole blood for assay at other laboratories. Send fixed or unfixed slides. Slides should be kept cool during transit. Allow films to dry at least 1 hour prior to fixation.

Notes

Test performance can be monitored by including blood films prepared from a healthy individual. In the absence of tartaric acid, all leukocytes demonstrate granular sites of activity. When incubated with tartrate, an occasional granule may be observed in lymphocytes, but all other normally occurring cells are devoid of activity.

Procedure

Fixation

1. Fix slide in Citrate/Acetone Solution for 30 seconds at room temperature (18–26°C).
2. Rinse slides in deionized water and allow to air dry for at least 15 minutes.

Staining

1. Label two 100 mL beakers A and B and add the following, mixing after each addition:

	Beaker A	Beaker B
Water, prewarmed to 37°C	46.0 mL	44.0 mL
Acetate Solution	2.0 mL	2.0 mL
Naphthol AS-BI Phosphoric Acid	2.0 mL	2.0 mL
Tartrate Solution	—	2.0 mL

2. Add contents of 1 capsule Fast Garnet GBC Salt to each beaker. Stir on magnetic mixer 30–60 seconds.
3. Rapidly filter solutions through Whatman® No. 54 (or equivalent) paper into separate labeled staining jars.
4. Warm solutions in staining jars to 37°C in water bath. Check that temperature is at 37°C before adding slides.
5. Place at least 2 duplicate fixed blood films in each jar containing the filtered substrate solution.
6. Incubate slides 1 hour at 37°C in the dark.
7. Remove slides and wash in deionized water for 3 minutes.
8. Stain in Acid Hematoxylin Solution for 5 minutes.
9. Rinse slides in deionized water for 3 minutes and allow to air dry.
10. Evaluate microscopically using oil immersion (900X) lens. Use of solvent based mounting media is not recommended as they may cause fading of the stains.

Performance Characteristics

Stained blood films are usually evaluated subjectively for the presence or absence of tartrate-resistant enzyme. In the absence of tartaric acid, most leukocytes demonstrate granular sites of activity. When incubated with tartrate, an occasional granule may be observed in lymphocytes and some specialized macrophages.³ In blood smears a positive reaction is denoted by the presence of more than two cells with diffuse and intense activity, i.e., more than 40 granules.³ To evaluate Golgi staining, characteristic of thymus derived lymphocytes (T-cells), Sigma-Aldrich Procedure No. 181 is recommended.

Acid Phosphatase activity appears as purplish to dark red granules in the cytoplasm of most leukocytes.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
3861	Citrate Concentrate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3862	Tartrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3863	Acetate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3864	Naphthol AS-BI Phosphoric Acid Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
38615	Fast Garnet GBC Salt	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
2852	Acid Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

386A:

- H290 May be corrosive to metals.
- H312 Harmful in contact with skin.
- H314 Causes severe skin burns and eye damage.
- H331 Toxic if inhaled.
- H350 May cause cancer.
- H360 May damage the unborn child.
- P202 Do not handle until all safety precautions have been read and understood.
- P234 Keep only in original packaging.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
- P304 + P340 + P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297, 1962
- Janckila A, Li CY et al: The cytochemistry of the tartrate-resistant acid phosphatase-technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 2.0	2014
Rev. 3.0	2016
Rev. 4.0	2023

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Removed reference number 4. Added Warnings and Hazards. Added CH-REP information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Saure Phosphotase, Leukozyten (TRAP) – Kit

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Referenzen

- 1. Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- 2. Jancikla A, Li CY, et al: The cytochemistry of the tartrate-resistant acid phosphatase-technical considerations. Am J Clin Pathol 70:45, 1978
- 3. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Mode d'emploi

Kit de phosphatase acide leucocytaire (TRAP)

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MDSS GmbH