

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# SILu™Lite SigmaMAb Cetuximab Monoclonal Antibody Standard

recombinant, expressed in CHO cells

Catalog Number **MSQC18** Storage Temperature –20 °C

#### **Product Description**

SILu™Lite SigmaMAb Cetuximab is a recombinant monoclonal antibody with a molecular mass of ~150 kDa expressed in CHO cells. SigmaMAb Cetuximab is designed to be used as a standard for optimization of bioanalytical assays of Cetuximab.

Each vial of SigmaMAb Cetuximab contains 500  $\mu$ g of lyophilized antibody from a solution of phosphate buffered saline. Vial content was determined by measuring A<sub>280</sub> and using an extinction coefficient (E<sup>0.1%</sup>) of 1.4.

## Sequence Information

## SigmaMAb Cetuximab Heavy Chain:

QVQLKQSGPGLVQPSQSLSITCTVSGFSLTNYGVHWVRQSP GKGLEWLGVIWSGGNTDYNTPFTSRLSINKDNSKSQVFFKM NSLQSNDTAIYYCARALTYYDYEFAYWGQGTLVTVSAASTK GPSVFPLAPSSKSTSGGTAALGCLVKDYFPEPVTVSWNSGA LTSGVHTFPAVLQSSGLYSLSSVVTVPSSSLGTQTYICNVN HKPSNTKVDKRVEPKSCDKTHTCPPCPAPELLGGPSVFLFP PKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVH NAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNK ALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTC LVKGFYPSDIAVEWESNGQPENNYKTTPPVLDSDGSFFLYS KLTVDKSRWQQGNVFSCSVMHEALHNHYTQKSLSLSPG

## SigmaMAb Cetuximab Light Chain:

DILLTQSPVILSVSPGERVSFSCRASQSIGTNIHWYQQRTN GSPRLLIKYASESISGIPSRFSGSGSGTDFTLSINSVESED IADYYCQQNNNWPTTFGAGTKLELKRTVAAPSVFIFPPSDE QLKSGTASVVCLLNNFYPREAKVQWKVDNALQSGNSQESVT EQDSKDSTYSLSSTLTLSKADYEKHKVYACEVTHQGLSSPV TKSFNRGEC

#### **Precautions and Disclaimer**

This product is for R&D use only. Not for drug, household or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

#### **Preparation Instructions**

Reconstitute the contents of the vial by adding 500  $\mu\text{L}$  of ultrapure water or phosphate buffer, and mixing vigorously for a 1 mg/mL solution.

If the lyophilized powder does not dissolve completely, make the solution slightly acidic by adding 0.1% formic acid until complete dissolution is achieved. The resulting acidic solution should be neutralized to pH 6–7 by addition of a base or dilution into suitable buffer. Note: **Avoid PBS for reconstitution.** 


## Storage/Stability


Store the lyophilized product at -20 °C.

SILu is a trademark of Sigma-Aldrich Co. LLC.

PJ,MAM 10/19-1

Appendices
Figure 1.
Mass Spectra





Deconvoluted mass spectra of partially reduced (a) light chain and (b) heavy chain SigmaMAb Cetuximab. The reduction was performed in non-denaturing conditions, where the interchain disulfide bonds (which are more susceptible to reduction) will break and produce the light chain and heavy chains, while the intrachain disulfide bonds within each individual domain may remain intact.

**Table 1.**The calculated molecular mass of light chains, heavy chains of fully reduced SigmaMAb Cetuximab with the most abundant glycoforms in this product

Description	Composition	Modification*	Average Mass (Da)**	Disulfide bond***
Light chain reduced	$C_{1025}H_{1599}N_{281}O_{338}S_5$	NA	22,426.73	2 intra-chain
Heavy chain reduced	C2202H3393N579O673S15 C2314H3577N587O751S15 C2320H3587N587O756S15 C2326H3597N587O766S15 C2332H3607N587O766S15 C2337H3614N588O769S15 C2343H3624N588O774S15 C2348H3631N589O777S15 C2354H3641N589O782S15	PyroGlu PyroGlu, G0FG0F PyroGlu, G0FG1F PyroGlu, G1FG1F PyroGlu, G1FG2F PyroGlu, G1FG1FS PyroGlu, G1FG2FS PyroGlu, G1FG2FS PyroGlu, G1FG1FS2 PyroGlu, G1FG2FS2	49226.07 52116.75 52278.89 52441.03 52603.17 52732.29 52894.43 53023.54 53185.68	4 intra-chain

G0F: GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>Fuc G1F: GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>GalFuc G2F: GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub> Gal<sub>2</sub>Fuc

G1FS: GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>GalFuc Neu5Ac1 G2FS: GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>Gal<sub>2</sub>Fuc Neu5Ac1

\* C-terminal Lys removed from the sequence and accounted in the table

\*\* Masses based on NIST Physical Reference Data

<sup>\*\*\*</sup> Intra disulfide bonds remain intact after partial reduction