

Data Sheet

HL-1 Cardiac Muscle Cell Line

SCC065

Pack Size: $\geq 1 \times 10^6$ viable cells/vial

Store in liquid nitrogen.

FOR RESEARCH USE ONLY

Not for use in diagnostic procedures. Not for human or animal consumption.

Background

The HL-1 Cardiac Muscle Cell Line is an immortalized mouse cardiomyocyte cell line able to continuously divide and spontaneously contract while maintaining a differentiated cardiac phenotype. HL-1 can be serially passaged without losing its differentiated cardiac myocyte phenotype, including morphological, biochemical, and electrophysiological properties. HL-1 is used in a variety of model systems to address questions of cardiac biology at the cellular and molecular levels.

Spontaneous contraction (for example, beating) of HL-1 may not be observable by bright-field microscopic inspection as the beating is extremely subtle. The beating phenotype may be observed by using a fluorescent calcium flux dye, Fluo-8®. Please refer to the protocol for the Fluo-8® functional assay.

Source

This product contains genetically modified organisms. The HL-1 line was derived from AT-1 subcutaneous tumor excised from an adult female C57BL/6J mouse. The parental AT-1 line was originally derived from an atrial tumor growing in a transgenic mouse in which expression of the SV40 large T-antigen was targeted to atrial cardiomyocytes via the atrial natriuretic factor (ANF) promoter.

Quality Control Testing

- Each vial contains $\geq 1 \times 10^6$ viable cells.
- Cells tested negative for infectious diseases by the Mouse Essential CLEAR panel PCR panel by Charles River Animal Diagnostic Services.
- Cells are negative for mycoplasma contamination.
- Functional Assay: Fluo-8® calcium dye fluorescent based assay for detecting changes in intracellular calcium.

Storage and Handling

The HL-1 Cardiac Muscle Cell Line should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

Representative Data




Figure 1. Cell morphology and immunocytochemistry staining of HL-1 cardiac muscle cell line. Bright-field images of HL-1 cells (**A**) one and (**B**) three days after thawing in a T25 flask. (**C**) HL-1 cells express α -actinin, (**D**) Titin, (**E**) Myosin (MF-20), and (**F**) ANF.

Protocols

Gelatin/Fibronectin ECM Coating of Flasks

1. Add 10 mL of 0.1% Gelatin Solution (ES-006) to 40 mL of Ultrapure Water in a sterile 50 mL conical tube. Final concentration = 0.02% Gelatin Solution
2. Add 250 μ L Fibronectin (1 mg/mL, F-1141) to the 50 mL 0.02% Gelatin Solution. Final concentration = 5 μ g/mL Fibronectin
3. Sterilize filter using 0.2 μ m filter (SE1M179M6 or SCGP00525).
4. Aliquot 6 mL of Gelatin/Fibronectin ECM Solution into sterile 15 mL centrifuge tubes. Label and store at -20°C .
5. Coat flasks with the Gelatin/Fibronectin ECM mixture (3 mL for T25, 6 mL for T75 or 15 mL for T225 flasks) at $2-8^{\circ}\text{C}$ overnight or at 37°C for a minimum of 1 hour. Flasks may be coated 5-6 days in advance and stored at $2-8^{\circ}\text{C}$.

Preparation of Medium

Preparation of 100X (10 mM) Stock Solution of (+)- Norepinephrine [(+)-arterenol] (A0937)

1. Prepare 50 mL of 30 mM L-Ascorbic Acid (A7506) by dissolving 0.264 g L-Ascorbic acid in 50 mL of distilled water.
2. Add 160 mg Norepinephrine to the 50 mL of 30 mM L-Ascorbic acid solution. Sterile filter using 0.2 μ m filter.
3. Aliquot the 100X Norepinephrine stock solution to sterile microtubes and store at -20°C for later use.
4. Wrap 100X Norepinephrine stock solution aliquots in aluminum foil to protect them from light. Minimize freeze-thaw cycles. Aliquots stored at -20°C may be used for up to 1 month.

Preparation of HL-1 Expansion Medium (10% FBS; 500 mL). Scale according to the volumes required.

Component	Quantity	Final Concentration	Cat. No.
Claycomb Basal Medium	435 mL		51800C
HL-1 Qualified FBS	50 mL	10%	TMS-016-B
Norepinephrine (10 mM, 100X)	5 mL	0.1 mM	A0937
L-Glutamine, 200 mM	5 mL	2 mM	TMS-002-C
Penicillin/Streptomycin, 100X	5 mL	1X	TMS-AB2-C

Wrap HL-1 Expansion Medium in aluminum foil as the supplemented medium is extremely light sensitive. Store at 4 °C. HL-1 Expansion Medium should be used within two weeks after preparation.

Thawing the Cells

The HL-1 cells should always be maintained at high cell density. Frozen vials should be thawed into a T25 flask. Culture medium should be exchanged daily. Cells should be passaged at 1:3 split ratio and only when they are 100% confluent. Do not split cells beyond 1:3 as this may cause cells to dedifferentiate.

1. Do not thaw the cells until the Gelatin/Fibronectin coated tissue culture flasks and recommended expansion medium are on hand.
2. Remove the vial of frozen HL-1 cells from liquid nitrogen and thaw in a 37 °C water bath. Closely monitor the cells to ensure complete, but not overthawing. Maximum cell viability depends on the rapid and complete thawing of frozen cells.

Important: Do not vortex the cells.

3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol.
4. In a laminar flow hood, use a 1- or 2-mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
5. Using a 10 mL pipette, slowly add dropwise 9 mL of HL-1 Cell Expansion Medium to the 15 mL conical tube.

Important: Do not add the entire volume of medium all at once to the cells. This may result in decreased cell viability due to osmotic shock.

6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.

Important: Do not vortex the cells.

7. Centrifuge the tube at 300 x *g* for 2-3 minutes to pellet the cells.
8. Decant as much of the supernatant as possible.
9. Resuspend the pelleted cells into 10-15 mL of HL-1 Expansion Medium.
10. Transfer the cell suspension to a Gelatin/Fibronectin coated T25 tissue culture flask.
11. Incubate the cells at 37 °C in a humidified incubator with 5% CO₂.

12. The next day, exchange the medium with 10-15 mL of fresh HL-1 Cell Expansion Medium. Exchange with fresh medium every day.

13. When the cells are 100% confluent, they can be dissociated with 0.05% Trypsin in 0.02% EDTA (T3924) and further passaged or, alternatively, frozen for later use.

Note: It is absolutely critical to use Trypsin-EDTA solution.

14. Cells should be passaged at 1:3 ratio into the appropriate Gelatin/Fibronectin ECM coated flasks. Do not passage cells beyond 1:3 split ratio as this may cause cells to dedifferentiate.

Cryopreservation of Cells

The HL-1 Cardiac Cell Line may be frozen in 95% HL-1 Qualified FBS (TMS-016-B) containing 5% DMSO using a Nalgene® slow freeze Mr. Frosty® container.

Protocol for Detection of Spontaneous Calcium Fluctuations in Confluent HL-1 Cells

Please have on hand a 100% confluent T25 flask of HL-1 cells. Please refer to protocol above on culture of HL-1 cells.

1. Prepare a Gelatin/Fibronectin coated 24-well tissue culture plate with 0.5 mL per well of ECM. Please refer to Gelatin/Fibronectin ECM coating of flasks in the protocol section above.
2. Prepare fresh HL-1 Expansion Medium. Please refer to protocol above for the formulation.
3. Dissociate the confluent layer of HL-1 cells with 0.5% Trypsin in 0.02% EDTA.
Note: It is absolutely critical to use Trypsin EDTA (T3924).
4. Count and dilute 3 million cells into 24 mL of fresh HL-1 Expansion Medium. Dispense 1 mL per well of the cell suspension into each well of the coated 24-well plate. This corresponds to approximately 125,000 cells per well. Swirl the plate to ensure even plating of the cells.
5. Allow HL-1 cells to grow to 100% confluency with media changes every day. It will typically take 3-4 days for the cells to be 100% confluent. Alternatively, the cells can be plated at a higher cell density but must be monitored to ensure that cells reach confluency but are not overgrown. It is important that fresh media is exchanged daily.

On the day of the assay (for example, when cells reach 100% confluency):

6. Prepare a 100X Fluo-8® AM stock solution (500 μ M): Add 100 μ L of DMSO to 1 vial (50 μ g) of Fluo-8® AM (AAT Bioquest, 21081). Final concentration = 500 μ M. Store at -20 °C in a tinfoil wrapped vial until ready to use.
7. Aspirate the growth medium from the 24-well plate containing confluent HL-1 cells. Rinse each well to be stained with 1 mL of Hepes Buffered Tyrode's Solution (ThermoFisher, 50151910) OR Hanks' Buffer with 20 mM Hepes (HHBS) (AAT Bioquest, 20011).
8. Prepare a 5 μ M Fluo-8® working solution by diluting 15 μ L of the 500 μ M Fluo-8® AM stock solution into 1.5 mL of Tyrode's Solution.
9. Aspirate the wells to be stained and replace with 0.5 mL per well of 5 μ M Fluo-8® solution.
10. Incubate at 37 °C for 1 hour.
11. After 1 hour, aspirate the staining solution. Rinse each well with 1 mL per well of Tyrode's Solution.
12. Aspirate and add 1 mL per well Tyrode's Solution to each well.
13. Immediately view the cells using a fluorescent microscope equipped with a 20X objective and a GFP or FITC filter set. With the fluorescent microscope, spontaneous contraction of HL-1 cells can be readily observed and videotaped.

References

1. Proc. Natl. Acad. Sci. USA. 95(6): 2979-2984.
2. Am. J. Physiol. Heart Circ. Physiol 286(3): H823-9.

Academic Use Agreement

Subject to local law

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means HL-1 Cardiac Muscle Cell Line (SCC065).

"Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

Genetically Modified Organisms (GMO)

This product contains genetically modified organisms.
Este producto contiene organismos genéticamente modificados.
Questo prodotto contiene degli organismi geneticamente modificati.
Dieses Produkt enthält genetisch modifizierte Organismen.
Ce produit contient des organismes génétiquement modifiés.
Dit product bevat genetisch gewijzigde organismen.
Tämä tuote sisältää geneettisesti muutettuja organismeja.
Denna produkt innehåller genetiskt ändrade organismer.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2017-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

SCC065CA Ver 5.0, Rev 14MAY2024, DP

**MILLIPORE
SIGMA**