1.01728.1000
1.01728.9010 **REF**

Microscopy

OSTEOSOFT®

mild decalcifier-solution for histology

For professional use only

IVD In Vitro Diagnostic Medical Device

Intended purpose

This ready-to-use product, "OSTEOSOFT® - mild decalcifier-solution for histology" has been specifically developed for human cell diagnostics in humans and is used for the demineralization and decalcification of bones, punched bone-marrow specimens, and other hard tissue of human origin. Using the auxiliary reagents from our portfolio creates the conditions that enable authorized and qualified investigators to make a correct diagnosis at the end of the diagnostic process. In this regard, auxiliary IVD reagents serve *inter alia* to process human specimen material (e.g. fixing, decalcifying, dehydrating, clarifying, paraffin-embedding, mounting, microscoping, archiving). When used together with the corresponding staining solutions, this enables the visualization of cellular structures that are otherwise low in contrast, thus rendering them evaluable under the optical microscope. Further examinations may be necessary to reach a definitive diagnosis.

Principle

Decalcification methods are necessary for optical microscopic examinations of hard tissue in routine, their purpose being to enable the materials to be cut into thin sections. Usually this decalcification process is carried out using acidic solutions based on HCl. While this method is swift, it results in the tissue undergoing substantial morphological changes as well as the destruction of antigens, meaning for example that it is no longer possible to perform an immunohistological analysis of the tissue decalcified in this manner.

A far more gentle method for decalcifying tissue specimens is by eliminating calcium ions from the tissue by the chelating principle, e.g. using EDTA. This is the principle on which decalcification with OSTEOSOFT® is based.

The decalcifying solution is poured onto the material to be decalcified in excess, thus demineralizing (decalcifying) the tissue in question. The decalcification time is dependent on the size and structural density of the respective tissue.

The composition of the decalcifying solution also has a major impact on the process: OSTEOSOFT® contains chelating agents, which bind calcium ions in the tissue. This preserves the antigen structures in the tissue, thus enabling immunohistological, molecular-biological, and enzymatic (e.g. NASDCL) methods to be carried out. A drawback of this gentle and antigen-preserving decalcification method, however, is that the process lasts much longer compared to the decalcification method based on the acid principle.

The decalcifying solution OSTEOSOFT® is stained yellow to make it easier to identify. The dye used is inert in terms of any effect on the tissue to be decalcified.

Sample material

Sensitive, hard tissue (e.g. punched iliac crest specimens) and calcified tissue (e.g. blood vessels) for the gentle preparation of paraffin sections in histology

Reagents

Cat. No. 1.01728
OSTEOSOFT®
mild decalcifier-solution for histology

1 l, 10 l Titripac®

Alternatively:

Cat. No. 1.01736 OSTEOMOLL®
Solution for rapid decalcification (with fixation)
for histology

1 l, 2.5 l

Sample preparation

The sampling must be performed by qualified personnel.

Sample material to be decalcified with OSTEOSOFT® must be fixed previously (e.g. with Formaldehyde solution 4%, Cat. No. 1.00496).

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

Reagent preparation

The OSTEOSOFT® - mild decalcifier-solution for histology used is ready-to-use, dilution of the solution is not necessary.

Procedure

Place all the tissue to be decalcified into a vessel (glass or plastic) containing an excess of the ready-to-use OSTEOSOFT® decalcifier-solution.

The decalcification time and the quantity of OSTEOSOFT® required are dependent on the size, type, and density of the respective material:

Decalcification of porous bone fragments and other hard tissue

Immerse a piece of bone measuring e.g. 15 x 9 x 3 mm in approx. 50 ml of OSTEOSOFT®; here the decalcification process lasts approx. 18 - 24 hours.

Determining the end-point of the decalcification process

The end-point of the decalcification process (the softness of the tissue) is determined by pricking a needle into the material.

Accelerated decalcification of punched iliac crest specimens using the ultrasonic bath

"OSTEOSOFT® - mild decalcifier-solution for histology" can be used to decalcify punched iliac crest specimens in combination with a suitable ultrasonic bath. Compared to the conventional procedure, the use of this method can considerably shorten the time required to decalcify punched iliac crest specimens, at the same time preserving the tissue and the stainability of the cells and tissue structures.

The ultrasonic bath used must be suited for the purpose and must have the following features:

- Ultrasonic bath with automatic cooling for *in vitro* diagnostic sample preparation (e.g. bone decalcification)
- Setting of the running time for at least 7:00 hours
- Temperature controller
- Unsupervised continuous operation
- Selectable ultrasonic power
- Ultrasonic bath with a lid

Program for the decalcification of punched iliac crest specimens with "OSTEOSOFT® - mild decalcifier-solution for histology"

Instrument	Ultrasonic unit with cooling
Ultrasonic power	45 watts
Temperature	25 °C
Running time	7:00 h

The compartments of the instrument must be filled with OSTEOSOFT® decalcifying solution according to the instructions provided by its supplier. The operating manual for the ultrasonic bath must be strictly followed.

The decalcified material is then taken for histoprocessing by the usual methods.

Result

Decalcified material is cartilaginous or rubber-like in its consistency and exhibits only a weak resistance.

Notes on use

Immunohistological procedures can be conducted without difficulty after decalcification with OSTEOSOFT®, since the antigen structures of the tissue are preserved by this decalcification method.

Enzymatic procedures such as the detection of specific esterase with naphthol AS-D chloroacetate and silver impregnation of reticulin fibers can be carried out without difficulty.

Molecular-biological methods can be carried out, since the corresponding macromolecules are not destroyed by this decalcification process.

Decalcification of compact, calcified tissue

In cases where it is not intended to use the specimen material in immunohistological, molecular-biological, or enzymatic tests, it should be decalcified with OSTEOMOLL® - Solution for rapid decalcification (with fixation) for histology, Cat. No. 1.01736, on an acid basis. OSTEOMOLL® should, however, not be used for sensitive material such as punched iliac crest specimens.

Depending on the size of the specimen, compact, calcified material is decalcified after approx. 6 - 72 hours in OSTEOMOLL®.

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using cryostats or histoprocessors, please follow the instructions for use supplied by the supplier of the system and software.

Analytical performance characteristics

The present auxiliary reagent "OSTEOSOFT®" aids in the microscopic examination of biological structures as described in the "Intended purpose" of this IFU. The use of the product is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, histoprocessing, decisions regarding suitable controls and more.

The analytical performance of the product is confirmed by testing each production batch.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100%:

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Decalcification method				
Performance test	20/20	20/20	8/8	8/8

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

The results of this Performance Evaluation confirms that the product is suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used.

This product is an auxiliary reagent that, when used together with other IVD products such as staining solutions, renders human specimen material evaluable for diagnostic purposes.

Further tests must be selected and implemented according to recognized methods.

Suitable controls should be conducted with each application in order to avoid an incorrect result.

Storage

Store the OSTEOSOFT® - mild decalcifier-solution for histology at +15 °C to +25 °C.

Shelf-life

The OSTEOSOFT® - mild decalcifier-solution for histology can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

The bottles must be kept tightly closed at all times.

Any diminution in the color intensity of OSTEOSOFT® has no effect on the decalcification capacity.

Capacity

The quantity of OSTEOSOFT® required is dependent on the size, type, and density of the respective material.

OSTEOSOFT® is used only once; a fresh solution must be prepared for each new tissue specimen.

Additional instructions

For professional use only.

In order to avoid errors, the application must be carried out by qualified personnel only.

National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

If necessary use a standard centrifuge suitable for medical diagnostic laboratory.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.

Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No. 1.00496	Formaldehyde solution 4%, buffered, pH 6.9 (approx. 10% Formalin solution) for histology	350 ml and 700 ml (in bottle with wide neck), 5 l, 10 l, 10 l Titripac®
Cat. No. 1.01736	OSTEOMOLL® Solution for rapid decalcification (with fixation) for histology	1 l, 2.5 l
Cat. No. 1.11609	Histosec™ pastilles solidification point 56-58°C embedding agent for histology	1 kg, 10 kg (4x 2.5 kg), 25 kg
Cat. No. 1.15161	Histosec™ pastilles (without DMSO) solidification point 56-58°C embedding agent for histology	10 kg (4x 2.5 kg), 25 kg

Hazard classification

Cat No. 1.01728

Please observe the hazard classification printed on the label and the information given in the safety data sheet.

The safety data sheet is available on the website and on request.

Main product components

Cat. No. 1.01728

$C_{10}H_{16}N_2O_8$ 177 g/l
pH 7.0 - 7.3

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and / or its authorised representative and to your national authority.

Literature

1. Romeis - Mikroskopische Technik, Editors: Maria Mulisch, Ulrich Welsch, 2015, Springer Spektrum, 19. Auflage
2. Theory and Practice of Histological Techniques, John D Bancroft, Marilyn Gamble, 2008, Churchill Livingstone ELSEVIER, sixth Edition
3. Histological and Histochemical Methods, Theory and practice, J.A. Kiernan, 2015, Scion Publishing Ltd, 5th Edition
4. Laboratory Manual of Histochemistry, Linda L. Vacca, 1985, Raven Press

H318: Causes serious eye damage.

H373: May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure if inhaled.

P260: Do not breathe dust/ fume/ gas/ mist/ vapors/ spray.

P280: Wear eye protection/ face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P314: Get medical advice/ attention if you feel unwell.

P501: Dispose of contents/ container to an approved waste disposal plant.

Revision History

Version	Modification Comment
2024-Jul-01	Initial version with the introduction of Revision History

Consult instructions
for use

Manufacturer

Catalog number

Batch code

Caution, consult
accompanying documents

Use by
YYYY-MM-DD

Temperature
limitation

Status: 2024-Jul-01

MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

EMD Millipore Corporation, 400 Summit Drive, Burlington MA 01803,
USA, Tel. +1-978-715-4321
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville, Ontario,
L6H 6J8, Canada, Phone: +1 800-565-1400
www.sigmaldrich.com

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EMD Millipore Corporation, 400 Summit Drive, Burlington MA 01803, USA, Tel. +1-978-715-4321
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville, Ontario, L6H 6J8, Canada, Phone: +1 800-565-1400