1.00327.1000
1.00327.5000 **REF**

Microscopy

Hydrochloric acid in ethanol

for microscopy

For professional use only

IVD In Vitro Diagnostic Medical Device

Intended purpose

This "Hydrochloric acid in ethanol - for microscopy" is used for human-medical cell diagnosis and serves the purpose of the bacteriological and histological investigation of sample material of human origin. It is a ready-to-use solution for decolorizing that when used together with other *in vitro* diagnostic products from our portfolio makes bacterial target structures evaluable for diagnostic purposes (acid-fast bacteria (AFB)) by fixing, embedding, staining, counterstaining, mounting in bacteriological and histological specimen materials, for example smears of enriched bacterial cultures or histological sections of e.g. the lung.

Unstained structures are relatively low in contrast and are extremely difficult to distinguish under the light microscope. The images created using the staining solutions help the authorized and qualified investigator to better define the form and structure in such cases. Further tests must be carried out according to recognized, valid methods to reach a definitive diagnosis.

Principle

Acid-fast bacteria are difficult to stain because of the high proportion of lipid and wax in their cell walls. The most efficient staining method is the hot staining according to Ziehl-Neelsen. In this method, carbolfuchsin solution is applied to the specimen and heated. This heating process accelerates the rate at which the fuchsin dye is absorbed and thus also that of the formation of the mycolatefuchsin complex in the cell wall. At the Kinyoun method, a modification of the AFB-Color carbol fuchsin solution renders heating unnecessary (cold staining).

Once the acid-fast bacteria have absorbed the fuchsin dye, it is virtually impossible to decolorize them again, even when they are intensively treated with a decolorizing solution such as e.g. hydrochloric acid in ethanol. Accordingly, acid-fast bacteria are termed as acid- and alcohol-fast for staining, and are stained red in the microscopic visualization. Correspondingly, all non-acid-fast microorganisms are counterstained with an appropriate dye. In this instruction for use malachite green is used for counterstaining. Pretreatment of the specimens with Sputofluor™ dissolves the bacteria from the surrounding viscid sputum and cell material.

Sample material

Smears of bacteriological material that have been air-dried, heat-fixed, and pretreated with Sputofluor™ like sputum, smears from fine needle aspiration biopsies (FNAB), rinses, imprints, effusions, pus, exudates, liquid and solid cultures

Sections of formalin-fixed tissue of human origin, embedded in paraffin (3 - 5 µm thick paraffin sections)

Reagents

Cat. No. 1.00327
Hydrochloric acid in ethanol 1 l, 5 l
for microscopy

Also required:

Cat. No. 1.08000 Sputofluor™ 1 l
for microscopy
Cat. No. 1.08512 AFB-Color carbol fuchsin solution 500 ml, 2.5 l
for the microscopic investigation of
acid-fast bacteria (AFB) (cold staining)
Cat. No. 1.10630 AFB-Color malachite green (oxalate) solution 500 ml
for the microscopic investigation of acid-fast
bacteria (AFB) (cold staining)

Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

When using the corresponding auxiliary reagents, the corresponding instructions for use must be observed.

Sputum

The acid-fast bacteria should be pretreated with Sputofluor™ to dissolve them from mucus and cellular structures. In this process, the active ingredient hypochlorite dissolves the organic material by oxidation and gently releases the acid-fast bacteria so that they can be processed further.

Reagent preparation: Preparation of Sputofluor™ solution 15%

For preparation of approx. 100 ml solution mix:

Sputofluor™	15 ml
Distilled water	85 ml

Preparing sample material in centrifuge tubes:	
Sample	1 part (min. 2 ml)
Sputofluor™ solution (15% in distilled water)	3 parts
Shake vigorously	10 min
Centrifuge at 3000 - 4800 rpm	20 min
Decant supernatant Prepare smears of the sediment Air-dry	

Punction and lavage material, sediments

After carrying out the appropriate enrichment measures smear out samples on the microscopic slides and allow to air dry.

Fixation of smear samples

Fixation is carried out over the flame of a Bunsen burner (2 - 3 times, avoiding excessive heating).

It is also possible to fix the smears in an oven at 100 - 110 °C for 20 min.

Impairment of staining must be expected if a higher temperature or longer heating is employed.

Reagent preparation

The Hydrochloric acid in ethanol - for microscopy used for decolorizing is ready-to-use, dilution of the solution is not necessary and merely produces a deterioration of the staining result and its stability.

The staining solutions, solid dyes, and test kits may have to be individually prepared. The preparation procedures are described in the respective instructions for use.

Procedure

The stated times should be adhered to in order to guarantee an optimal staining result.

Example of the cold staining of smear specimens

Staining on the staining rack

Slide with fixed smear		
AFB-Color carbol fuchsin solution	cover completely and stain	20 min
Running tap water	rinse until no further clouds of dye are produced	

Hydrochloric acid in ethanol	cover completely and leave to react	15 - 30 sec*
Running tap water	rinse immediately	
AFB-Color malachite green (oxalate) solution	cover completely and leave to react	1 min
Running tap water	rinse carefully	
Air-dry (e.g. over night or at 50 °C in the drying cabinet)		

* depending on thickness of specimen

Covering with non-aqueous mounting media (e.g. Neo-Mount™, Entellan™ new, or DPX new) and a cover glass is recommended for the storage of bacteriological specimens for several months. For this purpose, the stained specimens must be dried very well.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Alternatively, specimens can be hot-stained by the Ziehl-Neelsen method (see instructions for use for Ziehl-Neelsen carbol fuchsin solution, Cat. No. 1.09215).

Result

Acid-fast bacteria	red
Background	light green

Evaluation

A positive result means “acid-fast bacteria detected” and a negative result “acid-fast bacteria not detected”. A positive result does not mean that a taxonomic classification by microscopy is possible. If acid-fast bacteria are detected, further analyses must be performed in specially equipped laboratories. The vitality (active, inactive) of the bacteria can also not be determined.

Trouble-shooting

Fixation of smear samples

A sufficient degree of heat-fixing using a Bunsen burner or in a heating cabinet is essential to prevent the infectious potential of the specimens and further proliferation of the bacteria.

No staining of acid-fast bacteria

The critical step of this staining process is the decolorizing step, which can be influenced by the thickness of the specimen smear. In addition, a fresh solution of hydrochloric acid in ethanol is highly reactive, meaning that the result should be evaluated with caution. The incubation times stated in this protocol should be kept accurately in the decolorizing step, since otherwise false-negative results may ensue.

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory. When using automatic staining systems, please follow the instructions for use supplied by the supplier of the system and software. Remove surplus immersion oil before filing.

Analytical performance characteristics

“Hydrochloric acid in ethanol” decolorizes and thereby visualizes biological structures, as described in the “Result” chapter of this IFU. The use of the product is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, decisions regarding suitable controls and more.

The analytical performance of the product is confirmed by testing each production batch.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100%:

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
AFB staining				
Acid-fast bacteria	20/20	20/20	7/7	7/7
Background	20/20	20/20	7/7	7/7

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

The results of this Performance Evaluation confirms that the product is suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used. This method can be supplementarily used in human diagnostics. Further tests must be selected and implemented according to recognized methods.

Suitable controls (e.g. ISOSLIDE™ AFB, Cat. No. 1.02560.0001) should be conducted with each application in order to avoid an incorrect result.

Storage

Store the Hydrochloric acid in ethanol - for microscopy at +15 °C to +25 °C.

Shelf-life

The Hydrochloric acid in ethanol - for microscopy can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

The bottles must be kept tightly closed at all times.

Capacity

up to 250 stainings / 500 ml

Additional instructions

For professional use only. In order to avoid errors, the application must be carried out by qualified personnel only. National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used. If necessary use a standard centrifuge suitable for medical diagnostic laboratory.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines. Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link “Hints for Disposal of Microscopy Products” at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No.	1.00579	DPX new non-aqueous mounting medium for microscopy	500 ml
Cat. No.	1.00869	Entellan™ new for cover slipper for microscopy	500 ml
Cat. No.	1.02560	ISOSLIDE™ AFB Control Slides with reference tissue for the detection of acid-fast bacteria in histological tissue	25 tests
Cat. No.	1.04699	Immersion oil for microscopy	100-ml dropping bottle, 100 ml, 500 ml
Cat. No.	1.07961	Entellan™ new rapid mounting medium for microscopy	100 ml, 500 ml, 1 l
Cat. No.	1.08000	SputofluoI™ for microscopy	1 l
Cat. No.	1.08298	Xylene (isomeric mixture) for histology	4 l
Cat. No.	1.08512	AFB-Color carbol fuchsin solution for the microscopic investigation of acid-fast bacteria (AFB) (cold staining)	500 ml, 2.5 l
Cat. No.	1.08562	Aquatex® (aqueous mounting agent) for microscopy	50-ml dropping bottle
Cat. No.	1.09016	Neo-Mount™ anhydrous mounting medium for microscopy	100-ml dropping bottle, 500 ml
Cat. No.	1.09215	Ziehl-Neelsen carbolfuchsin solution for microscopy	100 ml, 500 ml, 2,5 l
Cat. No.	1.09843	Neo-Clear™ (xylene substitute) for microscopy	5 l
Cat. No.	1.10630	AFB-Color malachite green (oxalate) solution for the microscopic investigation of acid-fast bacteria (AFB) (cold staining)	500 ml

Hazard classification

Cat. No. 1.00327
Please observe the hazard classification printed on the label and the information given in the safety data sheet. The safety data sheet is available on the website and on request.

Main components of the product

Cat. No.	1.00327	
HCl		0.75%
contains C ₂ H ₅ OH		

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Literature

- 1. Romeis - Mikroskopische Technik, Editors: Mulisch, Maria, Welsch, Ulrich, 2015, Springer-Verlag Berlin Heidelberg, 19. Auflage
- 2. Theory and Practice of Histological Techniques, John D Bancroft and Marilyn Gamble, 6th Edition
- 3. Conn's Biological Stains: A Handbook of Dyes, Stains and Fluorochromes for Use in Biology and Medicine, 10th Edition, (ed. Horobin, R.W. and Kiernan, J.A). Bios, 2002
- 4. Kurzlehrbuch Medizinische Mikrobiologie und Infektiologie, Editor: Uwe Groß, Thieme 2009, 2. Auflage
- 5. Histological and Histochemical Methods, Theory and practise, J.A. Kiernan, Scion, 5th Edition

H225: Highly flammable liquid and vapor.

H290: May be corrosive to metals.

H319: Causes serious eye irritation.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P234: Keep only in original packaging.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Revision History

Version	Modificaton Comment
2024-Aug-12	Initial version with the introduction of Revision History

Consult instructions
for use

Manufacturer

Catalog number

Batch code

Caution, consult
accompanying documents

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© 2024 Merck KGaA, Darmstadt, Deutschland und/oder Tochterunternehmen. Alle Rechte vorbehalten. Merck und Sigma-Aldrich sind Marken der Merck KGaA, Darmstadt, Deutschland. Alle anderen Marken sind Eigentum der jeweiligen Inhaber. Ausführliche Informationen zu Markennamen sind über öffentlich zugängliche Informationsquellen erhältlich.