Instructions for Use

Leukocyte Peroxidase
(Myeloperoxidase) Kit

Procedure No. 390

Intended Use

Sigma-Aldrich Leukocyte Peroxidase (Myeloperoxidase) kit is intended for use as a general laboratory histochemical stain for the demonstration of leukocyte peroxidase. The Leukocyte Peroxidase (Myeloperoxidase) kit is for professional "In Vitro Diagnostic Use" only. This manual, qualitative procedure identifies leukocyte peroxidase in blood or bone marrow films of human specimens. Histological visualization of leukocyte peroxidase is a unique technique currently commonly used in medicine.

Classic methods of cytochemical localization of myeloperoxidase (MP) have involved use of benzidine¹ or diaminobenzidine.² In 1977, Hanker et al.³ described the use of p-phenylenediamine and catechol to detect injected horseradish peroxidase. That indicator system is the basis for the Sigma-Aldrich procedure when myeloperoxidase is detected by means of the following reaction:

Reagents

TRIZMAL™ 6.3 BUFFER Buffer Concentrate (Cat. No. 903C-50ML)

TRIZMA® maleate, 200 mmol/L. Chloroform added as preservative.

Peroxidase Indicator Reagent (Cat. No. 3901-10VL)

p-Phenylenediamine diHCl (1 part) and catechol (2 parts). Danger. Toxic if swallowed or in contact with skin. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. Harmful if inhaled. Very toxic to aquatic life with long lasting effects.

Acid Hematoxylin Solution (Cat. No. 2852-100ML)

Hematoxylin, certified, 1 g/L, C.I. 75290, pH 3.3 at 25°C.

Special Materials Required but Not Provided

- Formaldehyde, 37%, Solution
- Ethanol, 95% (v/v)
- Phosphate Buffered Saline, pH 7.4 (Cat. No. P3813)
- Hydrogen Peroxide, 30%

Storage and Stability

TRIZMAL™ 6.3 Buffer Concentrate and Acid Hematoxylin Solution should be stored at room temperature (18–26°C).

Peroxidase Indicator Reagent should be stored refrigerated (2–8°C).

Acid Hematoxylin Solution should not be returned to original container after use in Coplin jar.

Hydrogen Peroxide, 3% in Phosphate Buffered Saline Solution should be stored in refrigerator (2–8°C). Discard if turbidity develops.

Reagents are stable until expiration date.

Deterioration

Discard TRIZMAL™ 6.3 Buffer Concentrate if turbidity develops.

Discard Acid Hematoxylin Solution when the time required for suitable staining exceeds the time recommended in the procedure by more than 5 minutes.

Preparation

TRIZMAL™ 6.3 Dilute Buffer is prepared by mixing 1 volume of TRIZMAL™ 6.3 Buffer Concentrate with 9 volumes of deionized water. Use once and discard.

Fixative Solution is prepared by mixing 5 mL of 37% formaldehyde with 45 mL of 95% ethanol. Prepare fresh daily. Store tightly capped.

Hydrogen Peroxide, 3%, in Phosphate Buffered Saline Solution, is prepared by adding 1 part Hydrogen Peroxide, 30%, to 9 parts Phosphate Buffered Saline Solution pH 7.4. Should be prepared fresh.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Freshly prepared whole blood or bone marrow films should be used for the assay. Blood may be collected in heparin or EDTA. Exposure to light should be minimized as leukocyte peroxidase is photolabile. Unfixed films are reported to be stable for 3 weeks when kept in the dark.¹ Allow films to air dry for 10 minutes, protected from light, prior to fixation.

Notes

It is recommended that blood films prepared from healthy donors be processed along with patient samples as a system performance check.

Although myeloperoxidase is generally considered a marker for cells of myelocytic lineage, it is imperative to recognize that monocytoid cells may also display weak peroxidase activity.

Procedure

- Fix films at room temperature for 30 seconds in Fixative Solution.
- Wash slides in gently running tap water for 2 minutes and allow to air dry in the dark for 10 minutes.
- Prewarm 50 mL TRIZMAL™ 6.3 Dilute Buffer in 37°C water bath.
- Immediately before use, add 1 vial Peroxidase Indicator Reagent and 200 µl (0.2 mL) of 3% Hydrogen Peroxide to the prewarmed TRIZMAL™ 6.3 Dilute Buffer. Mix thoroughly. Discard after use.
- Place washed, fixed slides (Step 2) in Peroxidase Indicator Reagent Solution (Step 4) for 30 minutes in the dark in 37°C water bath.
- Following incubation, wash slides in gently running tap water for 15–30 seconds and allow to air dry.
- Counterstain slides in Acid Hematoxylin Solution for 10 minutes.
- Rinse slides in running deionized water for 15–30 seconds. Air dry and examine slides microscopically.

Performance Characteristics

Blood films prepared from normal donors were stained for myeloperoxidase according to this procedure and by a benzidine method.¹ Neutrophils showed brown-black granulation with this procedure and blue granulation with the benzidine procedure. In both cases, monocytes stained less intensely and lymphocytes did not show myeloperoxidase activity.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
903C	TRIZMAL™ 6.3 Buffer Concentrate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3901	Peroxidase Indicator Reagent	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
2852	Acid Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

390A:

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

H410 Very toxic to aquatic life with long lasting effects.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/ doctor if you feel unwell.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, pp 149–424.

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 3.0	2014
Rev. 4.0	2016
Rev. 5.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 EN ISO 15223 for symbols. Added adverse event contact information. Added Warnings and Hazards. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The initial M, TRIZMA, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- 1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- 2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- 3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- 4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- 5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, pp 149–424.

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Références

1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, pp 149-424.

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Bibliografia

- 1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- 2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- 3. Harker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- 4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- 5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, pp 149-424.

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Referencias

1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Sangre 26:215, 1965.
2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, pp 149–424.

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Referenser

- 1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- 2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- 3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- 4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- 5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, sid 149–424.

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Referências

- 1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- 2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- 3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- 4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- 5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Editado por D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febriger, Philadelphia, 1981, pp 149-424.

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- Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febriger, Philadelphia, 1981, pp 149–424.

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- Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febrieger, Philadelphia, 1981, pp 149–424.

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Referanslar

- Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febiger, Philadelphia, 1981, pp 149-424.

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Referenties

- 1. Kaplow SL: Simplified myeloperoxidase stain using benzidine dihydrochloride. Blood 26:215, 1965.
- 2. Graham RC, Karnovsky MJ: The early stage of absorption of injected horseradish peroxidase in the proximal tubule of mouse kidney; Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- 3. Hanker JS, Yates PE, Metz CB, Rustioni A: A new specific sensitive and non-carcinogenic reagent for the demonstration of horseradish peroxidase. Histochem 9:789, 1977.
- 4. Bennett JM, Catovsky D, Daniel M, Randrin G, Galton DAG, Gralnick HR, Sultan C: Proposals for the classification of the acute leukemias. French-American-British (FAB) Co-operative Group. Brit J Haematol 33:451, 1976.
- 5. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. In Atlas of Blood Cells-Function and Pathology, Vol. 1. Edited by D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Lea and Febriger, Philadelphia, 1981, pp 149-424.

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MDSS GmbH