



# **Cell Comb<sup>TM</sup> Scratch Assay**

Catalog No. 17-10191

(Patent Pending)  
**FOR RESEARCH USE ONLY**  
Not for use in diagnostic procedures.

This Page Left Intentionally Blank

---


## Introduction

Cell migration can be studied by a variety of methods, but the scratch assay remains a highly popular method due to the simplicity of the required materials, experimental setup, data collection and interpretation<sup>1,2</sup>. The scratch assay is typically initiated by scratching a confluent cell monolayer with a pipette tip to create a narrow wound-like gap. Shortly after wounding, the cells at the edge of the wound initiate a program to migrate into the gap, a process that continues until the gap has been completely repopulated with cells. Extent of wound closure is typically observed through light microscopy and protein expression patterns that occur during the wound healing process can be characterized by immunofluorescence.

Advances in understanding the repair mechanisms of wounded cell monolayers have been facilitated by the development of methods for performing the assay in multiwell plates. Such methods include delivering wounds to existing monolayers in the wells, or occlusion of the center of the well during monolayer formation to create a gap<sup>3-5</sup>. Each method then involves quantification of the extent of cell migration into the gap.

However, these methods are not optimal for biochemical analysis of the molecular events mediating wound repair. For example, the small scale of a basic pipette tip-derived wound provides insufficient and inadequate material for biochemical analysis. Using the same pipette tip to scratch a cell monolayer in a larger plate is tedious and irreproducible, and the proportion of migrating cells to quiescent cells is low. Several methods have been described to scale up the scratch assay by creating multiple wounds in cell monolayers but these methods require specialized tools<sup>6,7</sup>.

EMD Millipore has developed the Cell Comb™ Scratch Assay to address the need for an easy-to-use tool for creating multiple scratch wounds. The Cell Comb™ has been optimized to apply a high density field of scratches to maximize the area of wound edges, while leaving sufficient numbers of undamaged cells to migrate into the gap. This form of high density wounding creates a high proportion of migrating cells to quiescent monolayer cells, which permits sensitive detection of the biochemical events occurring, specifically in the migrating cell population.



*For Research Use Only; Not for use in diagnostic procedures  
Patent Pending*

---

## Kit Components

1. Cell Combs (Part No. CS210543): Quantity of 6 individually packaged, disposable combs.
2. Rectangular Cell Culture Plates (Part No. CS210544): Quantity of 6 individually packaged, cell culture-treated 86 mm x 128 mm plates.

The Cell Combs and Cell Culture Plates have been subjected to E-beam irradiation in order to minimize the possibility of contamination.

---

## Materials Not Supplied

- |                                                                                                                                         |                                                                      |
|-----------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------|
| 1. Sterile cell culture hood                                                                                                            | 6. Hemocytometer (e.g. Scepter™ Handheld Automated Cell Counter)     |
| 2. Pipettes, liquid aspirators, etc. for handling of cells and liquid reagents                                                          | 7. Trypan blue or equivalent viability stain                         |
| 3. Sterile plastic ware (cell culture flasks, centrifuge tubes, pipettes, pipette tips, etc. for handling of cells and liquid reagents) | 8. Low speed centrifuge for cell harvesting                          |
| 4. Sterile Dulbecco's phosphate-buffered saline (DPBS), without calcium or magnesium                                                    | 9. CO <sub>2</sub> tissue culture incubator                          |
| 5. Cell type of interest, with appropriate growth medium and cell detachment buffer (e.g., 0.25% trypsin)                               | 10. Cell lysis buffer, compatible with biochemical assay of interest |
|                                                                                                                                         | 11. Cell scraper                                                     |

---

## Storage & Usage

Store Cell Combs and Cell Culture Plates at room temperature. Use within **1 year** from date of receipt.

**Before scratching a cell monolayer, practice running the comb across a clean plate.**

---

## Assay Protocol

1. Plate 10 mL cells of interest, suspended in growth media at a density that will yield a confluent monolayer within 1-2 days after plating, into the provided rectangular plates. Put the plates into an incubator with appropriate settings (usually 37°C, 5% CO<sub>2</sub>).  
*Note: Use caution to keep plates level while transporting in and out of the incubator, as rectangular plates are more prone to spillage of media than circular culture dishes.*
2. When cells are confluent, remove media from the plate.
3. In a laminar flow hood, use scissors to open the pouch containing the Cell Comb. Remove Cell Comb™ from pouch, while handling only the top folded edge, and grasp the Cell Comb™ at the top folded edge, as shown in **Figure 1-A**.
4. Align the guide groove with the long edge of the plate, and draw the teeth of the comb across the monolayer, as shown in **Figure 1-B**. Apply sufficient force to apply even scratches on the monolayer, but **avoid excessive force that may bend the teeth**.


5. (Optional) To create scratches in 2 directions, align the guide groove with the short edge of the plate, and draw the teeth of the comb across the monolayer, perpendicularly to the first set of scratches, as shown in **Figure 1-C**
6. Discard comb after scratching one plate. Combs are designed for single use on one plate.  
**Figure 2** shows an example of the pattern of scratches made with the Cell Comb™.
7. Wash the scratched monolayer 1-2 times with media to remove detached cells, and leave 10 mL media in the plate.
8. Perform steps 2-8 for each additional plate of cells. It is recommended to include an unscratched plate of cells as a control.
9. Return the plates to incubator for desired time.

**Figure 3** depicts the time course of several key biochemical events, including Rac activation and phosphorylation of Erk and FAK.

10. After desired time period, remove the media and gently wash the cells with PBS. Collect or lyse cells in a manner appropriate for the intended biochemical analysis.

---

## Example Data




**Figure 1. Usage Schematics**

**A)** Removing the comb from the pouch.


**B)** Holding comb and drawing across a plate in the long direction (one-direction scratches).

**C)** Holding comb and drawing across a plate in the short direction (two-direction scratches).



**Figure 2. Schematic and photos of scratch pattern on a cell monolayer.**

*Top*, drawing of scratch pattern in one-direction (*left*) and two-directions (*right*). *Middle*, 2X magnification view of cell monolayers with multiple wounds applied in one direction (*left*) or two directions (*right*) with the Cell Comb™. *Bottom*, phase contrast image (10X) of scratched cell monolayers immediately after wounding.



**Figure 3. Biochemical analysis of cytoskeletal signaling after wounding with the Cell Comb™.** A confluent monolayer of NIH3T3 cells in the Rectangular Cell Culture Plate was left intact or wounded with a Cell Comb™, and lysed after the indicated time. Lysates were analyzed by Western blot for phospho-FAK (Y397) (EMD Millipore Catalogue No. 05-1140), total p130Cas (EMD Millipore Catalogue No. 06-500), phospho-Erk (T202/Y204 and T185/Y187) (EMD Millipore Catalogue No. 05-797R), total Erk (EMD Millipore Catalogue No. 05-1152), and total GAPDH (EMD Millipore Catalogue No. MAB374).

Major increases in phosphorylated Erk are observed to occur in a bimodal fashion at 30 min and 4 hours after wounding; such a bimodal increase has been previously described<sup>6</sup>. At 2-4 hours post-wounding, increases in phosphorylated FAK and total pCas are observed.

In addition, lysates were analyzed for the GTP-bound form of Rac1 GTPase by pulldown assay with immobilized PAK1 (EMD Millipore Catalogue No. 17-441). In quiescent cells, activated Rac1 is undetectable, whereas within 30 minutes of wounding, a large increase in activated Rac1 is evident. In sum, major increases in cytoskeletal signaling events, such as Rac1 activation and Erk phosphorylation, are detected by biochemical analysis in cells wounded with the Cell Comb™ Scratch Assay.

---

## References

1. Liang, C.-C. *et al.* (2007) In vitro scratch assay: a convenient and inexpensive method for analysis of cell migration in vitro. *Nature Protoc.* 2: 329-333.
2. Cory, G. (2011) Scratch-wound assay. *Meth. Mol. Biol.* 769: 25-30.
3. Yarrow, J.C. *et al.* (2004) A high-throughput cell migration assay using scratch wound healing, a comparison of image-based readout methods. *BMC Biotechnol.* 4: 21.
4. Simpson K.J. *et al.* (2008) Identification of genes that regulate epithelial cell migration using an siRNA screening approach. *Nat. Cell Biol.* 9: 1027-1038.
5. Gough, W. *et al.* (2011) A quantitative, facile, and high-throughput image-based cell migration method is a robust alternative to the scratch assay. *J. Biomol. Screen.* 16: 155-163.
6. Turchi, L. *et al.* (2002) Dynamic characterization of the molecular events during in vitro epidermal wound healing. *J. Invest. Dermatol.* 119: 56-63.
7. Lauder, H. *et al.* (1998) Quantification of the repair process involved in the repair of a cell monolayer using an in vitro model of mechanical injury. *Angiogenesis.* 2: 67-80.

---

## Warranty

**Millipore Corporation** ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the Company all monies paid for such applicable Product. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY COMPANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

Copyright ©2008-2012 EMD Millipore Corporation, a division of Merck KGaA, Darmstadt, Germany. All rights reserved. EMD, EMD Millipore, Millipore, the M mark, Upstate, Chemicon, Linco and all other trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany.

This Page Left Intentionally Blank

This Page Left Intentionally Blank

