

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Trizma® base

Catalog Number **T4661**Store at Room Temperature

CAS RN 77-86-1

Synonyms: Tris base, Tris(hydroxyamino)methane, 2-Amino-2-(hydroxymethyl)-1,3-propanediol, Tris(hydroxymethyl)aminomethane, THAM, Trometamol Note: The "Tris" described in this document is **not** the "Tris" used to flame-proof fabric. That compound, Tris(2,3-dibromopropyl)phosphate, has been reported to be a cancer suspect agent.

Molecular Formula: C₄H₁₁NO₃ Molecular Weight: 121.14

pK_a (20 °C):¹ 8.3

Melting Point: 171-172 °C

Product Description

Trizma[®] is the registered trademark of Sigma-Aldrich applied to various compounds of Tris(hydroxymethyl)-aminomethane that are prepared by Sigma-Aldrich.

This product meets all specifications for Catalog Number T1503, and all analytical and packaging requirements per the current United States Pharmacopoeia (USP). This product is not pyrogen-tested. Certificates of Analysis are available with lot-specific information.

Tris(hydroxyamino)methane, or "Tris" for short, is an established basimetric standard and buffer used in biochemistry and molecular biology.

It may be used by itself as a buffer or as a component of mixed buffer formulations.

These different buffer formulations include:

- Tris-EDTA (TE) buffer
- Tris magnesium buffer
- Tris-acetate-EDTA (TAE) buffer
- Tris-borate-EDTA (TBE) buffer
- Tris-buffered saline (TBS)
- Tris-buffered saline with dextrose (TBS-D)
- Tris-glycine buffer
- Tris-phosphate EDTA buffer
- Tris-SDS buffer
- Tris-sucrose
- Tris-Tricine-SDS buffer

When preparing Trizma solutions at a given pH and temperature, it is necessary to choose the proper mixture of Trizma free base and a corresponding Trizma salt to give the desired final pH at the desired temperature.

Trizma has a significant temperature coefficient, which affects the pH of the solution. For a given concentration the following changes are observed:

- From 5 °C to 25 °C, the pH decreases an average of 0.03 pH units per °C.
- From 25 °C to 37 °C, the pH decreases an average of 0.025 pH units per °C.

This product has been used in various studies and application fields, including:

- Transport properties of mucus³
- Enzymatic assays from cell extracts⁴
- Collagen analysis⁵
- Mitochondria preparation⁶

Precautions and Disclaimer

For research or further manufacturing use only. Not intended for direct use in humans or animals. TSCA-exempt use only in US. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

This product is soluble in water (666 mg/mL), yielding a clear, colorless solution.

Storage/Stability

Trizma solutions can be autoclaved.

References

- Gomori, G., Preparation of Buffers for Use in Enzyme Studies. *Methods Enzymol.*, 1, 138-146 (1955).
- Sambrook, J., and Russell, D.W., Molecular Cloning: A Laboratory Manual (3rd ed.). Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), pp. 5.8, 5.30, 5.43, 5.60, 5.76, 10.25-10.26, 12.75, 12.84, 12.87, 13.52, 16.29-16.31, A1.2-1.3, A1.7-A1.8, A1.17-A1.18, A1.22, A8.42-A8.43 (2001).
- 3. Odziomek, M. et al., Chem. Process Eng., **38(2)**, 217-229 (2017).

- 4. Çağlayan, M. et al., Bio. Protoc., **5(11)**, e1493 (2015).
- Griffiths, M.M. et al., "Collagen-Induced Arthristis in Rats", in *Methods in Molecular Biology: Arthritis* Research – Methods and Protocols, Volume 2 (A.P. Cope, ed.). Humana Press (Totowa, NJ), Chapter 15, pp. 201-214 (2007).
- 6. Guo, R. et al., Cell, 170(6), 1247-1257 (2017).

Trizma is a registered trademark of Sigma-Aldrich Co., LLC.

GCY,MBK,RXR,MAM 12/18-1