

Quick Start Guide

GenElute™ Mammalian Genomic DNA Kit

G1N10, G1N70, G1N350

Reagents to Prepare

- Dilute the Wash Solution Concentrate with 10 mL (10 preps), 80 mL (70 preps) or 360 mL (350 preps) of 95–100% ethanol. After each use, tightly cap the diluted Wash Solution to prevent ethanol evaporation.
- Dissolve one bottle of Proteinase K in water to obtain a 20 mg/mL stock solution. This solution can be stored for several days at 2–8 °C or for long term storage, aliquots may be stored at –20 °C until needed.

Protocol

All spins at $\geq 6,500 \times g$ unless otherwise noted.

Release DNA from Cells

1. Cultured Cells
 - a. Pellet up to 5×10^6 cells by spinning at $300 \times g$. Discard medium.
 - b. Resuspend cells in 200 μ L Resuspension Solution. **Optional:** Add 20 μ L RNase A. Mix & incubate at room temperature for 2 min.
 - c. Add 20 μ L Proteinase K, followed by 200 μ L Lysis Solution C. Vortex thoroughly. Incubate at 70°C for 10 min.
 - d. Continue to **Bind DNA to Column**.

Release DNA from Tissues

1. Mammalian Tissues
 - a. Mince up 25 mg tissue on ice. Transfer to microcentrifuge tube.
 - b. Add 180 μ L Lysis Solution T to tissue & 20 μ L Proteinase K. Vortex or invert to mix. Incubate at 55°C until fully digested (2–4 h). **Optional:** After digestion add 20 μ L RNase A. Mix and incubate at room temperature for 2 min.
 - c. Continue to **Bind DNA to Column**.

Bind DNA to Column

2. Column Preparation: Add 500 μ L Column Preparation Solution to assembled column in collection tube and spin at $12,000 \times g$ for 1 minute. Discard flowthrough.
3. Add 200 μ L ethanol to cell lysate. Vortex or invert to mix.
4. Transfer to column and spin for 1 min. Discard collection tube with flowthrough.

Wash to remove contaminants

5. Transfer column to new collection tube. Add 500 μ L diluted Wash Solution to column. Spin 1 min and discard flowthrough. Repeat once.
6. Spin 3 min to completely dry the column.

Elute purified DNA

7. Transfer column to new collection tube.
8. Add 200 μ L Elution Solution to column. Spin 1 min. Repeat elution in the same or new tube if desired.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/Offices.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/TechService.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/Terms.

MilliporeSigma, Millipore, GenElute and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.