

Product Information

10X PCR Buffer

Catalog Number **P2192**

Product Summary

DNase, RNase: None detected
Suitable for use in the Polymerase Chain Reaction (PCR) when diluted to 1X concentration (see comments section).

Composition of 10X concentrate

100 mM Trizma®-HCl, pH 8.3 at 25 °C
500 mM KCl
15 mM MgCl₂
0.01% (w/v) gelatin

Comments

Information about licenses to PCR can be obtained from The Perkin-Elmer Corporation or Roche Molecular Systems, Inc.

AmpliTaQ® is a registered trademark of Roche Molecular Systems, Inc.

PCR Suitability

10X PCR Buffer was tested at a final concentration of 1X (10 mM Trizma®-HCl, pH 8.3 at 25 °C, 50 mM KCl, 1.5 mM MgCl₂, 0.001% (w/v) gelatin), in a reaction containing each dNTP at 200 µM, primers defining an approximately 500 base pair region of λDNA at 1.0 µM each, λ DNA template at 1 ng/100 µl, and mPliTaQ® DNA Polymerase at 2.5 Units/100 µl. The reaction underwent 25 cycles of 94 °C to denature the double stranded DNA, 55 °C to anneal the DNA segments, and 72 °C to extend the DNA segments. A single band of approximately 500 base pairs was visualized following electrophoresis of the reaction product in a 1.5% agarose gel.

Endonuclease-exonuclease

One µg of λ Hind III fragments was incubated for 16 hours at 37 °C with 10X PCR Buffer at a final concentration of 1X in a 50 µl reaction mixture containing 30 mM Trizma®-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl₂. No degradation of the DNA fragments was detected by agarose gel electrophoresis. Detection limit: Degradation of 10% of the DNA substrate is detectable.

Endonuclease (Nickase)

One µg of pBR322 DNA was incubated for 16 hours at 37 °C with 10X PCR Buffer at a final concentration of 1X in a 50 µl reaction mixture containing 30 mM Trizma®-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl₂. No conversion of the covalently closed circular DNA to the nicked or linear form was observed by agarose gel electrophoresis. Detection limit: Conversion of 1% of the DNA substrate is detectable.

RNase

Two µg of transfer RNA were incubated for 16 hours at 37 °C with 10X PCR Buffer at a final concentration of 1X in a 50 µl reaction mixture containing 30 mM Trizma®-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl₂. No degradation of the tRNA was detected by polyacrylamide gel electrophoresis. Detection limit: Degradation of 10% of the tRNA substrate is detectable.

PHC 11/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.