1.14846.0001

MQuant®

Phosphate Test

P

for the determination of orthophosphate

1. Method

Determination with color-disk comparator

In sulfuric solution orthophosphate ions react with molybdate ions to form molybdophosphoric acid. Ascorbic acid reduces this to phosphomolybdenum blue (PMB). The phosphate concentration is measured **semiquantitatively** by visual comparison of the color of the measurement solution with the color fields of a color disk.

2. Measuring range and number of determinations

Measuring range / color-scale graduation ¹⁾	Number of determinations
0.2 - 0.4 - 0.6 - 0.8 - 1.0 - 1.5 - 2.0 - 2.5 - 3.0 mg/l PO ₄ -P	200
0.6 - 1.2 - 1.8 - 2.5 - 3.1 - 4.6 - 6.1 - 7.7 - 9.2 mg/l PO ₄ ³⁻	
0.46 - 0.9 - 1.4 - 1.8 - 2.3 - 3.4 - 4.6 - 5.7 - 6.9 mg/l P ₂ O ₅	

¹⁾ for conversion factors see section 8

3. Applications

This test measures only orthophosphate.

Sample material:

Groundwater and surface water, seawater
Drinking water and mineral water
Waters from aquaculture
Industrial and process water
Wastewater
Soils and fertilizers after appropriate sample pretreatment
Food after appropriate sample pretreatment
Steels after appropriate sample pretreatment

4. Influence of foreign substances

This was checked individually in solutions containing 1 and 0 mg/l PO₄-P. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag ⁺	1000	F ⁻	50	Pb ²⁺	25
AsO ₄ ³⁻	0.2	Fe ³⁺	1000	S ²⁻	2.5
Ca ²⁺	1000	Hg ²⁺	10	SiO ₃ ²⁻	1000
Cd ²⁺	1000	Mg ²⁺	1000	SO ₃ ²⁻	1000
CN ⁻	1000	Mn ²⁺	1000	Zn ²⁺	1000
Cr ³⁺	1000	NH ₄ ⁺	1000	EDTA	1000
Cr ₂ O ₇ ²⁻	5	Ni ²⁺	500	Na-acetate	1%
Cu ²⁺	250	NO ₂ ⁻	1000	NaCl	5%
				NaNO ₃	10%
				Na ₂ SO ₄	10%

Reducing agents interfere with the determination.

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

2 bottles of reagent PO₄-1
1 bottle of reagent PO₄-2
1 graduated 6-ml plastic syringe
2 test tubes with screw caps
1 color-disk comparator

Other reagents and accessories:

MQuant® Phosphate Test, Cat. No. 1.10428
measuring range 10 - 500 mg/l PO₄³⁻ (3.3 - 163 mg/l PO₄-P)
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
Phosphate standard solution Certipur®, 1000 mg/l PO₄³⁻, Cat. No. 1.19898
Hydrochloric acid 25% for analysis EMSURE®, Cat. No. 1.00316
Sodium hydroxide solution 1 mol/l (approx. 4%) Titripur®, Cat. No. 1.09137

MQuant® Flat-bottomed tubes with screw caps for MQuant® tests with color disk comparator (12 pcs), Cat. No. 1.17988

Refill pack:

Cat. No. 1.18465

Phosphate Test

Refill pack for 1.14846

(Reagents **without technical accessories** for the number of determinations stated in section 2)

6. Preparation

- Analyze immediately after sampling.
- Check the phosphate content with the MQuant® Phosphate Test. Samples containing more than 3.0 mg/l PO₄-P must be diluted with distilled water.
- The pH must be within the range 0 - 10.** Adjust, if necessary, with sulfuric acid.
- Filter strongly turbid samples.

7. Procedure

	Measurement sample right-hand tube (A) behind the color disk	Blank left-hand tube (B) behind the color disk	
Pretreated sample (10 - 35 °C)	6 ml	6 ml	Inject into the test tube with the syringe. Add, close the tube, and shake vigorously until the reagent is completely dissolved.
Reagent PO ₄ -1	5 drops ¹⁾	- and mix.	
Reagent PO ₄ -2	1 level blue microspoon (in the cap of the PO ₄ -2 bottle)	-	

Leave to stand for 1 min (reaction time).

Hold the comparator to the light, keeping it upright, and rotate the disk until the closest possible color match is achieved between the two large windows. Read off the result in mg/l PO₄-P shown in the small window.

¹⁾ Hold the bottle vertically while adding the reagent!

Notes on the measurement:

- The color of the measurement solution remains stable for at least 60 min after the end of the reaction time stated above.
- Turbidity in the measurement solution makes the color comparison more difficult.
- If the color of the measurement solution is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 3.0 mg/l PO₄-P is obtained. Concerning the result of the analysis, the dilution (see also section 6) must be taken into account:

$$\text{Result of analysis} = \text{measurement value} \times \text{dilution factor}$$

8. Conversions

required given	mg/l PO ₄ -P	mg/l PO ₄ ³⁻	mg/l P ₂ O ₅
1 mg/l PO ₄ -P	1	3.07	2.29
1 mg/l PO ₄ ³⁻	0.326	1	0.747
1 mg/l P ₂ O ₅	0.436	1.34	1

9. Method control

To check test reagents, measurement device, and handling: Dilute the phosphate standard solution with distilled water to 0.8 mg/l PO₄-P (2.5 mg/l PO₄³⁻) and analyze as described in section 7. Additional notes see under www.qa-test-kits.com.

10. Notes

- Reclose the reagent bottles immediately after use.
- Use only phosphate-free detergents to rinse the test tubes and the syringe. Otherwise fill with hydrochloric acid (approx. 10%) and leave to stand for several hours.
- Cleanse the test tubes from time to time as follows: Fill with sodium hydroxide solution (approx. 0.4%) and leave to stand for max. 1 hour.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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© 2024 Merck KGaA, Darmstadt, Allemagne et/ou ses sociétés affiliées. Tous droits réservés. Merck, Supelco, Sigma-Aldrich et MQuant sont des marques de Merck KGaA, Darmstadt, Allemagne. Toutes les autres marques citées appartiennent à leurs propriétaires respectifs. Des informations détaillées sur les marques sont disponibles via des ressources accessibles au public.