Instructions for Use

Leukocyte α -Naphthyl Acetate (Non-Specific Esterase) and Leukocyte Naphthol AS-D Chloroacetate (Specific Esterase) Kit

Procedure No. 91

Intended Use

Sigma-Aldrich Naphthol AS-D Chloroacetate Esterase and α -Naphthyl Acetate Esterase kits are intended as a general laboratory use histological stain. Naphthol AS-D chloroacetate esterase and α -naphthyl acetate esterase reagents are for professional "In Vitro Diagnostic Use" only. This manual, qualitative, histochemical procedure demonstrates naphthol AS-D chloroacetate esterase and α -naphthyl acetate esterase activity in leukocytes of human blood, bone marrow films, or tissue touch preparations. Histological visualization of leucocyte cellular esterases is a unique technique currently commonly used in medicine.

Cellular esterases are ubiquitous, apparently representing a series of different enzymes acting upon select substrates. Under defined reaction conditions, it may be possible to determine hemopoietic cell types, using specific esterase substrates. The described methods provide hematologists and hematopathologists means of distinguishing granulocytes from monocytes.¹⁻⁷

To perform the test, blood, bone marrow films or tissue touch preparations are incubated with either naphthol AS-D chloroacetate (NCAE) or α -naphthyl acetate (NAE) in the presence of freshly formed diazonium salt. Enzymatic hydrolysis of ester linkages liberates free naphthol compounds. These couple with the diazonium salt, forming highly colored deposits at sites of enzyme activity.

Most recent procedures, including those provided by Sigma-Aldrich, employ stable diazonium salts. These are formed by reacting an arylamine with sodium nitrite in an acid medium.⁸ The resulting diazonium chloride (usually unstable) can then be treated with compounds such as zinc chloride, zinc sulfate or naphthalene-1-6-disulfonate, forming stable salts. These stabilizers may exert marked inhibition on some enzymatic systems, whereas the diazonium chlorides are less inhibitory.⁹ For this reason, Sigma-Aldrich now provides stable solutions for Fast Red Violet LB Base, Fast Blue BB Base and Sodium Nitrite for esterase cytochemistry. To further simplify these methods, stable solutions of Naphthol AS-D Chloroacetate and α -Naphthyl Acetate are included. The availability of these stable solutions allows the customer to adjust working reagent volumes according to needs, eliminating waste.

Reagents

Naphthol AS-D Chloroacetate Solution (Cat. No. 911-10ML)

Naphthol AS-D chloroacetate, 8 mg/mL, and stabilizer. Warning. Causes eye irritation. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Fast Red Violet LB Base Solution (Cat. No. 912-10ML)

Fast red violet LB base, 15 mg/mL, in 0.4 mol/L hydrochloric acid with stabilizer. Danger. Causes severe skin burns and eye damage. Wear protective gloves/eye protection/face protection.

TRIZMAL™ 6.3 Concentrate (Cat. No. 913-50ML)

TRIZMA® maleate, 1 mol/L, with surfactant. pH 6.3 $\pm$ 0.15 at 25°C. Danger. Causes skin irritation. May cause allergic skin reaction. Causes serious eye damage. Wear protective gloves/eye protection/face protection.

Sodium Nitrite Solution (Cat. No. 914-10ML)

Sodium nitrite, 0.1 mol/L

Citrate Solution (Cat. No. 915-50ML)

Citric acid, 18 mmol/L, sodium citrate, 9 mmol/L, sodium chloride, 12 mmol/L, with surfactant. pH 3.6 $\pm$ 0.1 at 25°C

α -Naphthyl Acetate Solution (Cat. No. 916-10ML)

α -Naphthyl acetate, 12.5 mg/mL, in methanol solution with stabilizers. Danger. Flammable liquid and vapour. Toxic if swallowed. Harmful in contact with skin or if inhaled. Causes mild skin irritation. Causes serious eye irritation. Causes damage to organs.

Fast Blue BB Base Solution (Cat. No. 917-10ML)

Fast blue BB base, 15 mg/mL, C.I. 37175 in 0.4 mol/L hydrochloric acid with stabilizers. Warning. Harmful if swallowed. Causes mild skin irritation. Causes serious eye irritation. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

TRIZMAL™ 7.6 Concentrate (Cat. No. 918-50ML)

TRIZMA® maleate, 1 mol/L, with surfactant. pH 7.6 $\pm$ 0.15 at 25°C. Danger. Causes skin irritation. Causes serious eye damage. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause respiratory irritation. Harmful to aquatic life.

Hematoxylin Solution Gill No. 3 (Cat. No. GHS3-50ML)

Certified hematoxylin, 6.0 g/L, C.I. 75290, sodium iodate, 0.6 g/L, and aluminum sulfate, 52.8 g/L, with stabilizers. Danger. Harmful if swallowed. Causes serious eye damage. Wear protective gloves/eye protection/face protection.

Sodium Fluoride Solution (Cat. No. 919-25ML)

Sodium fluoride, 20 g/L. Warning. Harmful if swallowed. Causes mild skin irritation.

Special Materials Required but Not Provided

- Acetone, ACS Reagent
- Formaldehyde, 37%, ACS
- Sodium Fluoride Solution, Cat. No. 919-25ML (Required for α -Naphthyl Acetate Esterase with Fluoride Inhibition Procedure.)

Storage and Stability

Store Hematoxylin Solution Gill No. 3 at room temperature (18–26°C) protected from light. Store other reagents in refrigerator (2–8°C).

TRIZMAL™ 6.3 Concentrate, TRIZMAL™ 7.6 Concentrate and Citrate Solution are suitable for use in the absence of microbial growth. Other reagents are stable until expiration date shown on the labels.

Deterioration

Discard TRIZMAL™ Concentrate and Citrate Solution if microbial growth is evident. Discard Hematoxylin Solution Gill No. 3 if solution turns brown (over-oxidized from air) or purple (loss of acidity).

Preparation

Warm all reagents to room temperature (18–26°C) before use. Esterase reagents are provided ready for use.

Citrate-Acetone-Formaldehyde Fixative: To 25 mL Citrate Solution, add 65 mL Acetone and 8 mL 37% Formaldehyde. Place in glass bottle and cap tightly. Store in refrigerator (2–8°C). Bring to room temperature (18–26°C) prior to use. Stable up to 4 weeks if stored tightly capped in refrigerator.

Precautions

The IVDs included in these kits are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Blood, bone marrow films, tissue touch preparation and cytocentrifuge preparations may be used with both α -naphthyl acetate esterase and naphthol AS-D chloroacetate esterase. Either EDTA or heparin will serve as an anticoagulant.⁹ Frozen and paraffin embedded tissues may be used with naphthol AS-D chloroacetate esterase. α -Naphthyl acetate esterase may be used successfully on frozen tissue sections.¹⁰ Blood or bone marrow films may be stored fixed at room temperature (18–26°C) for several weeks or unfixed for several days without appreciable change in activity.^{4,9} Do not ship whole blood for assay at other laboratories. Send fixed or unfixed slides. Slides should be kept cool during transit. Allow films to dry at least 1 hour prior to fixation.

Limitations of the Procedure

The described procedures are performed at 37°C. If reagents are not at this temperature, weak or negative reactions may be obtained. It is recommended that temperatures be checked with an accurate thermometer. Controlled temperature water baths are more efficient than warm air incubators and should be used for enzyme cytochemical methods. Heat transfer through glass is more rapid than through plastic, thus, glass Coplin jars should be employed.

Many enzyme systems are sensitive to minute traces of detergent. Washing glassware with dilute bleach followed by rinsing in copious quantities of deionized water will prevent detergent effect upon cellular enzymes.

The fixative system described in the "Reagents" section contains formaldehyde. If fixative is not completely removed by vigorous rinsing, any amount of aldehyde, either concentrated on air-dried slides or trace amounts from wet slides, added to the incubation system may result in enzyme inhibition. To prevent loss of blood film during the rinse process, direct the water stream onto the slide above the feather edge. Rinse both sides of the slide.

Results are based on a certain degree of subjective interpretation. Individual laboratories should establish their own normal ranges.

Notes

Users of Sigma Kits 390A and 91C should exercise caution when using TRIZMA® Buffer Concentrates Catalog Nos. 913 and 903C as they are not interchangeable. Use of incorrect buffer will result in a negative reaction.

Procedure

The described procedures are performed at 37°C.

Naphthol AS-D Chloroacetate Esterase Procedure

1. Prewarm sufficient deionized water for substrate use to 37°C. Check temperature before use.
2. Immediately prior to fixation, add 1 mL Sodium Nitrite Solution to 1 mL Fast Red Violet LB Base Solution in a test tube. Mix gently by inversion and allow to stand 2 minutes. Active evolution of gas bubbles should be avoided.
3. Add solution from Step 2 to 40 mL prewarmed deionized water from Step 1.
4. Add 5 mL TRIZMAL™ 6.3 Buffer Concentrate. (See "Note")
5. Add 1 mL Naphthol AS-D Chloroacetate Solution. The solution should turn red. Mix well and pour into Coplin jar.
6. Bring Citrate-Acetone-Formaldehyde (CAF) Solution to room temperature (18–26°C). Fix slides by immersing in CAF solution for 30 seconds.
7. Rinse slides thoroughly in running deionized water for 45–60 seconds, then place in solution from Step 5. Do not allow slides to dry.
8. Incubate for 15 minutes, at 37°C protected from light.
9. After 15 minutes, remove slides and rinse thoroughly in deionized water for at least 2 minutes.
10. Counterstain 2 minutes in Hematoxylin Solution, Gill No. 3.
11. Rinse in tap water and air dry.
12. Evaluate microscopically. If coverslipping is required use only an aqueous mounting media.

Notes

- For use with Columbia jars, divide reagent volumes by 5.
- If substrate (See Step 5) appears turbid bring to room temperature (18–26°C) and mix well.
- If slides have been prefixed and stored, skip fixation (Steps 6 and 7) and begin staining of dry, prefixed slides at Step 8.

α-Naphthyl Acetate Esterase Procedure

1. Prewarm sufficient deionized water for substrate use to 37°C. Check temperature before use.
2. Immediately prior to fixation, add 1 mL Sodium Nitrite Solution to 1 mL Fast Blue BB Base Solution in a test tube. Mix by inversion and allow to stand at least 2 minutes. The color will change from dirty brown to deep yellow. Active evolution of gas bubbles should be avoided.
3. Add solution from Step 2 to 40 mL prewarmed deionized water from Step 1.
4. Add 5 mL TRIZMAL™ 7.6 Buffer Concentrate.
5. Add 1 mL α-Naphthyl Acetate Solution. The solution should turn greenish. Mix well and pour into Coplin jar.
6. Bring Citrate-Acetone-Formaldehyde (CAF) Solution to room temperature (18–26°C). Fix slides by immersing in CAF solution for 30 seconds. Agitate slides vigorously for the last 5 seconds.
7. Rinse slides thoroughly in running deionized water for 45–60 seconds, then place in solution from Step 5. Do not allow slides to dry.
8. Incubate for 30 minutes at 37°C protected from light.
9. After 30 minutes, remove slides and rinse thoroughly for at least 2 minutes in running deionized water.
10. Counterstain 2 minutes in Hematoxylin Solution, Gill No. 3.
11. Rinse in tap water and air dry.
12. Evaluate microscopically. If coverslipping is required use only an aqueous mounting media.

Notes

- For use with Columbia jars, divide reagent volumes by 5.
- If substrate (See Step 5) appears turbid bring to room temperature (18–26°C) and mix well.
- If slides have been prefixed and stored, skip fixation (Steps 6 and 7) and begin staining of dry, prefixed slides at Step 8.

Double Staining Esterase Procedure

1. Perform α-Naphthyl Acetate Esterase test as described in Procedure. Do not counterstain.
2. Rinse slide 5 minutes in deionized water.
3. Perform Naphthol AS-D Chloroacetate Esterase test as described in Procedure Steps 1-12. Omit Step 6.

Note

Fast Blue BB Base Solution may be substituted for Fast Red Violet LB Base Solution if blue granulation is preferred for Naphthol AS-D Chloroacetate Esterase.

α-Naphthyl Acetate Esterase with Fluoride Inhibition Procedure

Although α-naphthyl acetate esterase is found primarily in cells of monocytic lineage when performed as described, it should be recognized that megakaryocytes and erythroid precursors are positive for this enzyme.¹¹ Lymphocytes and some mature granulocytes also show occasional positivity.⁴ To differentiate these cells conclusively from monocytes, sodium fluoride is incorporated with the incubation system. The monocyte enzyme is inactivated in the presence of this compound.¹² The following procedure may be used to perform the fluoride inhibition test.

1. To 2 mL Fast Blue BB Base Solution add 2 mL Sodium Nitrite Solution. Mix gently by inversion. Allow to stand 2 minutes.
2. Label 2 beakers A and B, and add the following:

	Beaker A	Beaker B
Prewarmed 37°C deionized water	40 mL	40 mL
Diazotized Fast Blue BB from Step 1	2 mL	2 mL
TRIZMAL™ 7.6 Concentrate	5 mL	5 mL
α-Naphthyl Acetate Solution	1 mL	1 mL
Sodium Fluoride Solution	-	1 mL

3. Mix well and pour into Coplin jars labeled A and B.
4. Proceed as described in Steps 6-12 of α-Naphthyl Acetate Esterase Procedure.

Performance Characteristics

Method of Scoring

Scan the film and select a thin area with few erythrocytes. Sites of Naphthol AS-D Chloroacetate Esterase activity will appear as bright red granulation, α-Naphthyl Acetate Esterase as black granulation. Rate from 0 to 4+ on the basis of quantity and intensity of individual dyes within the cytoplasm of the respective cell types. Characteristics of scoring are based on somewhat subjective interpretation. A suggested scoring format is presented in Table 1. Conclusions center on relative presence or absence of staining.

Table 1. Characteristics of Scoring

Cell Rating	Staining Intensity	Interpretation
0+	None	–
1+	Faint to Moderate	±
2+	Moderate to Strong	+
3+	Strong	+
4+	Brilliant	+

Expected Observations

Naphthol AS-D Chloroacetate Esterase

(Fast Red Violet LB) - Enzyme is usually considered specific for cells of granulocytic lineage. Sites of activity show bright red granulation. Activity is weak or absent in monocytes and lymphocytes.

α-Naphthyl Acetate Esterase

(Fast Blue BB) - Enzyme is detected primarily in monocytes, macrophages and histiocytes, and is virtually absent in granulocytes. Monocytes should show black granulation. Lymphocytes may occasionally exhibit enzyme activity.

Double Staining Esterase

Specimens taken through the double staining procedure will demonstrate the granulocytes with red granulation and monocytes with black granulation.

α-Naphthyl Acetate Esterase with Fluoride Inhibition

All cells of monocytic lineage will be negative for enzyme activity, with the exception of differentiated histiocytes or specialized macrophages in tissue which may also be resistant to sodium fluoride.¹⁰

Quality Control

The reagent system should be monitored by the use of positive and negative control slides.

Positive control slides may be prepared from specific cell lines known to be positive.

Alternately, anti-coagulated blood from normal specimens (preferably with increased monocyte count if using α-naphthyl acetate esterase procedure) may also be used; however, they will provide less intense staining and will have fewer positive cells.

Known negative patient slides may be used as a negative control. If unavailable, staining a specimen in an incubation mixture with the substrate omitted will give the desired results. However, use of the former is highly recommended.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
911	Naphthol AS-D Chloroacetate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
912	Fast Red Violet LB Base Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
913	TRIZMAL™ 6.3 Concentrate	Osmolality	3 of 3	3 of 3	3 of 3	3 of 3
		pH	3 of 3	3 of 3	3 of 3	3 of 3
914	Sodium Nitrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
915	Citrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
916	α-Naphthyl Acetate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
917	Fast Blue BB Base Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
918	TRIZMAL™ 7.6 Concentrate	Osmolality	3 of 3	3 of 3	3 of 3	3 of 3
		pH	3 of 3	3 of 3	3 of 3	3 of 3
919	Sodium Fluoride Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
GHS3	Hematoxylin Solution, Gill No. 3	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

91A:

H226 Flammable liquid and vapor.

H290 May be corrosive to metals.

H301 Toxic if swallowed.

H312 + H332 Harmful in contact with skin or if inhaled.

H319 Causes serious eye irritation.

H370 Causes damage to organs (Eyes, Central nervous system).

H412 Harmful to aquatic life with long lasting effects.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell.

91C:

- H290 May be corrosive to metals.
- H302 Harmful if swallowed.
- H315 Causes skin irritation.
- H317 May cause an allergic skin reaction.
- H318 Causes serious eye damage.
- H335 May cause respiratory irritation.
- H360FD May damage fertility. May damage the unborn child.
- H373 May cause damage to organs (Kidney) through prolonged or repeated exposure if swallowed.
- P202 Do not handle until all safety precautions have been read and understood.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P301 + P312 IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell.
- P302 + P352 IF ON SKIN: Wash with plenty of water.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P308 + P313 IF exposed or concerned: Get medical advice/ attention.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
7. Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127-130
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 4.0	2014
Rev. 5.0	2016
Rev. 6.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
7. Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127–130
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Bibliografia

1. Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974.
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
7. Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113.
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127-130.
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38.
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111.
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973.

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Referencias

- Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
- Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
- Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
- Brown BA: In Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
- Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
- Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973

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Referencer

1. Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
7. Moloney WC, McPherson K, Fliegeman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127-130
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Referenser

- Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
- Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
- Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, sid. 88–113
- Brown BA: In Hematology: Principles and Procedures. Leas och Febriger, Philadelphia, 1984, sid. 127–130
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, sid. 24, 38
- Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, sid. 34, 111
- Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973

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- Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
- Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
- Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127–130
- Sun T: Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
- Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
- Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973

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Reference

1. Beckmann J, Neth R, Gaedlicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
7. Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127–130
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Referanser

- Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
- Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
- Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
- Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127-130
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
- Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
- Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Referanslar

- Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
- Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
- Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naftol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
- Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127-130
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
- Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
- Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Lähteet

- Beckmann J, Neth R, Gaedicke G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
- Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
- Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127–130
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
- Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
- Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Referenties

1. Beckmann J, Neth R, Gaedick G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127-130
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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Reference

1. Beckmann J, Neth R, Gaedick G, et al: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971
5. Yam LT, Li CY, Wolfe NJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957
7. Moloney WC, McPherson K, Fliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960
8. Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
9. Brown BA: In Hematology: Principles and Procedures. Leas and Febriger, Philadelphia, 1984, pp 127–130
10. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 24, 38
11. Hayhoe FGJ, Flemans RJ: In Color Atlas of Hematological Cytology. John Wiley & Sons, New York, 1982, pp 34, 111
12. Li CY, Lam KW, Lam LT: Esterase in human leukocytes. J Histochem Cytochem 21:1, 1973

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MDSS GmbH