1.14548.0001

Spectroquant® Sulfate Cell Test

SO₄²⁻

USEPA equivalent for wastewater

1. Method

Sulfate ions react with barium ions to form slightly soluble barium sulfate. The resulting turbidity is measured in the photometer (turbidimetric method).

The method is analogous to EPA 375.4, APHA 4500-SO₄²⁻ E, and ASTM D516-16.

2. Measuring range and number of determinations

Measuring range	Number of determinations
5 - 250 mg/l SO ₄ ²⁻	25

For programming data for selected photometers / spectrophotometers see www.sigmaaldrich.com/photometry.

3. Applications

Sample material:

Groundwater and surface water, seawater
Drinking water and mineral water
Wastewater
Nutrient solutions for fertilization
Soils after appropriate sample pretreatment
Water for concrete in the construction industry

4. Influence of foreign substances

This was checked individually in solutions containing 125 and 0 mg/l SO₄²⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag⁺	2	Hg²⁺	1000	S₂O₃²⁻	25
Al³⁺	1000	Mg²⁺	1000	Zn²⁺	1000
Ca²⁺	1000	Mn²⁺	1000		
Cd²⁺	1000	NH₄⁺	1000	Anionic surfactants	5
CN⁻	1000	Ni²⁺	1000	Cationic surfactants	500
CO₃²⁻	1000	NO₂⁻	1000	Nonionic surfactants	1000
Cr³⁺	100	Pb²⁺	100	Na-acetate	10 %
Cr₂O₇²⁻	50	PO₄³⁻	1000	NaCl	10 %
Cu²⁺	1000	S²⁻	10	NaNO₃	0.1 %
F⁻	1000	SiO₃²⁻	1000		
Fe³⁺	1000	SO₃²⁻	50		

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent SO₄-1K
25 reaction cells
1 sheet of round self-adhesive stickers for numbering the cells

Other reagents and accessories:

MQuant® Sulfate Test, Cat. No. 1.10019,
measuring range <200 - >1600 mg/l SO₄²⁻
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sodium hydroxide solution 1 mol/l Titripur®, Cat. No. 1.09137
Hydrochloric acid 1 mol/l Titripur®, Cat. No. 1.09057
Membrane filters 0.45 µm
Spectroquant® CombiCheck 10, Cat. No. 1.14676
Sulfate standard solution 40.0 mg/l SO₄²⁻, Cat. No. 1.25050
Sulfate standard solution 125 mg/l SO₄²⁻, Cat. No. 1.25051

Pipette for a pipetting volume of 5.0 ml

6. Preparation

- Analyze immediately after sampling.
- Check the sulfate content with the MQuant® Sulfate Test. Samples containing more than 250 mg/l SO₄²⁻ must be diluted with distilled water.
- The pH must be within the range 2 - 10.** Adjust, if necessary, with sodium hydroxide solution or hydrochloric acid.
- Filter turbid samples through a 0.45 µm membrane filter.

7. Procedure

Pretreated sample (20 - 40 °C)	5.0 ml	Pipette into a reaction cell and mix.
Reagent SO ₄ -1K	1 level green microspoon (in the cap of the SO ₄ -1K bottle)	Add, close the cell tightly, and shake vigorously until the reagent is completely dissolved.
Leave to stand for exactly 2 min (reaction time), then measure the sample in the photometer.		

Notes on the measurement:

- For turbidimetric measurement the cells must be clean. Wipe, if necessary, with a clean dry cloth.
- The pH of the measurement solution must be within the range 1.0 - 2.0.

8. Analytical quality assurance

recommended before each measurement series
To check the photometric measurement system (test reagent, measurement device, handling) and the mode of working, the sulfate standard solutions (see section 5) or Spectroquant® CombiCheck 10 can be used. Besides a **standard solution** with 100 mg/l SO₄²⁻, CombiCheck 10 also contains an **addition solution** for determining sample-dependent interferences (**matrix effects**).
Additional notes see under www.qa-test-kits.com.
For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Notes

- Reclose the reagent bottle immediately after use.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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