

Technical Bulletin

Gel Filtration Markers Kit for Protein Molecular Weights 12,000-200,000 Da

MWGF200

Product Description

Gel filtration chromatography is an established method for determining the size and molecular mass of proteins. Fractionation is based on the diffusion of molecules into the pores of the resin. Larger proteins do not enter the pores of the resin as readily, but instead pass through the fluid volume of the column faster than smaller proteins. These protein molecules elute from the column in order of decreasing molecular mass.

The molecular mass determination of unknown proteins is made by comparing the ratio of V_e/V_o for the protein in question to the V_e/V_o of protein standards of known molecular mass, where:

- V_e is the elution volume
- V_o is the void volume

The V_o of a given column is based on the volume of effluent required for the elution of a large molecule such as Blue Dextran (molecular mass of ~2,000 kDa, Cat. No. D4772). Plotting the logarithms of the known molecular masses of protein standards versus their respective V_e/V_o values produces a linear calibration curve (such as shown in **Figure 1**).

V_e/V_o is essentially independent of column size and protein concentration. However, V_e/V_o may be dependent on temperature for some proteins. Unreliable molecular masses may be obtained if the protein:

- Forms a complex with the gel
- Contains a large amount of carbohydrate
- Aggregates to larger complexes
- Or dissociates into subunits under the conditions used.¹

The molecular mass of an impure protein may be determined using this procedure if a specific detection test is available for the protein.

The procedure for determining molecular masses using gel filtration chromatography, as outlined in this bulletin, is a modification of published methods.^{1,2} The protein standards in this kit may be suitable for use in other chromatographic systems such as HPLC, although some buffer systems seem to alter the elution volumes of albumin (Cat. No. A8531) and carbonic anhydrase (Cat. No. C7025). The proteins in the MWGF200 Kit have a range of molecular masses from 12.4 kDa to 200 kDa.

Several theses³ and dissertations⁴⁻²¹ cite use of product MWGF200 in their protocols.

Reagent

- Cytochrome c from horse heart, Cat. No. C7150, 10 mg/vial. Approximate Molecular Mass: ~12.4 kDa
- Carbonic Anhydrase from bovine erythrocytes, Cat. No. C7025, 15 mg/vial. Approximate Molecular Mass: ~29 kDa
- Albumin, Bovine Serum, Cat. No. A8531, 50 mg/vial (contains ~0.3% dithiothreitol). Approximate Molecular Mass: ~66 kDa
- Alcohol Dehydrogenase from yeast, Cat. No. A8656, 25 mg/vial. Approximate Molecular Mass: ~150 kDa
- β -Amylase from sweet potato, Cat. No. A8781, 15 mg protein/vial (contains ~15% NaCl, 4% glucose, and 1% dithiothreitol). Approximate Molecular Mass: ~200 kDa
- Blue Dextran, Cat. No. D4772, 50 mg/vial. Approximate Molecular Mass: ~2,000 kDa

Storage/Stability

Store the MWGF200 product at -20 °C.

Procedure

Use of MWGF200 for Gel Filtration Chromatography

1. Buffer and Resin:

- It is recommended to use 50 mM Tris-HCl, pH 7.5, with 100 mM KCl as the equilibration buffer, with a 90 cm × 1.6 cm Sephacryl® S-200-HR (Cat. No. S200HR) column at 2-8 °C.
- For information on resin preparation, column packing, and equilibration, contact Technical Service.

2. Void Volume (V_0) Determination:

- Dissolve the Blue Dextran in equilibration buffer containing 5% glycerol at a concentration of 2 mg/mL.
- This concentration of Blue Dextran will give an A_{280} of ~1.0 in the peak fraction.
- Glycerol is added to increase the density of the solution, but its use is optional.
- The recommended sample volume is less than 2% of the total gel bed volume.
- Carefully apply the Blue Dextran sample to the column (avoid disturbing the gel bed surface) to determine V_0 and to check the column packing.
- Immediately after applying the sample, begin collecting fractions of 0.5-1.5% of the total gel bed volume. The flow rate should be ~7% of the column volume per hour.
- Skewing of the Blue Dextran band represents a fault in the column, although some tailing is normal.
- The leading peak indicates the void volume.
- Determine spectrophotometrically the elution volume for Blue Dextran (V_0 for the column) at 280 nm or 610 nm by measuring the volume of effluent collected from the point of sample application to the center of the effluent peak.
- Note:** Mixing Blue Dextran with kit standards or sample proteins is **not** recommended, since many proteins bind to Blue Dextran.

- Prepare the protein standards and Blue Dextran fresh.²² Occasionally some aggregated protein may appear at the void volume.

3. Elution Volume (V_e) Determination for Protein Standards:

- Dissolve individual protein standards in equilibration buffer containing 5% glycerol (see **Table 1**).
- If, upon reconstitution, any of the protein solutions contain insoluble material, then filter the protein solution through a 0.45 µm or 0.2 µm filter. The loss of protein from this filtration is negligible.
- For a 90 cm × 1.6 cm column, the application of 2.0 mL of individual samples at the recommended concentration (**Table 1**) gives an A_{280} of ~1 in the peak fraction.

Table 1: Recommended Protein Concentrations

Protein	Recommended concentration
Bovine Serum Albumin	10 mg/mL
β-Amylase	5 mg/mL
Alcohol dehydrogenase	4 mg/mL
Carbonic anhydrase	3 mg/mL
Cytochrome c	2 mg/mL

- The following pairs of proteins may be mixed and run together on the columns:
 - Cytochrome c and β-Amylase
 - Carbonic anhydrase and Alcohol dehydrogenase
- Apply the protein standards to the column, using the same sample volume and flow rate as used for the Blue Dextran sample.
- The elution of the standard proteins may be followed by absorbance readings at 280 nm.
- Determine the V_e for the protein standards by measuring the volume of effluent collected from the point of sample application to the center of the effluent peak.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
MWGF200bul Rev 10/21 FF,MAM,KMR,GCY

4. Standard Curve: Plot molecular mass vs. V_e/V_o for each respective protein standard on semi-log paper (see **Figure 1**).
5. Elution Volume (V_e) Determination for an Unknown Protein:
 - Apply the unknown sample to the column at an appropriate concentration using the same sample volume, fraction size, and flow rate as used for the Blue Dextran and the protein standards.
 - Determine the V_e of the unknown using the same methods applied to the standards.
 - Calculate the V_e/V_o for the unknown and determine its molecular mass from the standard curve.




Figure 1. Typical calibration curve obtained with proteins from the MWGF200 Kit run on Sephacryl® S-200-HR.

Precautions and Disclaimer

For Research Use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

References

1. Whitaker, J.R., *Anal. Chem.*, **35**(12), 1950-1953 (1963).
2. Andrews, P., *Biochem. J.*, **91**(2), 222-233 (1964).
3. Lee, Kelly, "Release of a Proposed Activator of Spreading Depolarization by Abrupt Hyperthermia". Queen's University (Kingston, Ontario), M.Sc. thesis, p. 34 (2016).
4. Martinez, Salette, "Investigating the Biochemical and Catalytic Properties of Nitrile Hydratases". Loyola University Chicago, Ph.D. dissertation, p. 88 (2014).
5. Burrows, Abigail S., "Detection of Food Proteins in Human Serum Using Mass Spectrometry Methods". University of Nebraska – Lincoln, Ph.D. dissertation, p. 215 (August 2020).
6. Male, Abigail, "Identification of Inhibitors of Protein-Protein Interactions Essential for Virulence of Pathogenic Bacteria". University of Southampton, Ph.D. dissertation, pp. 103, 118 (April 2014).
7. Bon Ramos, Adriana, "QueF and QueF-like: Diverse Chemistries in a Common Fold". Portland State University, Ph.D. dissertation, p. 33 (2016).
8. Gholami, Zahra, "PNA-protein conjugates for nano-scale modeling of protein aggregates". Nottingham Trent University, Ph.D. dissertation, p. 113 (November 2013).
9. Schwab, Ricarda, "Investigation of the interaction of FAT10 and VCP (p97)". Heinrich-Heine-Universität Konstanz, Dr. rer. nat. dissertation, p. 69 (2015).
10. Olivares, Philip, "Homeostasis of Primary and Secondary Metabolites". University of Illinois at Urbana-Champaign, Ph.D. dissertation, p. 58 (2019).
11. McMillan, Lana, "Oxidative Stress in Archaea: Proteome Responses and Characterization of an Inorganic Pyrophosphatase That Drives E1-Like Activation of Ubiquitin-Like Protein Modification". University of Florida, Ph.D. dissertation, p. 49 (2017).

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
MWGF200bul Rev 10/21 FF,MAM,KMR,GCY

MERCK

12. Magerl, Kathrin, "LOV domain signaling: A study of LOV-LOV interactions". Universität Regensburg, Dr. rer. nat. dissertation, pp. 107, 109 (2017).
13. Eichner, Ruth, "Functional characterization of Cereblon, the molecular target of immunomodulatory drugs". Technischen Universität München, Ph.D. dissertation, p. 85 (2016).
14. Bussmann, Michael, "cAMP dependent regulation of the TCA cycle by GlxR and oxidative stress response regulation by RosR in *Corynebacterium glutamicum*". Heinrich-Heine-Universität Düsseldorf, p. 69 (November 2009).
15. Gao, Song, "Structural and Functional Study of Human MxA Protein". Freien Universität Berlin, Dr. rer. nat. dissertation, pp. 40, 50 (February 2011).
16. Bjerregaard-Andersen, Kaare, "Structural and biophysical studies of the mammalian Na⁺-dependent Cl⁻-HCO₃⁻ exchanger NCBE and the bacterial enzyme isatin hydrolase". University of Oslo, Ph.D. dissertation, p. 30 (April 2014).
17. Piccoli, Ilaria, "Challenges of Conservation Agriculture on Silty Soils. Disentangling the Effects of Conservation Practices on Soil Organic Carbon Cycle and Soil Pore Network in Northeastern Italy". Università degli Studi di Padova, Ph.D. dissertation, p. 40 (2017).
18. Witmiss, Maria, "Untersuchungen zur Struktur und Funktion des Glycine Cleavage System in Cyanobakterien und Pflanzen" ("Investigations into the structure and function of the Glycine Cleavage System in cyanobacteria and plants"). Universität Rostock, Dr. rer. nat. dissertation, p. 23 (2019).
19. Meyer, Susanne Hilke, "Tth-IM60: eine Membraninsertase aus *Thermus thermophilus*" ("Tth-IM60: a membrane insert from *Thermus thermophilus*"). Universität Hohenheim, Dr. rer. nat. dissertation, p. 26 (2011).
20. Im, Ha Na, "Structural Studies of Csd6 Protein from *Helicobacter pylori* and Rv2258c Protein from *Mycobacterium tuberculosis*". Seoul National University, Ph.D. dissertation, p. 103 (2016).
21. Raasch, Katharina, "Der Oxoglutarat-Dehydrogenase-Komplex in *Corynebacterium glutamicum* und seine Interaktion mit OdhI" ("The oxoglutarate dehydrogenase complex in *Corynebacterium glutamicum* and its interaction with OdhI"). Forschungszentrum Jülich, Ph.D. dissertation, p. 33 (2012).
22. Marshall, J.J., *J. Chromatog.*, **53(2)**, 379-380 (1970).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
MWGF200bul Rev 10/21 FF,MAM,KMR,GCY

