1.08562.0050
1.04095.0250
1.08635.0100

REF

Microscopy

Aquatex®

(aqueous mounting agent) for microscopy

Glycerol

for fluorescence microscopy

Kaiser's glycerol gelatine, phenol-free

for microscopy

For professional use only

In Vitro Diagnostic Medical Device

Intended purpose

These ready-to-use aqueous mounting media are suited for mounting of water-containing sample material of human origin after these have been fixed and embedded as necessary, and then histologically, bacteriologically, hematologically (enzyme-cytochemically), or cytologically stained and, where applicable, counterstained with other *in vitro* diagnostic products from our portfolio, thus making them evaluable for further diagnostic procedures. Samples are mounted on slides to enable the specimen material to be examined by light microscopy, at the same time preserving it and thus enabling it to be re-examined many years later.

The appropriate aqueous mounting medium for the respective application is given in the corresponding instructions for use for our *In Vitro* Diagnostic staining solutions, solid dyes, and test kits.

Using the auxiliary reagents from our portfolio creates the conditions that enable authorized and qualified investigators to make a correct diagnosis at the end of the diagnostic process. In this regard, auxiliary IVD reagents serve *inter alia* to process human specimen material (e.g. fixing, decalcifying, dehydrating, clarifying, paraffin-embedding, mounting, microscoping, archiving). When used together with the corresponding staining solutions, this enables the visualization of cellular structures that are otherwise low in contrast, thus rendering them evaluable under the optical microscope. Further examinations may be necessary to reach a definitive diagnosis.

Principle

Mounting media are viscous, clear liquids with brilliant light-refraction properties. The aqueous mounting media in their dissolved form are then dropped onto the aqueous specimen of human origin, and the slide is covered air-tight with a cover glass. The evaporation of the water causes the mounting medium to harden, forming a solid, clear film under the cover glass, preserving the stained specimen material and thus enabling it to be kept for several years for re-analysis at a later date. As a result of the glass-similar refractive properties of the cover glass, the sample can now be observed under a microscope without any interference. Aqueous mounting media should be used especially with specimens prepared for enzyme and lipid determinations i.e. specimens that must not be dehydrated with an ascending series of alcohol.

Thanks to the practical user-friendly dropping bottle, the mounting medium can be easily and safely dropped onto the slide without smearing. The closure of the nozzle ensures that the viscosity of the medium remains constant, meaning that the mounting medium is immediately ready for use.

Sample material

Starting materials are

- fixed and stained cytological smears, e.g. sputum, fine needle aspiration biopsies (FNAB), rinses, imprints, effusions
- formalin-fixed, paraffin-embedded, histologically stained tissue specimens (3 - 5 µm thick paraffin sections) or stained cryosections
- air-dried, heat-fixed, and stained smears of bacteriological specimen material, e.g. liquid and solid enrichment cultures of bacteria from body fluids, exudates, pus
- enzyme cytochemically processed and stained blood or bone-marrow smears.

Reagents

Cat. No. 1.08562	Aquatex® (aqueous mounting agent) for microscopy	50-ml dropping bottle
Cat. No. 1.04095	Glycerol for fluorescence microscopy	250 ml
Cat. No. 1.08635	Kaiser's glycerol gelatine, phenol-free for microscopy	100-g dropping bottle

Specifications

Cat. No. 1.08562 - Aquatex® (aqueous mounting medium) for microscopy

is an aqueous mounting medium used primarily in connection with hematology specimens stained by means of enzyme-cytochemical reactions. It can also be used for mounting purposes in connection with immunohistological specimens or histological lipid stains in frozen sections, e.g. using Oil red O staining solution (Cat. No. 1.02419).

Refractive index (20 °C)	1.390 - 1.400
Viscosity (20 °C)	170 - 280 mPa*s
Fluorescence (at 340 nm as quinine)	≤ 500 ppb

Cat. No. 1.04095 - Glycerol for fluorescence microscopy

is used as a mounting medium for stained histology specimens (e.g. toluidine blue stain). It is primarily designed to preserve the fluorescence of durable specimens for future diagnosis.

Refractive index (20 °C)	1.472 - 1.474
Fluorescence (at 365 nm as quinine)	≤ 100 ppb
Fluorescence (at 313 nm as quinine)	≤ 100 ppb
Fluorescence (at 254 nm as quinine)	≤ 100 ppb
Purity (20°C)	≥ 99.5%
Density	1.26 g/cm³

Cat. No. 1.08635 - Kaiser's glycerol gelatine, phenol-free for microscopy

is an aqueous, phenol-free mounting medium consisting of gelatine and glycerol, and is used primarily in the mounting of stained histology specimens (e.g. Nile blue staining). The pH of the product lies in the neutral range. Prior to use, the mounting medium is melted in the heating cabinet at approx. 40 °C and is dropped onto the specimen from the practical dropping bottle.

Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

The specimen material is processed, stained (and counterstained where applicable), and mounted according to the instructions for use of our *In Vitro* Diagnostic staining solutions, solid dyes, and test kits.

Reagent preparation

Aquatex® and **Glycerol** are ready-to-use.

Kaiser's glycerol gelatine, phenol-free is melted at 40 °C in a heating cabinet.

A dilution of the listed mounting media is not necessary.

Procedure

The mounting medium is applied to the horizontal slide, using a glass rod or else directly dropping approx. 0.2 ml of one of the listed mounting media from the dropping bottle. As soon as a homogeneous distribution of the solution is guaranteed, gently add a clean cover glass, so that the space between the slide and cover glass is filled without air bubbles with mounting medium. Allow this setup to dry and harden overnight in a horizontal position.

If specimens are to be kept for a longer period of time, they must be covered, air-tight, with laquers (i.e. by colorless nail polish) or liquids that resinify or polymerise.

Note

Hematoxylin stains are unstable once mounted with **Glycerol** and **Kaiser's glycerol gelatine, phenol-free**.

The nuclei stains achieved with hematoxylin are swiftly destroyed.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Result

The use of these aqueous, ready-to-use mounting media results in completely airtight specimen slides, the structure and stain pattern of which remains preserved over the long term, enabling them to be microscopically re-analyzed at a later date.

Trouble-shooting

No color stability over longer storage times

- Care must be taken to maintain a minimum quality and dye concentration of the staining solutions as a measure to stabilize the stain of the specimen.

Air bubbles and inclusions

- The volume of the mounting medium applied to the specimen must be carefully monitored (to avoid too much or too little mounting medium).
- The drying times for the specimens must be observed. The specimens must be completely dehydrated before microscopy with immersion oil, i. e. always allow the specimens to dry completely and mount thoroughly.
- The evaporation of the solvent after mounting must be borne in mind, and specimen slides must be dried for at least overnight.

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using automated-mounting instruments, please follow the instructions for use supplied by the supplier of the system and software. Remove surplus immersion oil before filing.

Analytical performance characteristics

The present auxiliary reagents "Glycerol", "Aquatex®", and "Kaiser's glycerol gelatine, phenol-free" aid in the microscopic examination of biological structures as described in the "Intended purpose" of this IFU. The use of the products is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, histoprocessing (Cat. Nos. 108562 and 108635), decisions regarding suitable controls and more.

The analytical performance of the products is confirmed by testing each production batch.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100%:

Cat. No. 1.04095 - Glycerol

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Histological and physical methods				
Suitability for microscopy	7/7	7/7	6/6	6/6
Refractive index (n_{D}^{20})	7/7	7/7	N/A	N/A

Analytical performance results

Cat. No. 1.08562 - Aquatex®

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological and physical methods				
Suitability for microscopy	4/4	4/4	6/6	6/6
Refractive index (n_{D}^{20})	4/4	4/4	6/6	6/6

Analytical performance results

Cat. No. 1.08635 - Kaiser's glycerol gelatine, phenol-free

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Histological and physical methods				
Suitability for microscopy	3/3	3/3	9/9	9/9

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

The results of this Performance Evaluation confirms that the products are suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used.

These products are auxiliary reagents these, when used together with other IVD products such as staining solutions, render human specimen material evaluable for diagnostic purposes.

Further tests must be selected and implemented according to recognized methods.

Suitable controls should be conducted with each application in order to avoid an incorrect result.

Storage

Cat. Nos. 1.08562, 1.08635:

Store the listed mounting media at +15 °C to +25 °C.

Cat. No. 1.04095:

Store the Glycerol for fluorescence microscopy at +2 °C to +30 °C.

Shelf-life

Cat. Nos. 1.08562, 1.08635:

The listed mounting media can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

Cat. No. 1.04095:

The Glycerol for fluorescence microscopy can be used until the stated expiry date.

After opening the bottle the first time and subsequent storage of the tightly reclosed bottle at +2 °C to +30 °C, the medium can be used up to the printed expiry date.

Additional instructions

For professional use only.

In order to avoid errors, the application must be carried out by qualified personnel only.

National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.

Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No. 1.04699 Immersion oil for microscopy

100-ml dropping bottle,
100 ml, 500 ml

Hazard classification

Cat. Nos. 1.08562, 1.04095, 1.08635

Please observe the hazard classification printed on the label and the information given in the safety data sheet.

The safety data sheet is available on the website and on request.

Main components of the products

Cat. No. 1.08562

CAS-No 9003-39-8 33.22% (w/w)
(C₆H₅NO)_n 50%
 n_{D}^{20} ~1.4
pH ~7.0
1 l = 1.08 kg

Cat. No. 1.04095

C₃H₈O₃
M = 92.10 g/mol
1 l = 1.26 kg

Cat. No. 1.08635

C₃H₈O₃ 48.8% (w/w)
CAS-No 9000-70-8 7.3% (w/w)

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and / or its authorised representative and to your national authority.

Literature

- Romeis - Mikroskopische Technik, Editors: Maria Mulisch, Ulrich Welsch, 2015, Springer Spektrum, 19. Auflage
- Histological and Histochemical Methods, Theory and practice, J.A. Kiernan, 2015, Scion Publishing Ltd, 5th Edition

3. Laboratory Manual of Histochemistry, Linda L. Vacca, 1985, Raven Press
4. Theory and Practice of Histological Techniques, John D Bancroft, Marilyn Gamble, 2008, Churchill Livingstone ELSEVIER, sixth Edition

Cat. No. 1.08635

- H317: May cause an allergic skin reaction.
- P261: Avoid breathing dust.
- P272: Contaminated work clothing should not be allowed out of the workplace.
- P280: Wear protective gloves.
- P302 + P352: IF ON SKIN: Wash with plenty of water.
- P333 + P313: If skin irritation or rash occurs: Get medical advice/ attention.
- P362 + P364: Take off contaminated clothing and wash it before reuse.

Revision History

Version	Modification Comment
2024-Jul-01	Initial version with the introduction of Revision History

Consult instructions for use

Manufacturer

Catalog number

Batch code

Caution, consult accompanying documents

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© 2024 Merck KGaA, Darmstadt, Alemania y/o sus filiales. Todos los derechos reservados. Merck y Sigma-Aldrich son marcas comerciales de Merck KGaA, Darmstadt, Alemania. Todas las demás marcas comerciales son propiedad de sus respectivos propietarios. Tiene a su disposición información detallada sobre las marcas comerciales a través de recursos accesibles al público.
2. Histological and Histochemical Methods, Theory and practice, J. A. Kiernan, 2015, Scion Publishing Ltd, 5th Edition