1.02537.0001

Spectroquant® Sulfate Test

SO₄²⁻

1. Method

Sulfate ions react with barium ions to form slightly soluble barium sulfate. The resulting turbidity is measured in the photometer (turbidimetric method).

The method is analogous to EPA 375.4, APHA 4500-SO₄²⁻ E, and ASTM D516-16.

2. Measuring range and number of determinations

Measuring range	Number of determinations
5 - 300 mg/l SO ₄ ²⁻	100

For programming data for selected photometers / spectrophotometers see www.sigmaaldrich.com/photometry.

3. Applications

Sample material:

Groundwater and surface water, seawater
Drinking water and mineral water
Wastewater
Nutrient solutions for fertilization
Soils after appropriate sample pretreatment
Water for concrete in the construction industry

4. Influence of foreign substances

This was checked individually in solutions containing 150 and 0 mg/l SO₄²⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %			
Ag ⁺	2	Mn ²⁺	500
Al ³⁺	1000	NH ₄ ⁺	500
Ca ²⁺	1000	Ni ²⁺	100
Cd ²⁺	1000	NO ₂ ⁻	1000
CN ⁻	250	NO ₃ ⁻	100
CO ₃ ²⁻	250	Pb ²⁺	10
Cr ³⁺	20	PO ₄ ³⁻	100
Cr ₂ O ₇ ²⁻	5	S ²⁻	2
Cu ²⁺	1000	SiO ₃ ²⁻	1000
F ⁻	250	SO ₃ ²⁻	1
Fe ³⁺	50	S ₂ O ₃ ²⁻	0.5
Hg ²⁺	500	Zn ²⁺	500
Mg ²⁺	1000		
		EDTA	100
		Hydrazine	250
		Anionic surfactants¹⁾	0.5
		Cationic surfactants ²⁾	1000
		Nonionic surfactants ³⁾	1000
		Na-acetate	10%
		NaCl	10%

¹⁾ tested with Na-dodecyl sulfate

²⁾ tested with N-cetyl-N,N,N-trimethylammonium bromide

³⁾ tested with Triton® X-100

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent SO₄-1
1 bottle of reagent SO₄-2
1 AutoSelector

Other reagents and accessories:

MQuant® Sulfate Test, Cat. No. 1.10019,
measuring range <200 - >1600 mg/l SO₄²⁻
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sodium hydroxide solution 1 mol/l Titripur®, Cat. No. 1.09137
Hydrochloric acid 1 mol/l Titripur®, Cat. No. 1.09057
Membrane filters 0.45 µm
Spectroquant® CombiCheck 10, Cat. No. 1.14676
Sulfate standard solution, 40.0 mg/l SO₄²⁻, Cat. No. 1.25050
Sulfate standard solution, 125 mg/l SO₄²⁻, Cat. No. 1.25051

Pipettes for pipetting volumes of 0.50 and 5.0 ml
Rectangular cells 10 mm (2 pcs), Cat. No. 1.14946

6. Preparation

- Analyze immediately after sampling.
- The pH must be within the range 2 - 10.**
Adjust, if necessary, with sodium hydroxide solution or hydrochloric acid.
- Filter turbid samples through a 0.45 µm membrane filter.

7. Procedure

Reagent SO ₄ -1	0.50 ml	Pipette into a test tube.
Pretreated sample (20 - 40 °C)	5.0 ml	Add with pipette and mix.
Reagent SO ₄ -2	1 level blue microspoon (in the cap of the SO ₄ -2 bottle)	Add and shake vigorously until the reagent is completely dissolved.
Leave to stand for exactly 2 min (reaction time), then fill the measurement sample into the cell, and measure in the photometer.		

Notes on the measurement:

- Certain photometers may require a blank** (preparation as per measurement sample, but with distilled water instead of sample).
- For turbidimetric measurement the cells must be clean.
Wipe, if necessary, with a clean dry cloth.
- The pH of the measurement solution must be within the range 1.0 - 2.0.
- The turbidity of the measurement solution remains stable for only 2 min max. after the end of the reaction time stated above.**

8. Analytical quality assurance

recommended before each measurement series
To check the photometric measurement system (test reagent, measurement device, handling) and the mode of working, the sulfate standard solutions (see section 5) or Spectroquant® CombiCheck 10 can be used. Besides a **standard solution** with 100 mg/l SO₄²⁻, CombiCheck 10 also contains an **addition solution** for determining sample-dependent interferences (**matrix effects**).
Additional notes see under www.qa-test-kits.com.
For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Notes

- Reclose the reagent bottles immediately after use.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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