1.01812.0001

Spectroquant® Sulfate Test

1. Method

Sulfate ions react with barium ions to form slightly soluble barium sulfate. The resulting turbidity is measured in the photometer (turbidimetric method).

The method is analogous to EPA 375.4, APHA 4500-SO₄²⁻ E, and ASTM D516-16.

2. Measuring range and number of determinations

Cell mm	Measuring range mg/l SO ₄ ²⁻	Number of determinations
50	0.50 - 10.00	100
20	1.3 - 25.0	
10	2.5 - 50.0	

For programming data for selected photometers / spectrophotometers see www.sigmaaldrich.com/photometry.

3. Applications

Sample material:

Groundwater and surface water
Drinking water and mineral water
Wastewater
Nutrient solutions for fertilization
Soils after appropriate sample pretreatment
Water for concrete in the construction industry
This test is **not suited** for seawater.

4. Influence of foreign substances

This was checked individually in solutions containing 5 and 0 mg/l SO₄²⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag ⁺	1	Mn ²⁺	1000	EDTA	50
Al ³⁺	1000	NH ₄ ⁺	1000	Hydrazine	1000
Ca ²⁺	1000	Ni ²⁺	500	Anionic surfactants ¹⁾	0.5
Cd ²⁺	1000	NO ₂ ⁻	1000	Cationic surfactants ²⁾	250
CN ⁻	1000	NO ₃ ⁻	50	Nonionic surfactants ³⁾	250
CO ₃ ²⁻	1000	Pb ²⁺	2	Na -acetate	10%
Cr ³⁺	5	PO ₄ ³⁻	100	Na Cl	2%
Cr ₂ O ₇ ²⁻	2	S ²⁻	5		
Cu ²⁺	1000	SiO ₃ ²⁻	1000		
F ⁻	100	SO ₃ ²⁻	5		
Fe ³⁺	50	S ₂ O ₃ ²⁻	2		
Hg ²⁺	500	Zn ²⁺	1000		
Mg ²⁺	1000				

¹⁾ tested with Na-dodecyl sulfate

²⁾ tested with N-cetyl-N,N,N-trimethylammonium bromide

³⁾ tested with Triton® X-100

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent SO₄-1
1 bottle of reagent SO₄-2
1 AutoSelector

Other reagents and accessories:

MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sodium hydroxide solution 1 mol/l Titripur®, Cat. No. 1.09137
Hydrochloric acid 1 mol/l Titripur®, Cat. No. 1.09057
Membrane filters 0.45 µm
Sulfate standard solution Certipur®, 1000 mg/l SO₄²⁻, Cat. No. 1.19813

Pipettes for pipetting volumes of 0.50 and 10 ml
Rectangular cells 10, 20, and 50 mm (2 of each), Cat. Nos. 1.14946, 1.14947, and 1.14944

6. Preparation

- Analyze immediately after sampling.
- The pH must be within the range 2 - 10.**
Adjust, if necessary, with sodium hydroxide solution or hydrochloric acid.
- Filter turbid samples through a 0.45 µm membrane filter.

7. Procedure

Reagent SO ₄ -1	0.50 ml	Pipette into a test tube.
Pretreated sample (20 - 40 °C)	10 ml	Add with pipette and mix.
Reagent SO ₄ -2	1 level green microspoon (in the cap of the SO ₄ -2 bottle)	Add and shake vigorously until the reagent is completely dissolved.

Leave to stand for **exactly 2 min (reaction time)**, then fill the measurement sample into the cell, and measure in the photometer.

Notes on the measurement:

- Certain photometers may require a blank** (preparation as per measurement sample, but with distilled water instead of sample).
- When using the 50-mm cell**, perform the measurement against a separately prepared blank (preparation as per measurement sample, but with distilled water instead of sample).
Configure the photometer for blank measurement.
- For turbidimetric measurement the cells must be clean.
Wipe, if necessary, with a clean dry cloth.
- The pH of the measurement solution must be within the range 1.0 - 2.0.
- The turbidity of the measurement solution remains stable for only 2 min max. after the end of the reaction time stated above.**

8. Analytical quality assurance

recommended before each measurement series
To check the photometric measurement system (test reagents, measurement device, handling) and the mode of working, a dilute sulfate standard solution containing 5.00 mg/l SO₄²⁻ can be used.

Sample-dependent interferences (matrix effects) can be determined by means of standard addition.

Additional notes see under www.qa-test-kits.com.

For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Notes

- Reclose the reagent bottles immediately after use.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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