1.10537.0001

Carrez clarification

Cat.No. 1.10537.0001

Carrez clarification (5-fold concentrate)
Reagent kit for sample preparation in food analysis

Contents of package, storage, and stability

Carrez solution I: contains potassium hexacyanoferrate(II) trihydrate ($K_4[Fe(CN)_6] \times 3H_2O$)
 100 ml concentration: 15 g/100 ml

Carrez solution II: contains zinc sulfate heptahydrate ($ZnSO_4 \times 7H_2O$)
 100 ml concentration: 30 g/100 ml

The reagents may be used for at least 1 year if stored in a refrigerator at 4°C.

Carrez clarification is often used within the sample preparation in food analysis. The ready-to-use reagents of **Carrez clarification** may be used without further dilution. Depending on the volume and characteristics of the sample solution the reagents will be enough for about 100 preparations.

Carrez clarification is suitable for the **precipitation of proteins, the elimination of turbidity, and for breaking of emulsions**, which may interfere with subsequent determination steps.

On the addition of Carrez solution I and II to an aqueous sample a sparingly soluble precipitate consisting of $Zn_2[Fe(CN)_6]$ is formed to which high-molecular-weight compounds will be bound by adsorption and be precipitated, too. This sediment can be separated from the cleared solution by filtration or centrifugation.

Carrez clarification is used predominantly prior to the analysis of carbohydrates e. g. glucose, fructose, sucrose, lactose, lactulose, maltose etc. but also as sample preparation step during the analysis of some other important food ingredients like ethanol, glycerol, sorbitol, acetaldehyde, D-/L-lactate, nitrate.

But, the method is **unsuitable** for sample preparation prior to the determination of ascorbic acid (vitamin C), citrate and urea/ammonia.

Instructions for use

To a liquid sample (approx. 1 - 5 g in approx. 60 ml water), possibly pretreated (e. g. by homogenization or extraction with hot water) in an Erlenmeyer flask, **1 ml of Carrez solution I** is added and mixed carefully. Subsequently, **1 ml of Carrez solution II** is added and mixed again. If no precipitate appears, the clarification has to be repeated in a second sample using a five-fold volume of Carrez solutions (5 ml each).

Afterwards the pH of the solution is adjusted to 7.5 - 8.0 by addition of sodium hydroxide solution 1 mol/l (Cat. No. 1.09137) in order to precipitate excess zinc ions as $Zn(OH)_2$. The contents of the Erlenmeyer flask are transferred quantitatively to a 100 ml volumetric flask, and the flask filled to the mark with distilled water. Subsequently the solution is filtered through a fluted filter or centrifuged to separate the insoluble compounds. The clear filtrate is used for the following determination.

Avoiding interferences

In order to avoid any interferences it is very important to adjust the pH value to 7.5 - 8.0 after addition of the Carrez solutions. In specific cases we recommend to add some EDTA (Titriplex® III, EDTA- Na_2) (Cat. No. 1.08418) after the adjustment of the pH up to 1 mmol/l sample solution (approx. 0.37 mg/ml) making sure the excess zinc ions are complexated.

In case of doubt a possible adsorption of analyte to the precipitate should be checked by recovery tests.

Order informations – reagents

Cat. No.	Designation	Package size
1.10537	Carrez clarification Reagent kit for sample preparation in food analysis (5-fold concentrate)	1 pack
1.09137	Sodium hydroxide solution 1 mol/l (1N) Titripur®	1 l
1.08418	Titriplex® III GR for analyses ACS, ISO	100 g, 250 g, 1 kg

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MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, Supelco, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

EMD Millipore Corporation, 400 Summit Drive, Burlington MA 01803, USA, Tel. +1-978-715-4321
 MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville, Ontario, L6H 6J8, Canada, Phone: +1 800-565-1400
www.sigmaaldrich.com

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