Instructions for Use

Lymphocyte α -Naphthyl Butyrate Esterase Kit

Procedure No. 181


Intended Use

Lymphocyte α -Naphthyl Butyrate Esterase Kit is intended for use as a general laboratory histochemical stain for the detection of α -naphthyl butyrate esterase in blood, bone marrow films and tissue touch preparations. Lymphocyte α -Naphthyl Butyrate Esterase Kit is for "In Vitro Diagnostic Use." For professional use only. The data obtained from this manual, qualitative procedure serves only as an aid to diagnosis of human specimens and should be reviewed in conjunction with other clinical diagnostic tests or information.

The reaction listed below for α -Naphthyl Butyrate Esterase will also demonstrate nonspecific esterase activity in monocytes and macrophages. Several different methods are available for this purpose. The method listed below is not specific for monocytes and macrophages since the pH of the incubating media has been adjusted to demonstrate a focal (dot) staining pattern in some lymphocytes. The α -Naphthyl Butyrate Esterase is less sensitive than the α -Naphthyl Acetate Esterase procedure. For procedures more sensitive to monocytes Sigma offers two Esterase kits, Catalog Nos. 90A1 and 91A. These kits also include instructions for the demonstration of α -naphthyl acetate esterase with fluoride inhibition.

According to Sigma-Aldrich techniques, cytocentrifuge preparations or films are fixed in a citrate-acetone-formaldehyde solution. α -Naphthyl butyrate esterase is then visualized by the following simultaneous capture principles.

α -Naphthyl Butyrate Esterase

Reagents

α -Naphthyl Butyrate Solution (Cat. No. 1801-50ML)

α -Naphthyl Butyrate, 2.4 g/L, in methanol solution with solubilizers

Pararosaniline Solution (Cat. No. 1804-15ML)

Pararosaniline, 40 g/L, C.I. 42500, in 2 mol/L hydrochloric acid. Danger. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. May cause cancer. Avoid breathing dust/fume/gas/mist/vapours/spray.

Sodium Nitrite Tablets (Cat. No. 1809-10TAB)

Sodium nitrite, 250-260 mg per tablet. Danger. May intensify fire; oxidizer. Toxic if swallowed. Causes serious eye irritation. Very toxic to aquatic life. Keep/store away from clothing/combustible materials. Avoid release to the environment. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.

Citrate Solution (Cat. No. 915-50ML)

Citric acid, 18 mmol/L, sodium citrate, 9 mmol/L, sodium chloride, 12 mmol/L, and surfactant. The pH should be 3.6 ± 0.1 .

Phosphate Buffer (Cat. No. 1805-1VL)

Sodium phosphate, 59.6 mM, potassium phosphate, 7.1 mM.

Methylene Blue Solution (Cat. No. 1808-50ML)

Methylene blue, C.I. 52015, 1.4% (w/v) in 95% ethanol

Special Materials Required but Not Provided

- Acetone, ACS Reagent
- Formaldehyde, 37%
- Water bath capable of maintaining 37°C shielded from light

Storage and Stability

Store α -Naphthyl Butyrate Solution in freezer below 0°C. Warm solution to 37°C and mix well prior to use. Discard if reagents turn yellow or if precipitate forms.

Store Pararosaniline Solution, Phosphate Buffer, Sodium Nitrite Tablets and Methylene Blue Solution at room temperature (18–26°C) protected from light.

Store Sodium Nitrite Tablet Solution and Phosphate Buffer Solution refrigerated (2–8°C).

Reagent labels bear expiration date.

Citrate Solution is suitable for use in the absence of microbial growth.

Deterioration

Discard α -Naphthyl Butyrate Solution if reagent turn bright yellow or if precipitate forms.

Discard Pararosaniline Solution if solution does not turn amber upon addition of Sodium Nitrite Tablets Solution.

Phosphate Buffer Solution should be discarded if turbidity develops.

Preparation

Warm α -Naphthyl Butyrate Solution to 37°C. Sodium Nitrite Tablet Solution, 4 g/dL is prepared by dissolving 1 tablet in 6.25 mL deionized water. Store this solution tightly stoppered at 2–8°C. Warm to room temperature before use. Discard if turbidity develops.

Phosphate Buffer Solution is prepared by dissolving contents of Phosphate Buffer in 500 mL deionized water. Buffer has concentration of 0.067 mol/L, pH 7.7 at 25°C. Store in refrigerator (2–8°C). Discard if microbial growth is evident.

Note

Microbial growth may be retarded by filtration through a 0.22 micron Filter Unit.

Methylene Blue Counterstain is prepared by adding 5 mL Methylene Blue Solution to 45 mL deionized water. Mix well. Prepare fresh daily.

Citrate-Acetone-Formaldehyde Fixative: To 25 mL Citrate Solution add 65 mL Acetone and 8 mL of 37% Formaldehyde. Place in glass bottle and cap tightly. Store in refrigerator (2–8°C). Warm to 18–26°C prior to use. Stable if stored tightly capped in refrigerator.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Samples may be collected in either EDTA or heparin. After fixation slides may be stored at room temperature for at least 2 weeks. If mononuclear cells are to be isolated using HISTOPAQUE®-1077, separation should be performed within 4 hours although fair recovery has been noted after 24 hours. Blood films, particularly from leukopenic individuals, are not recommended for these procedures since evaluation of this material is quite time consuming. Cytocentrifuge or buffy coat preparations should be employed. Bone marrow films and tissue touch preparations pose no problems with respect to microscopic evaluation.

Notes

Cells from healthy donors should be included with each test.

Cells isolated on polysucrose-sodium diatrizoate gradients may be stored in liquid nitrogen (LN₂) for control purposes. To accomplish this, $1 \times 10^7 - 10^8$ cells are frozen at $1 \pm 0.3^\circ\text{C}$ per minute in a medium containing 50% fetal calf serum, 40% RPMI-1640 (or other appropriate tissue culture fluids) and 10% dimethylsulfoxide (DMSO). They may be stored in either the liquid or vapor phase of LN₂.

Slides should be evaluated using 1000X magnification. Didymium filters may enhance color particularly of the diffuse staining α -NB.

Procedure

1. Prewarm Phosphate Buffer Solution to 37°C.
2. Immediately prior to fixation, add 1.5 mL of Sodium Nitrite Tablet Solution to 1.5 mL Pararosaniline Solution. Mix gently by inversion and let stand at least 5 minutes, then add to 40 mL of prewarmed Phosphate Buffer Solution.
3. Add 5 mL α -Naphthyl Butyrate Solution.
4. Mix well and pour into Coplin jar. Solution will be amber. Formation of precipitate indicates reagent deterioration.
5. Fix slides for 10 seconds in Citrate-Acetone-Formaldehyde Fixative at room temperature (18–26°C). Rinse in deionized water for 45 seconds. Do not allow slides to dry.
6. Immediately after rinsing, place slides into solution from Step 4 and incubate 1 hour at 37°C. **Note:** If slides are not placed in incubation solution after fixation, they must be allowed to air dry for at least 45 minutes.
7. After 1 hour remove slides from Coplin jar and rinse 2–3 minutes in running tap water. Discard staining solution.
8. Allow slides to air dry at least 15 minutes before counterstaining.
9. Counterstain 5 minutes in Methylene Blue Counterstain.
10. Rinse in deionized water.

Performance Characteristics

The focal (dot) staining pattern observed in some lymphocytes is associated with mature T-cells bearing Fc receptors for IgM. This subset overlaps the T-helper cell population to some extent but is not an accurate measure of helper function. The diffuse or granular staining lymphocytes bear receptors to the Fc portion of IgG and overlap to some extent the T suppressor cell population.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
1801	α -Naphthyl Butyrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
1804	Pararosaniline Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
1809	Sodium Nitrite Tablets	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
915	Citrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
1805	Phosphate Buffer	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
1808	Methylene Blue Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

- H226 Flammable liquid and vapor.
- H272 May intensify fire, oxidizer.
- H290 May be corrosive to metals.
- H301 Toxic if swallowed.
- H312 + H332 Harmful in contact with skin or if inhaled.
- H315 Causes skin irritation.
- H319 Causes serious eye irritation.
- H335 May cause respiratory irritation.
- H350 May cause cancer.
- H370 Causes damage to organs (Eyes).
- H400 Very toxic to aquatic life.
- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P273 Avoid release to the environment.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.
- P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
- P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Zito DR, Johnson KL, Skipper ME, Normansell DE: Isolation of mononuclear cells from the whole blood: Comparison of two media. Lab Med 17:94, 1986

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 3.0	2014
Rev. 4.0	2016
Rev. 5.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards. Removed references to 181A kit. Added CH-REP information.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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MDSS GmbH