1.00798.0001
1.00798.0007

Spectroquant® Phosphate Test

P

for the determination of orthophosphate
USEPA equivalent for drinking water and wastewater

1. Method

In sulfuric solution orthophosphate ions react with molybdate ions to form molybdophosphoric acid. Ascorbic acid reduces this to phosphomolybdenum blue (PMB) that is determined photometrically.
The method is analogous to EPA 365.2+3, APHA 4500-P E, and DIN EN ISO 6878.

2. Measuring range and number of determinations

Cell mm	Measuring range			Number of determinations
	mg/l PO ₄ -P	mg/l PO ₄ ³⁻	mg/l P ₂ O ₅	
10	1.0 - 100.0	3 - 307	2 - 229	100

For programming data for selected photometers / spectrophotometers see
www.sigmaaldrich.com/photometry.

3. Applications

This test measures only orthophosphate.

Sample material:

Groundwater and surface water, seawater
Drinking water
Wastewater
Nutrient solutions for fertilization
Soils after appropriate sample pretreatment
Food after appropriate sample pretreatment

4. Influence of foreign substances

This was checked individually in solutions containing 50 and 0 mg/l PO₄-P. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag ⁺	1000	Cu ²⁺	500	Ni ²⁺	1000
AsO ₄ ³⁻	10	F ⁻	1000	NO ₂ ⁻	1000
Ca ²⁺	1000	Fe ³⁺	1000	Pb ²⁺	500
Cd ²⁺	1000	Hg ²⁺	100	SiO ₃ ²⁻	1000
CN ⁻	1000	Mg ²⁺	1000	Zn ²⁺	1000
Cr ³⁺	1000	Mn ²⁺	1000	EDTA	1000
Cr ₂ O ₇ ²⁻	500	NH ₄ ⁺	1000	Surfactants ¹⁾	100
				Na-Acetate	10%
				NaCl	15%
				NaNO ₃	20%
				Na ₂ SO ₄	10%

Reducing agents interfere with the determination.

¹⁾ tested with nonionic, cationic, and anionic surfactants

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!
The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent PO₄-1
1 bottle of reagent PO₄-2
2 blue dose-metering caps (each can be used for 50 dosages)
1 AutoSelector

Other reagents and accessories:

MQuant® Phosphate Test, Cat. No. 1.10428,
measuring range 10 - 500 mg/l PO₄³⁻ (3.3 - 163 mg/l PO₄-P)
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
Water for analysis EMSURE®, Cat. No. 1.16754
Phosphate standard solution Certipur®, 1000 mg/l PO₄³⁻, Cat. No. 1.19898
Hydrochloric acid 25% for analysis EMSURE®, Cat. No. 1.00316
Sodium hydroxide solution 1 mol/l Titripur® (approx. 4%),
Cat. No. 1.09137

Pipettes for pipetting volumes of 0.50 and 8.0 ml
Rectangular cells 10 mm (2 pcs), Cat. No. 1.14946

6. Preparation

- Use only phosphate-free detergents to rinse glassware. Otherwise fill with hydrochloric acid (approx. 10%) and leave to stand for several hours.

At the first use **replace the screw cap of the reagent bottle PO₄-2 by the blue dose-metering cap.**

Hold the reagent bottle **vertically** and, at each dosage, press the slide **all the way** into the dose-metering cap. **Before each dosage** ensure that the slide is **completely retracted**.

Reclose the reagent bottle with the screw cap at the end of the measurement series, since the function of the reagent is impaired by the absorption of atmospheric moisture.

- Analyze immediately after sampling.
- Check the phosphate content with the MQuant® Phosphate Test. Samples containing more than 100.0 mg/l PO₄-P (307 mg/l PO₄³⁻) must be diluted with distilled water.
- The pH must be within the range 0 - 10.**
Adjust, if necessary, with sulfuric acid.
- Filter turbid samples.

7. Procedure

Distilled water ¹⁾ (10 - 35 °C)	8.0 ml	Pipette into a test tube.
Pretreated sample (10 - 35 °C)	0.50 ml	Add with pipette and mix.
Reagent PO ₄ -1	0.50 ml	Add with pipette and mix.
Reagent PO ₄ -2	1 dose	Add and shake vigorously until the reagent is completely dissolved .

Leave to stand for 5 min (reaction time), then fill the sample into a 10-mm cell, and measure in the photometer.

¹⁾ It is recommended to use water for analysis EMSURE®, Cat. No. 1.16754.

Notes on the measurement:

- Certain photometers may require a blank** (preparation as per measurement sample, but with distilled water instead of sample). The blank is slightly yellow.
- For photometric measurement the cells must be clean. Wipe, if necessary, with a clean dry cloth.
- Measurement of turbid solutions yields false-high readings.
- The pH of the measurement solution must be within the range 0.80 - 0.95.
- The color of the measurement solution remains stable for at least 60 min after the end of the reaction time stated above.

8. Analytical quality assurance

recommended before each measurement series

To check the photometric measurement system (test reagents, measurement device, handling) and the mode of working, a dilute phosphate standard solution containing 50.0 mg/l PO₄-P (153.3 mg/l PO₄³⁻) can be used.

Sample-dependent interferences (matrix effects) can be determined by means of standard addition.

Additional notes see under www.qa-test-kits.com.

For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Notes

- Reclose the reagent bottles immediately after use.
- Cleanse the test tubes and the cells from time to time as follows:
Fill with sodium hydroxide solution (approx. 0.4%) and leave to stand for max. 1 hour.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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© 2024 Merck KGaA, Darmstadt, Allemagne et/ou ses sociétés affiliées. Tous droits réservés.