1.14942.0001

Spectroquant® Nitrate Test in seawater

NO₃⁻

USEPA equivalent for wastewater

1. Method

In strongly sulfuric solution, in the presence of chloride nitrate ions react with resorcinol to form a red-violet indophenol dye that is determined photometrically.

2. Measuring range and number of determinations

Cell mm	Measuring range		Number of determinations
	mg/l NO ₃ -N	mg/l NO ₃ ⁻	
10	0.2 - 17.0	0.9 - 75.3	50

For programming data for selected photometers / spectrophotometers see www.sigmaaldrich.com/photometry.

3. Applications

This test enables the measurement of very low nitrate contents **even in the presence of large quantities of chloride**.

Sample material:

Groundwater and water of precipitation
Surface water including seawater and brackish water
Wastewater and percolating water

4. Influence of foreign substances

This was checked individually in solutions containing 12 and 0 mg/l NO₃-N. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag ⁺	50	Cu ²⁺	100	Ni ²⁺	1000
Al ³⁺	1000	F ⁻	1000	NO ₂ ⁻	2.5 ¹⁾
Ca ²⁺	1000	Fe ³⁺	10	Pb ²⁺	250
Cd ²⁺	1000	Hg ²⁺	250	PO ₄ ³⁻	1000
CN ⁻	1000	Mg ²⁺	1000	SiO ₃ ²⁻	100
Cr ³⁺	100	Mn ²⁺	1000	Zn ²⁺	500
Cr ₂ O ₇ ²⁻	10	NH ₄ ⁺	500	Surfactants ²⁾	500
				COD	500
				(glucose)	500
				NaCl	20%
				Na ₂ SO ₄	20%

¹⁾ In cases of higher concentrations, eliminate nitrite ions acc. to section 6.

²⁾ tested with nonionic, cationic, and anionic surfactants

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent NO₃-1
1 bottle of reagent NO₃-2
1 bottle of reagent NO₃-3
1 AutoSelector

Other reagents and accessories:

MQuant® Nitrite Test, Cat. No. 1.10007,
measuring range 2 - 80 mg/l NO₂⁻ (0.6 - 24 mg/l NO₂-N)
Amidosulfuric acid for analysis EMSURE®, Cat. No. 1.00103
MQuant® pH-indicator strips pH 0 - 6.0, Cat. No. 1.09531
Sulfuric acid 25% for analysis EMSURE®, Cat. No. 1.00716
Spectroquant® CombiCheck 20, Cat. No. 1.14675
Nitrate standard solution, 0.500 mg/l NO₃-N, Cat. No. 1.25036
Nitrate standard solution, 2.50 mg/l NO₃-N, Cat. No. 1.25037
Nitrate standard solution, 15.0 mg/l NO₃-N, Cat. No. 1.25038
Nitrate standard solution, 1.00 mg/l NO₃⁻, Cat. No. 1.32240
Nitrate standard solution, 10.0 mg/l NO₃⁻, Cat. No. 1.32241
Nitrate standard solution, 50.0 mg/l NO₃⁻, Cat. No. 1.32242

Empty cells 16 mm with screw caps (25 pcs), Cat. No. 1.14724
Pipettes for pipetting volumes of 1.0, 1.5, and 5.0 ml
Rectangular cells 10 mm (2 pcs), Cat. No. 1.14946

6. Preparation

- Analyze immediately after sampling.
- Check the nitrite content with the MQuant® Nitrite Test.
If necessary, eliminate interfering nitrite ions (stated amounts apply for nitrite contents of up to 10 mg/l):
To 5 ml of sample add approx. 50 mg of amidosulfuric acid and dissolve.
The pH of this solution must be within the range 1 - 3. Adjust, if necessary, with sulfuric acid.
- Filter turbid samples.

7. Procedure

Reagent NO ₃ -1	5.0 ml	Pipette into a dry test tube ¹⁾ . Carefully allow from a pipette to run down the inside of the tilted test tube onto the reagent (Wear eye protection! The mixture becomes hot!). Immediately after the addition of the sample, carefully allow from a second pipette to run down the inside of the tilted test tube onto the solution (Wear eye protection! The mixture becomes hot!). Shake vigorously , holding only the upper part of the tube!
Pretreated sample (10 - 40 °C)	1.0 ml	
Reagent NO ₃ -2	1.5 ml	
Leave the hot reaction solution to stand for 15 min (reaction time A). Do not cool with cold water!		
Reagent NO ₃ -3	2 level grey microspoons (in the cap of the NO ₃ -3 bottle)	Add and shake vigorously until the reagent is completely dissolved .
Leave to stand for 60 min (reaction time B) , then fill the measurement sample into a 10-mm rectangular cell, and measure in the photometer.		

¹⁾ Empty cells Cat. No. 1.14724 are recommended that can be sealed with screw caps, thus enabling the sample to be mixed safely.

Notes on the measurement:

- Certain photometers may require a blank** (preparation as per measurement sample, but with distilled water instead of sample).
- For photometric measurement the cells must be clean.
Wipe, if necessary, with a clean dry cloth.
- Measurement of turbid solutions yields false-high readings.
- After the end of the reaction time B stated above, the color of the measurement solution remains virtually stable for 60 min. (The measurement value decreases by 3% during this time.)

8. Analytical quality assurance

it is recommended prior to each measurement series
To check the photometric measurement system (test reagents, measurement device, handling) and the mode of working, the nitrate standard solutions (see section 5) or Spectroquant® CombiCheck 20 can be used. Besides a **standard solution** with 9.0 mg/l NO₃-N, this article also contains an **addition solution** for determining sample-dependent interferences (**matrix effects**).
Additional notes see under www.qa-test-kits.com.
For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Notes

- Reclose the reagent bottles immediately after use.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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© 2024 Merck KGaA, Darmstadt, Alemania y/o sus filiales. Todos los derechos reservados. MilliporeSigma, Supelco, Sigma-Aldrich y Spectroquant son marcas comerciales de Merck KGaA, Darmstadt, Alemania. Todas las demás marcas comerciales son propiedad de sus respectivos propietarios. Tiene a su disposición información detallada sobre las marcas comerciales a través de recursos accesibles al público.