

## 51239 Oxytetra Selective Supplement (OGYE Selective Supplement)

Added to OGY Agar (Prod No 75310) or OGYE *ChromoSelect* Agar (Prod No 66481)  
Oxytetracycline suppresses the growth of accompanying bacterial flora during culturing yeasts and moulds.

Composition  
(per vial sufficient for 500 ml medium)  
Oxytetracycline 50.0 mg

Directions  
Aseptically add 5 ml of sterile distilled water to the content of the vial. Mix well and add to 500 ml of sterile, molten OGY Agar (Prod No 75310) or OGYE *ChromoSelect* Agar Base (Prod No 66481). Mix well and pour into sterile petri plates.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

