



QUALITY CONTROL RANGES
MILLIPLEX® MAP Human Kidney Injury Magnetic Panel 6
Catalog # HKI6MAG-99K

Control Catalog # HKI6-6099-6 Lot # HKI6-107 and HKI6-207

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
RBP4	Control 1	0.14 - 0.29	ng/mL
	Control 2	1.5 - 3.0	ng/mL
B-2-Microglobulin	Control 1	0.27 - 0.57	ng/mL
	Control 2	2.9 - 6.1	ng/mL
Clusterin	Control 1	17 - 35	ng/mL
	Control 2	167 - 348	ng/mL
Cystatin C	Control 1	0.29 - 0.60	ng/mL
	Control 2	3.5 – 7.4	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using assay buffer provided in the kit.

HKI6-6099-6: HKI6-107.207 (16 MAY 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.