

08449 Hemoglobin from bovine blood

An enrichment supplement, whose 2% aqueous solution is autoclavable. The solution contains flocculent dispersible precipitate. It is used for the preparation of Thayer Martin Medium.

contains hemoglobin powder
Loss on drying <10 %

Directions

Dissolve 20 g in 1 litre water and autoclave at 121 °C for 15 min. Cool to approx. 50 °C. Add 250 ml of this solution aseptically to 240 ml sterile, molten Thayer-Martin Medium Base (Cat. No. 88296)/Chocolate Agar Base along with the rehydrated contents of one bottle of each Thayer-Martin Supplement I and II (Cat. No. 44051 and 52083). Mix gently and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

