

15821 Klebsiella Selective Supplement

For selective isolation and detection of Klebsiella species.

Composition

(per vial, sufficient for 500 ml medium)

Carbenicillin 25 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten Klebsiella *ChromoSelect* Selective Agar Base (Cat. No. 90925). Mix well and pour into sterile petri plates

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

