

Separation of ingredients in analgesic agents

Plate	HPTLC pre-coated plate RP-18 WF 254s with concentrating zone (1.13124)
Eluent	Acetonitrile/water 50/50 (v/v)
Migration Distance	5 cm
Chamber	Normal chamber without chamber saturation
Compounds	1. Propyphenazone 0,1% 2. Phenazone 0,04% 3. 4-Aminobenzoic acid 0,03%
Application volume	5 µl
Detection	In-situ evaluation with TLC/HPTLC scanner (Camag) UV 265 nm

