

Product information

Anti-Actin, a-Smooth Muscle - FITC antibody, Mouse monoclonal

1

clone 1A4, purified from hybridoma cell culture

F3777

Product Description

Anti-Actin, a-Smooth Muscle-FITC is a purified mouse monoclonal antibody conjugated with fluorescein isothiocyanate (FITC) isomer I. Monoclonal Anti-Actin, a-Smooth Muscle (mouse IgG2a isotype) is derived from the 1A4 hybridoma produced by the fusion of mouse myeloma cells and splenocytes of immunized BALB/c mice. The NH₂ terminal synthetic decapeptide of a-smooth muscle actin coupled to keyhole limpet hemocyanin (KLH) was used as the immunogen.¹ The isotype is determined by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents, Cat. No. ISO2.

Anti-Actin, a-Smooth Muscle specifically recognizes the a-smooth muscle isoform of actin (42 kDa) by ELISA and immunoblotting. It does not react with the other major actin isoforms present in fibroblasts or epithelial cells (β and γ -cytoplasmic), striated muscle (a-sarcomeric), myocardium (a-myo-cardial), or y-smooth muscle isoform.

Anti-Actin, a-Smooth Muscle (also known as Anti-a-SM-1) recognizes the a-smooth muscle isoform of actin. The antibody reacts with normal and neoplastic, human vascular and visceral smooth muscle cells. It reacts with normal myoepithelial cells, pericytes, eye lens cells, hair follicle cells and certain stromal cells in the intestine, testis, lymphoid organs, liver, ovary and bone marrow. 1,2,3,4,5,6 The antibody also reacts with stromal myofibroblasts in hypertrophic scars, and in neoplastic tissues.7 a-Smooth muscle actin is transiently co-expressed with sarcomeric a-actin during myogenesis in chicken and rat embryos.^{8,9} It was also found in the ventricular conducting tract of adult mammalian heart. It is expressed in leiomyomas, leiomyosarcomas and leiomyoblastomas, as well as in a proportion of rhabdomyosarcomas. 10,11 The antibody cross reacts with actin expressing cells in human, bovine, goat, sheep, rabbit, cat, dog, mouse, rat, hamster, guinea pig, chicken, viper, lizard, frog, snail and crayfish tissues.

It can be used for staining acetone-fixed, frozen sections, smears, cytospins and EM preparations. The FITC conjugated monoclonal anti-a smooth muscle actin is especially useful for direct staining of tissues and cells.



Monoclonal Anti-Actin, a-Smooth Muscle-FITC may be used for:

- Identification of developing and adult smooth muscle pericytes and myoepithelial cells.
- Detection and characterization of smooth muscle tumors, glomus tumors and certain myoepithelial tumors, osteosarcomas and soft tissue tumors.
- Differentiation between glomus tumors and hemangiopericytomas, and between epitheliosis and intraductal breast carcinoma.
- 4. Studies on the expression of actins in cultured cells.
- Detection of a-smooth muscle actin positive cells in hepatic fibrosis, bone marrow fibrosis, experimental gliosis, atherosclerosis, pulmonary hypertension, and wound healing.

Reagents

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% BSA with 15 mM sodium azide as a preservative.

Precautions/Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight

turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Protect from prolonged exposure to light.

Product Profile

Immunohistochemistry: A minimum working dilution of 1:500 was determined using formalin-fixed, paraffin-embedded sections of human tonsil or appendix.

Note: In order to obtain best results in different techniques and preparations, it is recommended that each individual user determine their optimum working dilutions by titration assay.

References

- 1. Skalli, O., et al.., J. Cell Biol., 103, 2787 (1986).
- 2. Skalli, O., et al., Lab. Invest., 60, 275 (1989).
- 3. Skalli, O., et al., J. Histochem. Cytochem., 37, 315 (1989).
- 4. Peled, A., et al., Blood, 78, 304 (1991).
- 5. Schmitt-Graff, A., et al., Differentiation, 43, 115 (1990).
- 6. Jahoda, C., et al., J. Cell Sci., 99, 627 (1991).
- 7. Sappino, A., et al., Lab. Invest., 63, 144 (1990).
- Ruzicka, D., and Schwartz, R., J. Cell Biol., 107, 2575 (1988).
- 9. Babai, F., et al., Differentiation, 44, 132 (1990).
- 10. Jones, H., et al., J. Pathol., 162, 29 (1990).
- 11. Roholl, P., et al., Hu. Pathol., 21, 1269 (1990).



Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

