



QUALITY CONTROL RANGES
MILLIPLEX® MAP Canine Pituitary Magnetic Panel
Catalog # CANPIT-96K
Control Catalog # CANP-6096 Lot# CANP-105 and
CANP-7096 Lot# CANP-205

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
ACTH	CANP-6096	CANP-105	31 – 64	pg/mL
	CANP-7096	CANP-205	241 - 501	pg/mL
FSH	CANP-6096	CANP-105	28 - 59	ng/mL
	CANP-7096	CANP-205	242 - 503	ng/mL
GH	CANP-6096	CANP-105	0.15 - 0.31	ng/mL
	CANP-7096	CANP-205	1.3 - 2.6	ng/mL
TSH	CANP-6096	CANP-105	43 - 90	pg/mL
	CANP-7096	CANP-205	396 - 822	pg/mL
LH	CANP-6096	CANP-105	1.1 - 2.2	ng/mL
	CANP-7096	CANP-205	9.2 - 19	ng/mL
Prolactin	CANP-6096	CANP-105	15 - 32	ng/mL
	CANP-7096	CANP-205	126 - 261	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.

Quality Control values in culture media are not tested.

CANP-6096 Lot# CANP-105 and CANP-7096 Lot# CANP-205 (31 January 2025)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA