



QUALITY CONTROL RANGES

MILLIPLEX® MAP Human Kidney Injury Magnetic Bead Panel 1 Catalog # HKI1MAG-99K

Control Catalog # HKI1-6099-1 Lot # HKI1-105 and HKI1-205

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Calbindin	Control 1	1.4 – 2.9	ng/mL	IP-10	Control 1	0.030 – 0.062	ng/mL
	Control 2	14 – 29	ng/mL		Control 2	0.31 – 0.65	ng/mL
Osteoactivin	Control 1	0.10 – 0.58	ng/mL	Renin	Control 1	0.16 – 0.33	ng/mL
	Control 2	6.1 – 13	ng/mL		Control 2	1.5 – 3.2	ng/mL
GST α	Control 1	0.25 – 0.53	ng/mL	FABP-1	Control 1	44 – 91	ng/mL
	Control 2	2.7 – 5.5	ng/mL		Control 2	424 – 881	ng/mL
TIMP-1	Control 1	0.62 – 1.3	ng/mL	Collagen IV	Control 1	2.2 – 4.6	ng/mL
	Control 2	6.1 – 13	ng/mL		Control 2	23 – 48	ng/mL
KIM-1	Control 1	0.29 – 0.60	ng/mL	TFF-3	Control 1	3.0 – 6.3	ng/mL
	Control 2	3.1 – 6.5	ng/mL		Control 2	31 – 64	ng/mL

CAT# HKI1-6099-1, Lot# HKI1-105.205 (18 February 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.