

QC RANGE SHEET
Human Cardiovascular Disease Magnetic Bead Panel 4
Catalog # HCVD4MAG-67K
QC Catalog #: HCVMG4-6067-4



<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
sE-Selectin	Control 1	HCVM4-109	73 – 152	ng/mL
	Control 2	HCVM4-209	293 – 609	ng/mL
Follistatin	Control 1	HCVM4-109	0.95 – 2.0	ng/mL
	Control 2	HCVM4-209	3.9 – 8.1	ng/mL
Pecam-1	Control 1	HCVM4-109	6.5 – 14	ng/mL
	Control 2	HCVM4-209	24 – 49	ng/mL
Pentraxin-3	Control 1	HCVM4-109	3.6 – 7.4	ng/mL
	Control 2	HCVM4-209	14 – 30	ng/mL
Tissue Factor	Control 1	HCVM4-109	4.0 – 8.3	ng/mL
	Control 2	HCVM4-209	17 – 35	ng/mL
Thrombomodulin	Control 1	HCVM4-109	8.1 – 17	ng/mL
	Control 2	HCVM4-209	33 – 68	ng/mL
Troponin T	Control 1	HCVM4-109	0.15 – 0.31	ng/mL
	Control 2	HCVM4-209	0.63 – 1.3	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.
Quality Control values in culture media are not tested.

HCVMG4-6067-4: HCVM4-109.209 (17 November 2020)
 EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com