

Product Information

SIGMA QUALITY CONTROL TEST PROCEDURE

Enzymatic Assay of α -GALACTOSIDASE
(EC 3.2.1.22)
Sigma Prod. No. G-8507

PRINCIPLE:

PNP α -D-Galactopyranoside + H₂O $\xrightarrow{\alpha\text{-Galactosidase}}$ p-Nitrophenol + D-Galactose

Abbreviation used:

PNP α -D-Galactopyranoside = p-Nitrophenyl α -D-Galactopyranoside

CONDITIONS: T = 25°C, pH = 6.5, A_{405nm}, Light path = 1 cm

METHOD: Spectrophotometric Stop Rate Determination

REAGENTS:

- A. 100 mM Potassium Phosphate Monobasic Solution
(Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No. P-5379.)
- B. 100 mM Potassium Phosphate Dibasic Solution
(Prepare 100 ml in deionized water using Potassium Phosphate, Dibasic, Trihydrate, Sigma Prod. No. P-5504.)
- C. 100 mM Potassium Phosphate Buffer, pH 6.5 at 25°C
(Prepare 100 ml by adjusting 50 ml of Reagent A to pH 6.5 at 25°C by adding Reagent B.)
- D. 9.9 mM p-Nitrophenyl α -D-Galactopyranoside Solution (PNP-Gal)
(Prepare 4 ml in deionized water using p-Nitrophenyl α -D-Galactopyranoside, Sigma Prod. No. N-0877.)
- E. 200 mM Borate Buffer, pH 9.8 at 25°C
(Prepare 100 ml in deionized water using Boric Acid, Sigma Prod. No. B-0252. Adjust to pH 9.8 at 25°C with 1 M NaOH.)
- F. α -Galactosidase Enzyme Solution
(Immediately before use, prepare a solution containing 0.05 - 0.10 units/ml of α -Galactosidase in cold Reagent C.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Reagent C (Potassium Phosphate Buffer)	0.70	0.70
Reagent D (PNP-Gal)	0.20	0.20

Mix by swirling and equilibrate to 25°C. Then add:

Reagent F (Enzyme Solution)	0.10	-----
-----------------------------	------	-------

Immediately mix by swirling and incubate at 25°C for exactly 5 minutes. Then add:

Reagent E (Borate Buffer)	2.00	2.00
Reagent F (Enzyme Solution)	-----	0.10

Mix by swirling and record the $A_{405\text{nm}}$ for both the Test and Blank, using a suitably thermostatted spectrophotometer.

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(A_{405\text{nm}} \text{ Test} - A_{405\text{nm}} \text{ Blank})(3.0)(\text{df})}{(18.5)(5)(0.1)}$$

3.0 = Total volume of assay

df = Dilution factor

5 = Conversion factor for 5 minutes to 1 minute

18.5 = Millimolar extinction coefficient of p-Nitrophenol at 405 nm

0.1 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will hydrolyze 1.0 μmole of p-nitrophenyl α -D-galactoside to p-nitrophenol and D-galactose per minute at pH 6.5 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 1.00 ml reaction mix, the final concentrations are
80 mM potassium phosphate, 2.0 mM p-nitrophenyl α -D-galctopyranoside, and 0.005 - 0.01 units
 α -galactosidase.

REFERENCES:

(1968) *Eur. J. Biochem.* **8**, 395

NOTES:

1. This assay is based on the cited reference.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.