

51387 MRSA Selective Supplement

An antimicrobial supplement recommended for the selective isolation of Methicillin Resistant *Staphylococcus aureus* (MRSA) from clinical specimens.

Composition
(per vial, sufficient for 500 ml medium)
Methicillin 2 mg

Directions
Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten MeReSa ChromoSelect Agar Base (Cat. No. 90923), cooled to 45-50°C. Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

