

## Technical Data Sheet

### Anaerobic Jar

Ordering number: 1.13681.0001

#### Mode of Action

For cultivation of anaerobic and microaerophilic microorganisms in defined atmospheric conditions

- using Anaerocult® A (Cat. No. 1.32381.0001) and Anaerotest® (Cat. No. 1.32371.0001) for anaerobic microorganisms
- using Anaerocult® C (Cat. No. 1.32383.0001) for microaerophilic microorganisms.

The anaerobic jar is to be used together with Cat. No. 1.07040.0001 Petri dish rack (for up to 12 Petri dishes).

#### Ordering Information

| Product       | Cat. No.     |
|---------------|--------------|
| Anaerobic Jar | 1.13681.0001 |