

08582 Urea Solution 40%

Filter sterilized urea solution recommended for detection of urease activity.

Composition

(1 vial is sufficient for 100 ml medium)

Urea	2.0g
Distilled water	5.0ml

Directions

Warm up the refrigerated urea solution to room temperature and aseptically add 5 ml in 95 ml sterile molten Urea Agar Base (Prod No U1757) or Christensen's Urea Agar (Prod No 27048). Mix well and dispense in sterile tubes.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

