

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

SILu™Lite CRP, C-reactive protein, human recombinant, expressed in HEK cells MS protein standard

Catalog Number MSST0022 Storage Temperature -20 °C

Synonyms: pentraxin, PTX1

Product Description

SILu™Lite CRP is a recombinant human protein expressed in human 293 cells. It is a monomer consisting of 226 amino acids (including a C-terminal polyhistidine and FLAG® tags), with a calculated molecular mass of 25.4 kDa. SILu™Lite CRP is an analytical standard designed to be used as starting material for preparation of calibrators and controls in LC-MS applications.

CRP is a hepatically derived pentraxin. It exists mainly as a pentamer in a cyclic structure with 5 identical 23 kDa subunits in the blood. In addition to its role as an "acute phase" protein or first-line defense molecule against pathogenic organisms,2 baseline levels of CRP correlate positively to adverse cardiovascular events.2 CRP is, therefore, considered a "risk factor" for cardiovascular disease.3 In humans, treatment with statin therapy reduces levels of both LDL-C and CRP, and concurrently there is a reduction in the number of cardiovascular events.^{4,5} A recent study suggests patients with elevated basal levels of CRP are at an increased risk of diabetes, in addition to cardiovascular disease. 6 A study of over 700 nurses showed that those in the highest quartile of trans fat consumption had blood levels of CRP that were 73% higher than those in the lowest quartile.

Each vial contains 50-65 μg of SILu™Lite CRP standard, in a 0.1 mg/mL solution of 20 mM sodium phosphate, pH 8.0, 1 M NaCl, 1 mM EDTA, and 25% glycerol. Vial content was determined by the Bradford method using BSA as a calibrator. The correction factor from the Bradford method to Amino Acid Analysis is 130% for this protein.

Identity: Confirmed by peptide mapping

Purity: ≥95% (SDS-PAGE)

UniProt: P02741

Sequence Information

The C-terminal polyhistidine and FLAG tags are italicized.

QTDMSRK**AFVFPKESDTSYVSLK**APLTKPLKAFTVCL HFYTELSSTR**GYSIFSYATK**RQDNEILIFWSKDIGYSF TVGGSEILFEVPEVTVAPVHICTSWESASGIVEFWVD GKPRVRKSLKKGYTVGAEASIILGQEQDSFGGNFEGS QSLVGDIGNVNMWDFVLSPDEINTIYLGGPFSPNVLN WRALKYEVQGEVFTKPQLWPDYKDDDDKGHHHHHH HHGGQ

Transitions for three peptides (underlined) suggested for selected reaction monitoring analysis (SRM) are provided for download on the product display page at www.sigmaaldrich.com.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the product at -20 °C. The product is stable for at least 2 years as supplied. After initial thawing it is recommended to store the protein in working aliquots at -20 °C.

References

- Pepys, M.B., and Hirschfield, G.M., C-reactive protein: a critical update. *J. Clin. Invest.*, 111, 1805–1812 (2003).
- 2. Young, B. et al., C-reactive protein: A critical review. *Pathology*, **23**, 118-124 (1991).
- Tracy, R.P. et al., Lifetime smoking exposure affects the association of Creactive protein with cardiovascular disease risk factors and subclinical disease in elderly subjects. *Arterioscler. Thromb. Vasc. Biol.*, 17, 2167-2176 (1997).
- 4. Kinlay, S. Low-density lipoprotein-dependent and -independent effects of cholesterol-lowering therapies on C-reactive protein: a meta-analysis. *J. Am. Coll. Cardiol.*, **49** 2003–2009 (2007).

- 5. Nissen, S.E. et al., Statin therapy, LDL cholesterol, C-reactive protein, and coronary artery disease. *N. Engl. J. Med.*, **352**, 29–38 (2005).
- 6. Pradhan, A.D. et al., C-Reactive Protein, Interleukin 6, and Risk of Developing Type 2 Diabetes Mellitus. *J. Am. Med. Assoc.*, **286**, 327-334 (2001).
- 7. Lopez-Garcia, E. et al., Consumption of trans fatty acids is related to plasma biomarkers of inflammation and endothelial dysfunction. *J. Nutr.*, **135**, 562–566 (2005).

FLAG is a registered trademark and SILu is a trademark of Sigma-Aldrich Co. LLC.

NA, AI, MAM, KR 02/17-1