

90577 Bifido Selective Supplement B (Acetic acid Supplement)

A selective supplement for the selective isolation of Bifidobacterium species by colony count technique from milk products.

Composition

(per vial, sufficient for 500 ml medium)

Glacial acetic acid (1%) 5 ml

Directions

Aseptically add contents of one vial in 500 ml sterile, molten, cooled (45-50°C)

TOS-propionate agar medium (acc. ISO) (Cat. No. 43314) along with 1 vial of Mupirocin Supplement (Cat. No. 89914). Mix well and pour into sterile Petri plates.

Note:

After addition of Bifido Selective Supplement B (Cat. No. 90577) pH of the medium will shift to the acidic side, which does not affect the performance of the medium.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

