



QUALITY CONTROL RANGES
MILLIPLEX® MAP Human Cardiovascular Disease Magnetic Bead Panel 2
Kit Catalog# HCVD2MAG-67K
Control Catalog # HCVMG2-6067-2 Lot# HCV2-112.212

Analyte	QC Level	Expected Range	Units	Analyte	QC Level	Expected Range	Units
ADMTS13	Control 1	15 - 31	ng/mL	NGAL	Control 1	0.17 - 0.36	ng/mL
	Control 2	121 - 252	ng/mL		Control 2	1.3 - 2.6	ng/mL
D-dimer	Control 1	3.0 - 6.2	ng/mL	sVCAM-1	Control 1	1.6 - 3.4	ng/mL
	Control 2	24 - 50	ng/mL		Control 2	13 - 26	ng/mL
GDF-15	Control 1	0.007 - 0.02	ng/mL	SAA	Control 1	2.3 - 4.7	ng/mL
	Control 2	0.06 - 0.13	ng/mL		Control 2	19 - 39	ng/mL
Myoglobin	Control 1	0.29 - 0.60	ng/mL				
	Control 2	2.3 - 4.8	ng/mL				
sICAM-1	Control 1	1.1 - 2.3	ng/mL				
	Control 2	9.0 - 19	ng/mL				
MPO	Control 1	0.25 - 0.53	ng/mL				
	Control 2	2.4 - 5.0	ng/mL				
P-selectin	Control 1	2.7 - 5.6	ng/mL				
	Control 2	24 - 50	ng/mL				

HCVMG2-6067-2: HCV2-112.212 (18 October 2023)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA