

QUALITY CONTROL RANGES



MILLIPLEX[®] MAP Rat Kidney Toxicity Magnetic Bead Panel 1

Catalog # RKTXMAG-37K

Control Catalog # RKTXMG1-6037-1

Lot: RKTX1-108 and RKTX1-208

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Clusterin	Control I	38 - 80	ng/mL
	Control II	193 - 400	ng/mL
GSTα	Control I	61 - 190	ng/mL
	Control II	422 - 876	ng/mL
IP-10	Control I	0.54 - 1.1	ng/mL
	Control II	2.1 - 4.4	ng/mL
KIM-1	Control I	0.49 - 1.0	ng/mL
	Control II	2.5 - 5.2	ng/mL
OPN	Control I	0.49 - 1.0	ng/mL
	Control II	2.5 - 5.3	ng/mL
TIMP-1	Control I	5.5 - 11	ng/mL
	Control II	27 - 56	ng/mL
VEGF	Control I	0.51 - 1.1	ng/mL
	Control II	2.6 - 5.4	ng/mL

Note: The Quality Control Ranges are generated with using serum matrix provided in the kit. Quality Control values in culture media are not tested.

RKTXMG1-6037-1: RKTX1-108.208 (12 May 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.