

Product Information

Dulbecco's Modified Eagle's Medium (DME)/Ham's Nutrient Mixture F-12

D0547, D0697, D2906, D6421, D6434, D8900, D8062, D8437, D9785

The use of serum-free media has expanded as researchers investigate the conditions for successful culture of a variety of cell lines. Instead of serum, the media are supplemented with a defined combination of nutrients, growth factors, and hormones. Variation of the supplements and their concentrations are specific for the type of cell under study.

The media often used for these defined studies contain a 1:1 mixture of Dulbecco's Modified Eagle's Medium (DME) and Ham's F-12 Nutrient Mixture. HEPES buffer may be included in the formulations at a final concentration of 15 mM to compensate for the loss of buffering capacity incurred by eliminating serum.

COMPONENT	D0547 [powder] g/L	D0697 [1x] g/L	D2906 [powder] g/L	D6421 [1x] g/L	D6434 [1x] g/L	D8900 [powder] g/L	D8062 [1x] g/L	D8437 [1x] g/L	D9785 [powder] g/L
Inorganic Salts									
(NH ₄) ₆ Mo ₇ O ₂₄ • 4H ₂ O	0.00000618	—	—	—	—	—	—	—	—
NH ₄ BO ₃	0.00000058	—	—	—	—	—	—	—	—
CaCl ₂	0.1166	—	0.1166	—	—	0.1166	—	—	—
CaCl ₂ • 2H ₂ O	—	0.1545	—	0.1545	0.1545	—	0.1545	0.1545	—
CuSO ₄ • 5H ₂ O	0.0000013	0.0000013	0.0000013	0.0000013	0.0000013	0.0000013	0.0000013	0.0000013	0.0000013
Fe(NO ₃) ₃ • 9H ₂ O	0.00005	0.00005	0.00005	0.00005	0.00005	0.00005	0.00005	0.00005	0.00005
FeSO ₄ • 7H ₂ O	0.000417	0.000417	0.000417	0.000417	0.000417	0.000417	0.000417	0.000417	0.000417
MgCl ₂ • 6H ₂ O	0.0612	0.0612	0.0612	0.0612	0.0612	0.0612	0.0612	0.0612	—
MgSO ₄	0.04884	0.04884	0.04884	0.04884	0.04884	0.04884	0.04884	0.04884	—
MnSO ₄	0.000000151	—	—	—	—	—	—	—	—
NiCl ₂	0.00000012	—	—	—	—	—	—	—	—
KCl	0.3118	0.3118	0.3118	0.3118	0.3118	0.3118	0.3118	0.3118	0.3118
NaHCO ₃	—	1.2	—	1.2	1.2	—	1.2	1.2	—
NaCl	6.996	6.996	6.996	6.996	6.996	6.996	6.996	6.996	6.996
NaSiO ₃ • 9H ₂ O	0.0000142	—	—	—	—	—	—	—	—
Na ₂ SeO ₃	0.00000519	—	—	—	—	—	—	—	—
Na ₂ HPO ₄	0.07102	0.07102	0.07102	0.07102	0.07102	0.07102	0.07102	0.07102	0.07102
NaH ₂ PO ₄	0.0543	0.0543	0.0543	0.0543	0.0543	0.0543	0.0543	0.0543	0.0543

COMPONENT	D0547 [powder] g/L	D0697 [1x] g/L	D2906 [powder] g/L	D6421 [1x] g/L	D6434 [1x] g/L	D8900 [powder] g/L	D8062 [1x] g/L	D8437 [1x] g/L	D9785 [powder] g/L
SnCl ₂ • 7H ₂ O	0.00000011	—	—	—	—	—	—	—	—
ZnSO ₄ • 7H ₂ O	0.000432	0.000432	0.000432	0.000432	0.000432	0.000432	0.000432	0.000432	0.000432
Amino Acids									
L-Alanine	0.00445	0.0045	0.00445	0.0045	0.00445	0.00445	0.00445	0.0045	0.00445
L-Alanyl-L-Glutamine	—	0.54253	—	—	—	—	—	—	—
L-Arginine • HCl	0.1475	0.1475	0.1475	0.1475	0.1475	0.1475	0.1475	0.1475	0.1475
L-Asparagine • H ₂ O	0.0075	0.0075	0.0075	0.0075	0.0075	0.0075	0.0075	0.0075	0.0075
L-Aspartic Acid	0.00665	0.00665	0.00665	0.00665	0.00665	0.00665	0.00665	0.00665	0.00665
L-Cystine • 2HCl	0.03129	0.03129	0.03129	0.03129	0.03129	0.03129	0.03129	0.03129	0.03129
L-Cysteine • HCl • H ₂ O	0.01756	0.01756	0.01756	0.01756	0.01756	0.01756	0.01756	0.01756	0.01756
L-Glutamic Acid	0.00735	0.00735	0.00735	0.00735	0.00735	0.00735	0.00735	0.00735	0.00735
L-Glutamine	0.365	—	0.365	—	—	0.365	0.365	0.365	—
Glycine	0.01875	0.01875	0.01875	0.01875	0.01875	0.01875	0.01875	0.01875	0.01875
L-Histidine • HCl • H ₂ O	0.03148	0.03148	0.03148	0.03148	0.03148	0.03148	0.03148	0.03148	0.03148
L-Isoleucine	0.05447	0.05447	0.05447	0.05447	0.05447	0.05447	0.05447	0.05447	0.05447
L-Leucine	0.05905	0.05905	0.05905	0.05905	0.05905	0.05905	0.05905	0.05905	—
L-Lysine • HCl	0.09125	0.09125	0.09125	0.09125	0.09125	0.09125	0.09125	0.09125	—
L-Methionine	0.01724	0.01724	0.01724	0.01724	0.01724	0.01724	0.01724	0.01724	—
L-Phenylalanine	0.03548	0.03548	0.03548	0.03548	0.03548	0.03548	0.03548	0.03548	0.03548
L-Proline	0.01725	0.01725	0.01725	0.01725	0.01725	0.01725	0.01725	0.01725	0.01725
L-Serine	0.02625	0.02625	0.02625	0.02625	0.02625	0.02625	0.02625	0.02625	0.02625
L-Threonine	0.05345	0.05345	0.05345	0.05345	0.05345	0.05345	0.05345	0.05345	0.05345
L-Tryptophan	0.00902	0.00902	0.00902	0.00902	0.00902	0.00902	0.00902	0.00902	0.00902
L-Tyrosine • 2Na • 2H ₂ O	0.05579	0.05579	0.05579	0.05579	0.05579	0.05579	0.05579	0.05579	0.05579
L-Valine	0.05285	0.05285	0.05285	0.05285	0.05285	0.05285	0.05285	0.05285	0.05285

COMPONENT	D0547 [powder] g/L	D0697 [1x] g/L	D2906 [powder] g/L	D6421 [1x] g/L	D6434 [1x] g/L	D8900 [powder] g/L	D8062 [1x] g/L	D8437 [1x] g/L	D9785 [powder] g/L
Vitamins									
D-Biotin	0.0000035	0.0000035	0.0000035	0.0000035	0.0000035	0.0000035	0.0000035	0.0000035	0.0000035
Choline Chloride	0.00898	0.00898	0.00898	0.00898	0.00898	0.00898	0.00898	0.00898	0.00898
Folic Acid	0.00266	0.00266	0.00266	0.00265	0.00266	0.00266	0.00266	0.00266	0.00266
myo-Inositol	0.0126	0.0126	0.0126	0.0126	0.0126	0.0126	0.0126	0.0126	0.0126
Niacinamide	0.00202	0.00202	0.00202	0.00202	0.00202	0.00202	0.00202	0.00202	0.00202
D-Pantothenic Acid • ½Ca	0.00224	0.00224	0.00224	0.00224	0.00224	0.00224	0.00224	0.00224	0.00224
Pyridoxal • HCl	0.002	—	0.002	—	—	0.002	0.002	—	0.002
Pyridoxine • HCl	0.000031	0.002031	0.000031	0.002031	0.002031	0.000031	0.002031	0.002031	0.000031
Riboflavin	0.000219	0.000219	0.000219	0.000219	0.000219	0.000219	0.000219	0.000219	0.000219
Thiamine • HCl	0.00217	0.00217	0.00217	0.00217	0.00217	0.00217	0.00217	0.00217	0.00217
Vitamin B ₁₂	0.00068	0.00068	0.00068	0.00068	0.00068	0.00068	0.00068	0.00068	0.00068
Other									
D-Glucose	3.15	3.15	3.15	3.15	3.15	3.15	3.15	3.15	3.15
HEPES	—	3.5745	3.5745	3.5745	3.5745	3.5745	—	3.5745	3.5745
Hypoxanthine	0.0021	0.00244	0.0021	0.00244	0.00244	0.0021	0.00244	0.00244	0.0021
Linoleic Acid	0.000042	0.000042	0.000042	0.000042	0.000042	0.000042	0.000042	0.000042	0.000042
Phenol Red • Na	0.00863	0.00863	—	0.00863	—	0.00863	0.00863	0.00863	—
Putrescine • 2HCl	0.000081	0.000081	0.000081	0.000081	0.000081	0.000081	0.000081	0.000081	0.000081
Pyruvic Acid • Na	0.11	0.055	0.055	0.055	0.055	0.055	0.055	0.055	0.055
DL-Thioctic Acid	0.000105	0.000105	0.000105	0.000105	0.000105	0.000105	0.000105	0.000105	0.000105
Thymidine	0.000365	0.000365	0.000365	0.000365	0.000365	0.000365	0.000365	0.000365	0.000365
ADD									
NaHCO ₃	1.2	—	1.2	—	—	1.2	—	—	1.2
CaCl ₂ • 2H ₂ O	—	—	—	—	—	—	—	—	0.1545
L-Glutamine	—	—	—	0.365	0.365	—	—	—	0.365
L-Leucine	—	—	—	—	—	—	—	—	0.05905
L-Lysine • HCl	—	—	—	—	—	—	—	—	0.09125
MgCl ₂ • 6H ₂ O	—	—	—	—	—	—	—	—	0.0612
MgSO ₄	—	—	—	—	—	—	—	—	0.04884
L-Methionine	—	—	—	—	—	—	—	—	0.0172

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2024-2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

D0547pis Rev 06/25

**Millipore
Sigma**