



QUALITY CONTROL RANGES
MILLIPLEX® MAP Canine Pituitary Magnetic Panel
Catalog # CANPIT-96K
Control Catalog # CANP-6096 Lot# CANP-104 and
CANP-7096 Lot# CANP-204

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
ACTH	CANP-6096	CANP-104	31 - 65	pg/mL
	CANP-7096	CANP-204	249 - 516	pg/mL
FSH	CANP-6096	CANP-104	32 - 66	ng/mL
	CANP-7096	CANP-204	249 - 517	ng/mL
GH	CANP-6096	CANP-104	0.17 - 0.35	ng/mL
	CANP-7096	CANP-204	1.3 - 2.7	ng/mL
TSH	CANP-6096	CANP-104	46 - 96	pg/mL
	CANP-7096	CANP-204	413 - 858	pg/mL
LH	CANP-6096	CANP-104	1.1 - 2.4	ng/mL
	CANP-7096	CANP-204	9.3 - 19	ng/mL
Prolactin	CANP-6096	CANP-104	16 - 34	ng/mL
	CANP-7096	CANP-204	133 - 277	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.

Quality Control values in culture media are not tested.

CANP-6096 Lot# CANP-104 and CANP-7096 Lot# CANP-204 (06 August 2024)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA