

XXXXXX

7.76118.0006 01/2025

Supelco

sigmaaldrich.com

50 Tests

1.17952


Reflectoquant®

**Hydroxymethylfurfural (HMF) Test
for RQflex® 20****Hydroxymethylfurfural (HMF)-Test
für RQflex® 20**

Store cold and dry • (2 - 8 °C) • Kalt und trocken lagern

1.0 - 60.0 mg/l

Reaction time: 120 sec



Please consult Instruction for Use for further information and observe the RQflex® manual.



For IfU please scan code with My M Safety App or go to sigmaaldrich.com



©2025

Lot

min. shelf life (YYYY/MM/DD)

Merck KGaA, 64271 Darmstadt, Germany
EMD Millipore Corporation, Burlington MA 01803, USA
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
Ontario, L6H 6J8, Canada