

Technical Bulletin

Anti-D₁ Dopamine Receptor Antibody, Rat Monoclonal

Clone 1-1-F11 s.E6, purified from hybridoma cell culture

D2944

Product Description

Anti-D₁ Dopamine Receptor (rat IgG2a isotype) is derived from the rat hybridoma 1-1-F11 S.E6 produced by the fusion of mouse myeloma cells and splenocytes from rat immunized with recombinant fusion protein containing the C-terminal 97 amino acid of human D₁ dopamine receptor (Gene ID: 1812).¹

Anti-D₁ Dopamine Receptor recognizes human, monkey, and rat D₁ Dopamine Receptor. The antibody may be used in various immunochemical techniques including immunohistochemistry, immunoblotting (~ 50 kDa) and immunocytochemistry.²⁻⁵

Dopamine receptors belong to the family of seven transmembrane domain G protein coupled receptors. They were initially divided into two general categories on the basis of differences in receptor pharmacology and biochemical mechanisms of signal transduction. With the application of the techniques of molecular biology, two predominant dopamine receptors, D₁ and D₂, were cloned. Later other dopamine receptors with homology to either the D₁ or D₂ receptor were identified. Thus, at present, two families of vertebrate dopamine receptors (designated as D₁-like and D₂-like) are recognized. The D₁-like family consists of the D₁ and D₅ receptors, generally associated to a stimulatory function, while the D₂-like family consists of the D₂, D₃ and D₄ receptors, generally associated to an inhibitory function.⁶ Dopamine receptors are mainly localized in the striatum, limbic system, the brain cortex and the infundibulum. However, the presence of the dopamine receptors has been demonstrated in most areas of the central nervous system.⁶ The D₁ and D₂ receptors occur in sufficiently high concentrations that they can be studied in situ, while D₃, D₄, and D₅ receptors occur in such low concentrations that study of them in situ is difficult. Thus, the majority of study of these receptors has been accomplished using cell lines cloned to express these receptors.

Reagent

The product is supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

Antibody concentration: ~ 1 mg/mL

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For extended storage, freeze at -20 °C in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

Immunocytochemistry: A working concentration of 5-10 µg/mL is recommended using D₁ Dopamine receptor transfected HEK-293T cells.

Note: In order to obtain the best results in various techniques and preparations, we recommend determining optimal working concentration by titration.

References

1. Levey, A.L., et al., Proc. Natl. Acad. Sci. USA, 90, 8861-8865 (1993).
2. Hersch, S.M., et al., J. Neurosci., 15, 5222-5237 (1995).
3. Smiley, J.F., et al., Proc. Natl. Acad. Sci USA, 91, 5720-5724 (1994).
4. Yung, K.K.L., et al., Neuroscience, 65, 709-730 (1995).
5. Fiorentini, C., et al., J. Biol. Chem., 278, 20196-20202 (2003).
6. Pivonello, R., et al., Eur. J. Endocrinol., 156, S13-S21 (2007).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
D2944dat Rev07/22

